

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов измерительной системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадает во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

1. **BS-120**: BS-120, BS-130;

2. **BS-200**: BS-200, BS-220;

3. **BS-200E**: BS-200E, BS-220E;

4. **BS-230**: BS-230, BS-240;

5. **BS-240E**: BS-240E, BS-240Pro;

6. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), **BS-350E(V35.00)** (Serial Number starts with "W9-" and software version starts with "35.00");

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

7. **BS-380**: BS-380, BS-390;

8. **BS-400**: BS-400, BS-420;

9. **BS-480**: BS-480, BS-490;

10. **BS-600**: BS-600, BS-620;

11. **BS-600M**: BS-600M, BS-620M;

12. **BS-800**: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

13. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

14. **BS-2800M**: BS-2600M, BS-2800M.

LOT: 057523004

EXP: 2024-12-12

Abbreviated name	G6PD Control (L)						G6PD Control (H)							
	Model	Unit	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value±3SD)		1 SD		
G6PD	BS-120 ¹	U/L	1124	899	—	1349	75	BS-120 ¹	U/L	2329	1861	—	2797	156
		µkat/L	18.8	15.0	—	22.5	1.3		µkat/L	38.9	31.1	—	46.7	2.6
	BS-200 ²	U/L	1137	909	—	1365	76	BS-200 ²	U/L	2348	1877	—	2819	157
		µkat/L	19.0	15.2	—	22.8	1.3		µkat/L	39.2	31.3	—	47.1	2.6
	BS-200E ³	U/L	1026	819	—	1233	69	BS-200E ³	U/L	2127	1698	—	2556	143
		µkat/L	17.1	13.7	—	20.6	1.2		µkat/L	35.5	28.4	—	42.7	2.4
	BS-230 ⁴	U/L	1059	846	—	1272	71	BS-230 ⁴	U/L	2206	1762	—	2650	148
		µkat/L	17.7	14.1	—	21.2	1.2		µkat/L	36.8	29.4	—	44.3	2.5
	BS-240E ⁵	U/L	1112	887	—	1337	75	BS-240E ⁵	U/L	2315	1850	—	2780	155
		µkat/L	18.6	14.8	—	22.3	1.3		µkat/L	38.7	30.9	—	46.4	2.6
	BS-360E ⁶	U/L	1043	833	—	1253	70	BS-360E ⁶	U/L	2161	1726	—	2596	145
		µkat/L	17.4	13.9	—	20.9	1.2		µkat/L	36.1	28.8	—	43.4	2.4
	BS-380 ⁷	U/L	1062	849	—	1275	71	BS-380 ⁷	U/L	2221	1774	—	2668	149
		µkat/L	17.7	14.2	—	21.3	1.2		µkat/L	37.1	29.6	—	44.6	2.5
	BS-400 ⁸	U/L	1044	834	—	1254	70	BS-400 ⁸	U/L	2196	1755	—	2637	147
		µkat/L	17.4	13.9	—	20.9	1.2		µkat/L	36.7	29.3	—	44.0	2.5
	BS-480 ⁹	U/L	1097	878	—	1316	73	BS-480 ⁹	U/L	2254	1801	—	2707	151
		µkat/L	18.3	14.7	—	22.0	1.2		µkat/L	37.6	30.1	—	45.2	2.5

Abbreviated name	G6PD Control (L)						G6PD Control (H)							
	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD		
G6PD	BS-600 ¹⁰	U/L	1097	878	—	1316	73	BS-600 ¹⁰	U/L	2254	1801	—	2707	151
		μ kat/L	18.3	14.7	—	22.0	1.2		μ kat/L	37.6	30.1	—	45.2	2.5
	BS-600M ¹¹	U/L	1102	880	—	1324	74	BS-600M ¹¹	U/L	2290	1831	—	2749	153
		μ kat/L	18.4	14.7	—	22.1	1.2		μ kat/L	38.2	30.6	—	45.9	2.6
	BS-800 ¹²	U/L	1097	878	—	1316	73	BS-800 ¹²	U/L	2254	1801	—	2707	151
		μ kat/L	18.3	14.7	—	22.0	1.2		μ kat/L	37.6	30.1	—	45.2	2.5
	BS-2000 ¹³	U/L	1137	909	—	1365	76	BS-2000 ¹³	U/L	2301	1839	—	2763	154
		μ kat/L	19.0	15.2	—	22.8	1.3		μ kat/L	38.4	30.7	—	46.1	2.6
	BS-2800M ¹⁴	U/L	1139	911	—	1367	76	BS-2800M ¹⁴	U/L	2320	1855	—	2785	155
		μ kat/L	19.0	15.2	—	22.8	1.3		μ kat/L	38.7	31.0	—	46.5	2.6

English	Abbreviated name	G6PD Control	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)
Русский	сокращенное наименование	Контроль G6PD	модель	Прибор	Результат анализа	Диапазон(результат анализа \pm 3CO)
Português	Nome abreviado	Controle G6PD	Modelo	Unidade	Valores da análise	Faixa (Valores da análise \pm 3SD)
Español	nombre abreviado	Control G6PD	modelo	Unidad	Valor de ensayo	Rango (Valor de ensayo \pm 3SD)
Italiano	abbreviazione	Controllo G6PD	modelli	Unità	Valori di dosaggio	Intervallo(valore di concentrazione \pm 3 DS)
Türkçe	kısaltılmış ad	G6PD Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri \pm 3SD)

G6PD	
English	Glucose-6-Phosphate Dehydrogenase
Русский	глюкозо-6-фосфатдегидрогеназа
Português	Glicose-6-fosfato desidrogenase
Español	Glucosa-6-fosfato deshidrogenasa
Italiano	Glucosio-6-fosfato deidrogenasi
Türkçe	Glikoz 6 Fosfat Dehidrogenaz