

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов измерительной системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадает во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1. **BS-120**: BS-120, BS-130;

7. **BS-380**: BS-380, BS-390;

2. **BS-200**: BS-200, BS-220;

8. **BS-400**: BS-400, BS-420;

3. **BS-200E**: BS-200E, BS-220E;

9. **BS-480**: BS-480, BS-490;

4. **BS-230**: BS-230, BS-240;

10. **BS-600**: BS-600, BS-620;

5. **BS-240E**: BS-240E, BS-240Pro;

11. **BS-600M**: BS-600M, BS-620M;

6. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

12. **BS-800**: BS-800, BS-820, BS-800M, BS-820M,

BS-1800, BS-1800plus;

BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"),
BS-350E(V35.00)
(Serial Number starts with "W9-" and software version starts with "35.00");

13. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

14. **BS-2800M**: BS-2600M, BS-2800M.

LOT: 057523005

EXP: 2025-01-11

Abbreviated name	G6PD Control (L)							G6PD Control (H)						
	Model	Unit	Assay Value	Range (Assay Value±3SD)			1 SD	Model	Unit	Assay Value	Range (Assay Value±3SD)			1 SD
G6PD	BS-120 ¹	U/L	1144	913	—	1375	77	BS-120 ¹	U/L	2382	1902	—	2862	160
		µkat/L	19.1	15.2	—	23.0	1.3		µkat/L	39.8	31.8	—	47.8	2.7
	BS-200 ²	U/L	1183	946	—	1420	79	BS-200 ²	U/L	2455	1963	—	2947	164
		µkat/L	19.8	15.8	—	23.7	1.3		µkat/L	41.0	32.8	—	49.2	2.7
	BS-200E ³	U/L	1014	810	—	1218	68	BS-200E ³	U/L	2123	1697	—	2549	142
		µkat/L	16.9	13.5	—	20.3	1.1		µkat/L	35.5	28.3	—	42.6	2.4
	BS-230 ⁴	U/L	1061	848	—	1274	71	BS-230 ⁴	U/L	2211	1767	—	2655	148
		µkat/L	17.7	14.2	—	21.3	1.2		µkat/L	36.9	29.5	—	44.3	2.5
	BS-240E ⁵	U/L	1113	888	—	1338	75	BS-240E ⁵	U/L	2324	1856	—	2792	156
		µkat/L	18.6	14.8	—	22.3	1.3		µkat/L	38.8	31.0	—	46.6	2.6
	BS-360E ⁶	U/L	1024	817	—	1231	69	BS-360E ⁶	U/L	2153	1721	—	2585	144
		µkat/L	17.1	13.6	—	20.6	1.2		µkat/L	36.0	28.7	—	43.2	2.4
	BS-380 ⁷	U/L	1064	851	—	1277	71	BS-380 ⁷	U/L	2218	1771	—	2665	149
		µkat/L	17.8	14.2	—	21.3	1.2		µkat/L	37.0	29.6	—	44.5	2.5
	BS-400 ⁸	U/L	1038	828	—	1248	70	BS-400 ⁸	U/L	2198	1757	—	2639	147
		µkat/L	17.3	13.8	—	20.8	1.2		µkat/L	36.7	29.3	—	44.1	2.5
	BS-480 ⁹	U/L	1095	876	—	1314	73	BS-480 ⁹	U/L	2265	1809	—	2721	152
		µkat/L	18.3	14.6	—	21.9	1.2		µkat/L	37.8	30.2	—	45.4	2.5

Abbreviated name	G6PD Control (L)							G6PD Control (H)						
	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)			1 SD	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)			1 SD
G6PD	BS-600 ¹⁰	U/L	1095	876	—	1314	73	BS-600 ¹⁰	U/L	2265	1809	—	2721	152
		μ kat/L	18.3	14.6	—	21.9	1.2		μ kat/L	37.8	30.2	—	45.4	2.5
	BS-600M ¹¹	U/L	1101	879	—	1323	74	BS-600M ¹¹	U/L	2300	1838	—	2762	154
		μ kat/L	18.4	14.7	—	22.1	1.2		μ kat/L	38.4	30.7	—	46.1	2.6
	BS-800 ¹²	U/L	1095	876	—	1314	73	BS-800 ¹²	U/L	2265	1809	—	2721	152
		μ kat/L	18.3	14.6	—	21.9	1.2		μ kat/L	37.8	30.2	—	45.4	2.5
	BS-2000 ¹³	U/L	1132	904	—	1360	76	BS-2000 ¹³	U/L	2295	1833	—	2757	154
		μ kat/L	18.9	15.1	—	22.7	1.3		μ kat/L	38.3	30.6	—	46.0	2.6
	BS-2800M ¹⁴	U/L	1134	906	—	1362	76	BS-2800M ¹⁴	U/L	2320	1855	—	2785	155
		μ kat/L	18.9	15.1	—	22.7	1.3		μ kat/L	38.7	31.0	—	46.5	2.6

English	Abbreviated name	G6PD Control	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)
Русский	сокращенное наименование	Контроль G6PD	модель	Прибор	Результат анализа	Диапазон(результат анализа \pm 3CO)
Português	Nome abreviado	Controle G6PD	Modelo	Unidade	Valores da análise	Faixa (Valores da análise \pm 3SD)
Español	nombre abreviado	Control G6PD	modelo	Unidad	Valor de ensayo	Rango (Valor de ensayo \pm 3SD)
Italiano	abbreviazione	Controllo G6PD	modelli	Unità	Valori di dosaggio	Intervallo(valore di concentrazione \pm 3 DS)
Türkçe	kısaltılmış ad	G6PD Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri \pm 3SD)

G6PD	
English	Glucose-6-Phosphate Dehydrogenase
Русский	глюкозо-6-фосфатдегидрогеназа
Português	Glicose-6-fosfato desidrogenase
Español	Glucosa-6-fosfato deshidrogenasa
Italiano	Glucosio-6-fosfato deidrogenasi
Türkçe	Glikoz 6 Fosfat Dehidrogenaz