

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.BS-120: BS-120, BS-130, BS-180, BS-190;

8.BS-400: BS-400, BS-420;

2.BS-200: BS-200, BS-220, BS-330, BS-350;

9.BS-480: BS-480, BS-490;

3.BS-200E: BS-200E, BS-220E,;

10.BS-600: BS-600, BS-620;

4.BS-240: BS-230, BS-240;

11.BS-600M: BS-600M, BS-620M;

5.BS-240E: BS240E, BS240Pro;

12.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

6.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-

13.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

330E(V35.00)(Serial Number starts with "W8-" and software version

starts with "35.00"), BS-350E(V35.00) (Serial Number starts with "W9-"

14.BS-2800M:BS-2600M, BS-2800M;

and software version starts with "35.00");

15.Applicable models of the chemistry shall be subject to the parameter list and instructions.

7.BS-380: BS-380, BS-390;

LOT : 057523005

: 2025-01-11

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

G6PD

English	Glucose-6-Phosphate Dehydrogenase
Русский	глюкозо-6-фосфатдегидрогеназа
Português	Glicose-6-fosfato desidrogenase
Español	Glucosa-6-fosfato deshidrogenasa
Italiano	Glucosio-6-fosfato deidrogenasi
Türkçe	Glikoz 6 Fosfat Dehidrogenaz

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
G6PD	U/L	BS-120 ¹	1144	913	—	1375	77	BS-120 ¹	2382	1902	—	2862	160
		BS-200 ²	1183	946	—	1420	79	BS-200 ²	2455	1963	—	2947	164
		BS-200E ³	1014	810	—	1218	68	BS-200E ³	2123	1697	—	2549	142
		BS-240 ⁴	1061	848	—	1274	71	BS-240 ⁴	2211	1767	—	2655	148
		BS-240E ⁵	1113	888	—	1338	75	BS-240E ⁵	2324	1856	—	2792	156
		BS-360E ⁶	1024	817	—	1231	69	BS-360E ⁶	2153	1721	—	2585	144
		BS-380 ⁷	1064	851	—	1277	71	BS-380 ⁷	2218	1771	—	2665	149
		BS-400 ⁸	1038	828	—	1248	70	BS-400 ⁸	2198	1757	—	2639	147
		BS-480 ⁹	1095	876	—	1314	73	BS-480 ⁹	2265	1809	—	2721	152
		BS-600 ¹⁰	1095	876	—	1314	73	BS-600 ¹⁰	2265	1809	—	2721	152
		BS-600M ¹¹	1101	879	—	1323	74	BS-600M ¹¹	2300	1838	—	2762	154
		BS-800 ¹²	1095	876	—	1314	73	BS-800 ¹²	2265	1809	—	2721	152
		BS-2000 ¹³	1132	904	—	1360	76	BS-2000 ¹³	2295	1833	—	2757	154
		BS-2800M ¹⁴	1134	906	—	1362	76	BS-2800M ¹⁴	2320	1855	—	2785	155
µkat/L		BS-120 ¹	19.1	15.2	—	23.0	1.3	BS-120 ¹	39.8	31.8	—	47.8	2.7
		BS-200 ²	19.8	15.8	—	23.7	1.3	BS-200 ²	41.0	32.8	—	49.2	2.7
		BS-200E ³	16.9	13.5	—	20.3	1.1	BS-200E ³	35.5	28.3	—	42.6	2.4

G6PD	μkat/L	BS-240⁴	17.7	14.2	—	21.3	1.2	BS-240⁴	36.9	29.5	—	44.3	2.5
		BS-240E⁵	18.6	14.8	—	22.3	1.3	BS-240E⁵	38.8	31.0	—	46.6	2.6
		BS-360E⁶	17.1	13.6	—	20.6	1.2	BS-360E⁶	36.0	28.7	—	43.2	2.4
		BS-380⁷	17.8	14.2	—	21.3	1.2	BS-380⁷	37.0	29.6	—	44.5	2.5
		BS-400⁸	17.3	13.8	—	20.8	1.2	BS-400⁸	36.7	29.3	—	44.1	2.5
		BS-480⁹	18.3	14.6	—	21.9	1.2	BS-480⁹	37.8	30.2	—	45.4	2.5
		BS-600¹⁰	18.3	14.6	—	21.9	1.2	BS-600¹⁰	37.8	30.2	—	45.4	2.5
		BS-600M¹¹	18.4	14.7	—	22.1	1.2	BS-600M¹¹	38.4	30.7	—	46.1	2.6
		BS-800¹²	18.3	14.6	—	21.9	1.2	BS-800¹²	37.8	30.2	—	45.4	2.5
		BS-2000¹³	18.9	15.1	—	22.7	1.3	BS-2000¹³	38.3	30.6	—	46.0	2.6
		BS-2800M¹⁴	18.9	15.1	—	22.7	1.3	BS-2800M¹⁴	38.7	31.0	—	46.5	2.6

mindray

G6PD Control

For use on: BS-2000, BS-2800M,
BS-600M, BS-800

LOT 057523005

2025-01-11

