

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов измерительной системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадает во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

1. **BS-120**: BS-120, BS-130;

2. **BS-200**: BS-200, BS-220;

3. **BS-200E**: BS-200E, BS-220E;

4. **BS-230**: BS-230, BS-240;

5. **BS-240E**: BS-240E, BS-240Pro;

6. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), **BS-350E(V35.00)** (Serial Number starts with "W9-" and software version starts with "35.00");

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

7. **BS-380**: BS-380, BS-390;

8. **BS-400**: BS-400, BS-420;

9. **BS-480**: BS-480, BS-490;

10. **BS-600**: BS-600, BS-620;

11. **BS-600M**: BS-600M, BS-620M;

12. **BS-800**: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

13. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

14. **BS-2800M**: BS-2600M, BS-2800M.

LOT : 057523006

: 2025-02-10

Abbreviated name	G6PD Control (L)						G6PD Control (H)							
	Model	Unit	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value±3SD)		1 SD		
G6PD	BS-120 ¹	U/L	1176	939	—	1413	79	BS-120 ¹	U/L	2406	1923	—	2889	161
		µkat/L	19.6	15.7	—	23.6	1.3		µkat/L	40.2	32.1	—	48.2	2.7
	BS-200 ²	U/L	1171	937	—	1405	78	BS-200 ²	U/L	2410	1927	—	2893	161
		µkat/L	19.6	15.6	—	23.5	1.3		µkat/L	40.2	32.2	—	48.3	2.7
	BS-200E ³	U/L	1013	809	—	1217	68	BS-200E ³	U/L	2095	1675	—	2515	140
		µkat/L	16.9	13.5	—	20.3	1.1		µkat/L	35.0	28.0	—	42.0	2.3
	BS-230 ⁴	U/L	1050	840	—	1260	70	BS-230 ⁴	U/L	2168	1733	—	2603	145
		µkat/L	17.5	14.0	—	21.0	1.2		µkat/L	36.2	28.9	—	43.5	2.4
	BS-240E ⁵	U/L	1119	894	—	1344	75	BS-240E ⁵	U/L	2288	1829	—	2747	153
		µkat/L	18.7	14.9	—	22.4	1.3		µkat/L	38.2	30.5	—	45.9	2.6
	BS-360E ⁶	U/L	1042	832	—	1252	70	BS-360E ⁶	U/L	2148	1716	—	2580	144
		µkat/L	17.4	13.9	—	20.9	1.2		µkat/L	35.9	28.7	—	43.1	2.4
	BS-380 ⁷	U/L	1042	832	—	1252	70	BS-380 ⁷	U/L	2139	1710	—	2568	143
		µkat/L	17.4	13.9	—	20.9	1.2		µkat/L	35.7	28.6	—	42.9	2.4
	BS-400 ⁸	U/L	1037	830	—	1244	69	BS-400 ⁸	U/L	2167	1732	—	2602	145
		µkat/L	17.3	13.9	—	20.8	1.2		µkat/L	36.2	28.9	—	43.5	2.4
BS-480 ⁹	U/L	1091	872	—	1310	73	BS-480 ⁹	U/L	2234	1784	—	2684	150	
	µkat/L	18.2	14.6	—	21.9	1.2		µkat/L	37.3	29.8	—	44.8	2.5	

Abbreviated name	G6PD Control (L)						G6PD Control (H)							
	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD		
G6PD	BS-600 ¹⁰	U/L	1091	872	—	1310	73	BS-600 ¹⁰	U/L	2234	1784	—	2684	150
		μ kat/L	18.2	14.6	—	21.9	1.2		μ kat/L	37.3	29.8	—	44.8	2.5
	BS-600M ¹¹	U/L	1105	883	—	1327	74	BS-600M ¹¹	U/L	2261	1808	—	2714	151
		μ kat/L	18.5	14.7	—	22.2	1.2		μ kat/L	37.8	30.2	—	45.3	2.5
	BS-800 ¹²	U/L	1091	872	—	1310	73	BS-800 ¹²	U/L	2234	1784	—	2684	150
		μ kat/L	18.2	14.6	—	21.9	1.2		μ kat/L	37.3	29.8	—	44.8	2.5
	BS-2000 ¹³	U/L	1134	906	—	1362	76	BS-2000 ¹³	U/L	2275	1819	—	2731	152
		μ kat/L	18.9	15.1	—	22.7	1.3		μ kat/L	38.0	30.4	—	45.6	2.5
	BS-2800M ¹⁴	U/L	1153	922	—	1384	77	BS-2800M ¹⁴	U/L	2309	1844	—	2774	155
		μ kat/L	19.3	15.4	—	23.1	1.3		μ kat/L	38.6	30.8	—	46.3	2.6

English	Abbreviated name	G6PD Control	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)
Русский	сокращенное наименование	Контроль G6PD	модель	Прибор	Результат анализа	Диапазон(результат анализа \pm 3CO)
Português	Nome abreviado	Controle G6PD	Modelo	Unidade	Valores da análise	Faixa (Valores da análise \pm 3SD)
Español	nombre abreviado	Control G6PD	modelo	Unidad	Valor de ensayo	Rango (Valor de ensayo \pm 3SD)
Italiano	abbreviazione	Controllo G6PD	modelli	Unità	Valori di dosaggio	Intervallo(valore di concentrazione \pm 3 DS)
Türkçe	kısaltılmış ad	G6PD Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri \pm 3SD)

G6PD	
English	Glucose-6-Phosphate Dehydrogenase
Русский	глюкозо-6-фосфатдегидрогеназа
Português	Glicose-6-fosfato desidrogenase
Español	Glucosa-6-fosfato deshidrogenasa
Italiano	Glucosio-6-fosfato deidrogenasi
Türkçe	Glikoz 6 Fosfat Dehidrogenaz