

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.BS-120: BS-120, BS-130, BS-180, BS-190;

9.BS-480: BS-480, BS-490;

2.BS-200: BS-200, BS-220, BS-330, BS-350;

10.BS-600: BS-600, BS-620;

3.BS-200E: BS-200E, BS-220E,;

11.BS-600M: BS-600M, BS-620M;

4.BS-240: BS-230, BS-240;

12.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

5.BS-240E: BS240E, BS240Pro;

13.BS-1000M: BS-1000M, BS-1100M;

6.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version

14.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

starts with "35.00"), BS-350E(V35.00) (Serial Number starts with "W9-" and software version starts with "35.00");

15.BS-2800M:BS-2600M, BS-2800M;

7.BS-380: BS-380, BS-390;

16.For applicable models of the analyte, please refer to the reagent parameter sheet and instruction.

8.BS-400: BS-400, BS-420;

LOT : 057524003

: 2025-07-29

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

G6PD

English	Glucose-6-Phosphate Dehydrogenase
Русский	глюкозо-6-фосфатдегидрогеназа
Português	Glicose-6-fosfato desidrogenase
Español	Glucosa-6-fosfato deshidrogenasa
Italiano	Glucosio-6-fosfato deidrogenasi
Türkçe	Glikoz 6 Fosfat Dehidrogenaz

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
G6PD	U/L	BS-120 ¹	1129	901	—	1357	76	BS-120 ¹	2265	1809	—	2721	152
		BS-200 ²	1164	930	—	1398	78	BS-200 ²	2309	1844	—	2774	155
		BS-200E ³	1003	802	—	1204	67	BS-200E ³	2081	1664	—	2498	139
		BS-240 ⁴	1052	842	—	1262	70	BS-240 ⁴	2111	1688	—	2534	141
		BS-240E ⁵	1098	876	—	1320	74	BS-240E ⁵	2228	1781	—	2675	149
		BS-360E ⁶	1040	830	—	1250	70	BS-360E ⁶	2081	1664	—	2498	139
		BS-380 ⁷	1068	852	—	1284	72	BS-380 ⁷	2267	1811	—	2723	152
		BS-400 ⁸	1034	827	—	1241	69	BS-400 ⁸	2187	1746	—	2628	147
		BS-480 ⁹	1151	920	—	1382	77	BS-480 ⁹	2267	1811	—	2723	152
		BS-600 ¹⁰	1151	920	—	1382	77	BS-600 ¹⁰	2267	1811	—	2723	152
		BS-600M ¹¹	1100	878	—	1322	74	BS-600M ¹¹	2230	1783	—	2677	149
		BS-800 ¹²	1094	875	—	1313	73	BS-800 ¹²	2187	1746	—	2628	147
		BS-1000M ¹³	1029	822	—	1236	69	BS-1000M ¹³	2051	1640	—	2462	137
		BS-2000 ¹⁴	1118	893	—	1343	75	BS-2000 ¹⁴	2219	1772	—	2666	149
		BS-2800M ¹⁵	1151	920	—	1382	77	BS-2800M ¹⁵	2267	1811	—	2723	152
		BS-120 ¹	18.9	15.0	—	22.7	1.3	BS-120 ¹	37.8	30.2	—	45.4	2.5
		BS-200 ²	19.4	15.5	—	23.3	1.3	BS-200 ²	38.6	30.8	—	46.3	2.6


	BS-200E³	16.8	13.4	—	20.1	1.1	BS-200E³	34.8	27.8	—	41.7	2.3
	BS-240⁴	17.6	14.1	—	21.1	1.2	BS-240⁴	35.3	28.2	—	42.3	2.4
	BS-240E⁵	18.3	14.6	—	22.0	1.2	BS-240E⁵	37.2	29.7	—	44.7	2.5
	BS-360E⁶	17.4	13.9	—	20.9	1.2	BS-360E⁶	34.8	27.8	—	41.7	2.3
	BS-380⁷	17.8	14.2	—	21.4	1.2	BS-380⁷	37.9	30.2	—	45.5	2.5
	BS-400⁸	17.3	13.8	—	20.7	1.2	BS-400⁸	36.5	29.2	—	43.9	2.5
$\mu\text{kat/L}$	BS-480⁹	19.2	15.4	—	23.1	1.3	BS-480⁹	37.9	30.2	—	45.5	2.5
	BS-600¹⁰	19.2	15.4	—	23.1	1.3	BS-600¹⁰	37.9	30.2	—	45.5	2.5
	BS-600M¹¹	18.4	14.7	—	22.1	1.2	BS-600M¹¹	37.2	29.8	—	44.7	2.5
	BS-800¹²	18.3	14.6	—	21.9	1.2	BS-800¹²	36.5	29.2	—	43.9	2.5
	BS-1000M¹³	17.2	13.7	—	20.6	1.2	BS-1000M¹³	34.3	27.4	—	41.1	2.3
	BS-2000¹⁴	18.7	14.9	—	22.4	1.3	BS-2000¹⁴	37.1	29.6	—	44.5	2.5
	BS-2800M¹⁵	19.2	15.4	—	23.1	1.3	BS-2800M¹⁵	37.9	30.2	—	45.5	2.5

mindray

G6PD Control

For use on: BS-1000M, BS-2000,
BS-2800M, BS-600M,
BS-800

LOT 057524003

 2025-07-29

