

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

- 1.BS-120: BS-120, BS-130, BS-180, BS-190;
- 2.BS-200: BS-200, BS-220, BS-330, BS-350;
- 3.BS-200E: BS-200E, BS-220E,;
- 4.BS-240: BS-230, BS-240;
- 5.BS-240E: BS240E, BS240Pro;
- 6.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00");
- 7.BS-380: BS-380, BS-390;
- 8.BS-400: BS-400, BS-420;

- 9.BS-430: BS-430, BS-450,BS-460;
- 10.BS-480: BS-480, BS-490;
- 11.BS-600: BS-600, BS-620;
- 12.BS-600M: BS-600M, BS-620M;
- 13.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;
- 14.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;
- 15.BS-2800M:BS-2600M, BS-2800M;
- 16.Applicable models of the chemistry shall be subject to the parameter list and instructions.

LOT : 057623007

Exp. Date: 2025-02-05

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

β-HB

English	β-Hydroxybutyrate
Русский	β-гидроксибутирата
Português	β-Hidroxibutirato
Español	β-hidroxibutirato
Italiano	β-idrossibutirato
Türkçe	β-Hidroksibötrat

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
β-HB	mmol/L	BS-120 ¹	0.298	0.238	—	0.358	0.020	BS-120 ¹	1.16	0.92	—	1.40	0.08
		BS-200 ²	0.294	0.234	—	0.354	0.020	BS-200 ²	1.20	0.96	—	1.44	0.08
		BS-200E ³	0.295	0.235	—	0.355	0.020	BS-200E ³	1.17	0.93	—	1.41	0.08
		BS-240 ⁴	0.298	0.238	—	0.358	0.020	BS-240 ⁴	1.16	0.92	—	1.40	0.08
		BS-240E ⁵	0.297	0.237	—	0.357	0.020	BS-240E ⁵	1.17	0.93	—	1.41	0.08
		BS-360E ⁶	0.297	0.237	—	0.357	0.020	BS-360E ⁶	1.17	0.93	—	1.41	0.08
		BS-380 ⁷	0.292	0.232	—	0.352	0.020	BS-380 ⁷	1.17	0.93	—	1.41	0.08
		BS-400 ⁸	0.297	0.237	—	0.357	0.020	BS-400 ⁸	1.19	0.95	—	1.43	0.08
		BS-430 ⁹	0.297	0.237	—	0.357	0.020	BS-430 ⁹	1.17	0.93	—	1.41	0.08
		BS-480 ¹⁰	0.297	0.237	—	0.357	0.020	BS-480 ¹⁰	1.17	0.93	—	1.41	0.08
		BS-600 ¹¹	0.297	0.237	—	0.357	0.020	BS-600 ¹¹	1.17	0.93	—	1.41	0.08
		BS-600M ¹²	0.294	0.234	—	0.354	0.020	BS-600M ¹²	1.17	0.93	—	1.41	0.08
		BS-800 ¹³	0.297	0.237	—	0.357	0.020	BS-800 ¹³	1.17	0.93	—	1.41	0.08
		BS-2000 ¹⁴	0.295	0.235	—	0.355	0.020	BS-2000 ¹⁴	1.17	0.93	—	1.41	0.08
		BS-2800M ¹⁵	0.294	0.234	—	0.354	0.020	BS-2800M ¹⁵	1.17	0.93	—	1.41	0.08

mindray
β-HB Control

For use on: BS-2000, BS-2800M,
BS-600M, BS-800

LOT 057623007

EXP 2025-02-05

