

## Value sheet of Mindray BS Measurement System

Русский: Таблица результатов измерительной системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадает во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.**BS-120**: BS-120, BS-130;

7.**BS-380**: BS-380, BS-390;

2.**BS-200**: BS-200, BS-220;

8.**BS-400**: BS-400, BS-420;

3.**BS-200E**: BS-200E, BS-220E;

9.**BS-430**: BS-430, BS-450, BS-460;

4.**BS-230**: BS-230, BS-240;

10.**BS-480**: BS-480, BS-490;

5.**BS-240E**: BS-240E, BS-240Pro;

11.**BS-600**: BS-600, BS-620;

6.**BS-360E**: BS-360E,BS-370E,BS-350S,BS-360S, **BS-330E(V35.00)**

12.**BS-600M**: BS-600M, BS-620M;

(Serial Number starts with "W8-" and software version starts

13.**BS-800**: BS-800, BS-820, BS-800M, BS-820M,

BS-1800, BS-1800plus;

with "35.00"), **BS-350E(V35.00)**(Serial Number starts with

14.**BS-2000**: BS-2000, BS-2200, BS-2000M,

BS-2200M;

"W9-" and software version starts with "35.00");

15.**BS-2800M**: BS-2600M, BS-2800M.

**LOT**: 057722003

**EXP**: 2023-03-05

Abbreviated name	ACE Control (L)					ACE Control (H)							
	Model	Unit	Assay Value	Range (Assay Value $\pm$ 3SD)	1 SD	Model	Unit	Assay Value	Range (Assay Value $\pm$ 3SD)	1 SD			
ACE	<b>BS-120</b> <sup>1</sup>	U/L	/	/	-	/	/	U/L	/	/	-	/	/
		$\mu$ kat/L	/	/	-	/	/	$\mu$ kat/L	/	/	-	/	/
	<b>BS-200</b> <sup>2</sup>	U/L	/	/	-	/	/	U/L	/	/	-	/	/
		$\mu$ kat/L	/	/	-	/	/	$\mu$ kat/L	/	/	-	/	/
	<b>BS-200E</b> <sup>3</sup>	U/L	/	/	-	/	/	U/L	/	/	-	/	/
		$\mu$ kat/L	/	/	-	/	/	$\mu$ kat/L	/	/	-	/	/
	<b>BS-230</b> <sup>4</sup>	U/L	/	/	-	/	/	U/L	/	/	-	/	/
		$\mu$ kat/L	/	/	-	/	/	$\mu$ kat/L	/	/	-	/	/
	<b>BS-240E</b> <sup>5</sup>	U/L	/	/	-	/	/	U/L	/	/	-	/	/
		$\mu$ kat/L	/	/	-	/	/	$\mu$ kat/L	/	/	-	/	/
	<b>BS-360E</b> <sup>6</sup>	U/L	/	/	-	/	/	U/L	/	/	-	/	/
		$\mu$ kat/L	/	/	-	/	/	$\mu$ kat/L	/	/	-	/	/
	<b>BS-380</b> <sup>7</sup>	U/L	/	/	-	/	/	U/L	/	/	-	/	/
		$\mu$ kat/L	/	/	-	/	/	$\mu$ kat/L	/	/	-	/	/
	<b>BS-400</b> <sup>8</sup>	U/L	/	/	-	/	/	U/L	/	/	-	/	/
		$\mu$ kat/L	/	/	-	/	/	$\mu$ kat/L	/	/	-	/	/

Abbreviated name	ACE Control (L)							ACE Control (H)						
	Model	Unit	Assay Value	Range (Assay Value±3SD)			1 SD	Model	Unit	Assay Value	Range (Assay Value±3SD)			1 SD
ACE	BS-430 <sup>9</sup>	U/L	/	/	-	/	/	BS-430 <sup>9</sup>	U/L	/	/	-	/	/
		µkat/L	/	/	-	/	/		µkat/L	/	/	-	/	/
	BS-480 <sup>10</sup>	U/L	/	/	-	/	/	BS-480 <sup>10</sup>	U/L	/	/	-	/	/
		µkat/L	/	/	-	/	/		µkat/L	/	/	-	/	/
	BS-600 <sup>11</sup>	U/L	/	/	-	/	/	BS-600 <sup>11</sup>	U/L	/	/	-	/	/
		µkat/L	/	/	-	/	/		µkat/L	/	/	-	/	/
	BS-600M <sup>12</sup>	U/L	45.0	33.8	-	56.2	3.7	BS-600M <sup>12</sup>	U/L	120	90	-	150	10
		µkat/L	0.752	0.564	-	0.939	0.062		µkat/L	2.00	1.50	-	2.51	0.17
	BS-800 <sup>13</sup>	U/L	/	/	-	/	/	BS-800 <sup>13</sup>	U/L	/	/	-	/	/
		µkat/L	/	/	-	/	/		µkat/L	/	/	-	/	/
	BS-2000 <sup>14</sup>	U/L	/	/	-	/	/	BS-2000 <sup>14</sup>	U/L	/	/	-	/	/
		µkat/L	/	/	-	/	/		µkat/L	/	/	-	/	/
	BS-2800M <sup>15</sup>	U/L	/	/	-	/	/	BS-2800M <sup>15</sup>	U/L	/	/	-	/	/
		µkat/L	/	/	-	/	/		µkat/L	/	/	-	/	/

English	Abbreviated name	ACE Control	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	Контроль ACE	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Controle ACE	Modelo	Unidade	Valores da análise	Faixa (Valores da análise ±3SD)
Español	nombre abreviado	Control ACE	modelo	Unidad	Valor de ensayo	Rango (Valor de ensayo ±3SD)
Italiano	abbreviazione	Controllo ACE	modelli	Unità	Valori di dosaggio	Intervallo(valore di concentrazione ±3 DS)
Türkçe	kısaltılmış ad	ACE Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

ACE	
English	Angiotensin Converting Enzyme
Русский	ангиотензинпревращающего фермента
Português	Enzima de conversão da angiotensina
Español	Enzima convertidora de angiotensina
Italiano	Enzima convertitore dell'angiotensina
Türkçe	Anjiotensin Dönüştürücü Enzim