

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов измерительной системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадает во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.**BS-120**: BS-120, BS-130;

7.**BS-380**: BS-380, BS-390;

2.**BS-200**: BS-200, BS-220;

8.**BS-400**: BS-400, BS-420;

3.**BS-200E**: BS-200E, BS-220E;

9.**BS-430**: BS-430, BS-450, BS-460;

4.**BS-230**: BS-230, BS-240;

10.**BS-480**: BS-480, BS-490;

5.**BS-240E**: BS-240E, BS-240Pro;

11.**BS-600**: BS-600, BS-620;

6.**BS-360E**: BS-360E,BS-370E,BS-350S,BS-360S, **BS-330E(V35.00)**

12.**BS-600M**: BS-600M, BS-620M;

13.**BS-800**: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

(Serial Number starts with "W8-" and software version starts with "35.00"), **BS-350E(V35.00)**(Serial Number starts with

14.**BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

"W9-" and software version starts with "35.00");

15.**BS-2800M**: BS-2600M, BS-2800M.

LOT : 057723001

EXP : 2024-02-29

Abbreviated name	ACE Control (L)						ACE Control (H)							
	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD		
ACE	BS-120 ¹	U/L	53.7	40.3	—	67.1	4.5	BS-120 ¹	U/L	120	90	—	150	10
		μ kat/L	0.897	0.673	—	1.121	0.075		μ kat/L	2.00	1.50	—	2.51	0.17
	BS-200 ²	U/L	54.0	40.6	—	67.4	4.5	BS-200 ²	U/L	122	92	—	152	10
		μ kat/L	0.902	0.678	—	1.126	0.075		μ kat/L	2.04	1.54	—	2.54	0.17
	BS-200E ³	U/L	51.5	38.7	—	64.3	4.3	BS-200E ³	U/L	117	88	—	146	10
		μ kat/L	0.860	0.646	—	1.074	0.072		μ kat/L	1.95	1.47	—	2.44	0.17
	BS-230 ⁴	U/L	53.7	40.3	—	67.1	4.5	BS-230 ⁴	U/L	120	90	—	150	10
		μ kat/L	0.897	0.673	—	1.121	0.075		μ kat/L	2.00	1.50	—	2.51	0.17
	BS-240E ⁵	U/L	44.6	33.5	—	55.7	3.7	BS-240E ⁵	U/L	103	77	—	129	9
		μ kat/L	0.745	0.559	—	0.930	0.062		μ kat/L	1.72	1.29	—	2.15	0.15
	BS-360E ⁶	U/L	44.6	33.5	—	55.7	3.7	BS-360E ⁶	U/L	103	77	—	129	9
		μ kat/L	0.745	0.559	—	0.930	0.062		μ kat/L	1.72	1.29	—	2.15	0.15
BS-380 ⁷	U/L	49.6	37.2	—	62.0	4.1	BS-380 ⁷	U/L	110	83	—	137	9	
	μ kat/L	0.828	0.621	—	1.035	0.068		μ kat/L	1.84	1.39	—	2.29	0.15	
BS-400 ⁸	U/L	48.7	36.6	—	60.8	4.0	BS-400 ⁸	U/L	106	80	—	132	9	
	μ kat/L	0.813	0.611	—	1.015	0.067		μ kat/L	1.77	1.34	—	2.20	0.15	

Abbreviated name	ACE Control (L)						ACE Control (H)							
	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD		
ACE	BS-430 ⁹	U/L	44.6	33.5	—	55.7	3.7	BS-430 ⁹	U/L	103	77	—	129	9
		μ kat/L	0.745	0.559	—	0.930	0.062		μ kat/L	1.72	1.29	—	2.15	0.15
	BS-480 ¹⁰	U/L	44.6	33.5	—	55.7	3.7	BS-480 ¹⁰	U/L	103	77	—	129	9
		μ kat/L	0.745	0.559	—	0.930	0.062		μ kat/L	1.72	1.29	—	2.15	0.15
	BS-600 ¹¹	U/L	44.6	33.5	—	55.7	3.7	BS-600 ¹¹	U/L	103	77	—	129	9
		μ kat/L	0.745	0.559	—	0.930	0.062		μ kat/L	1.72	1.29	—	2.15	0.15
	BS-600M ¹²	U/L	49.3	37.0	—	61.6	4.1	BS-600M ¹²	U/L	111	83	—	139	9
		μ kat/L	0.823	0.618	—	1.029	0.068		μ kat/L	1.85	1.39	—	2.32	0.15
	BS-800 ¹³	U/L	44.6	33.5	—	55.7	3.7	BS-800 ¹³	U/L	103	77	—	129	9
		μ kat/L	0.745	0.559	—	0.930	0.062		μ kat/L	1.72	1.29	—	2.15	0.15
	BS-2000 ¹⁴	U/L	45.4	34.1	—	56.7	3.8	BS-2000 ¹⁴	U/L	110	83	—	137	9
		μ kat/L	0.758	0.569	—	0.947	0.063		μ kat/L	1.84	1.39	—	2.29	0.15
	BS-2800M ¹⁵	U/L	49.3	37.0	—	61.6	4.1	BS-2800M ¹⁵	U/L	113	85	—	141	9
		μ kat/L	0.823	0.618	—	1.029	0.068		μ kat/L	1.89	1.42	—	2.35	0.15

English	Abbreviated name	ACE Control	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)
Русский	сокращенное наименование	Контроль ACE	модель	Прибор	Результат анализа	Диапазон(результат анализа \pm 3CO)
Português	Nome abreviado	Controle ACE	Modelo	Unidade	Valores da análise	Faixa (Valores da análise \pm 3SD)
Español	nombre abreviado	Control ACE	modelo	Unidad	Valor de ensayo	Rango (Valor de ensayo \pm 3SD)
Italiano	abbreviazione	Controllo ACE	modelli	Unità	Valori di dosaggio	Intervallo(valore di concentrazione \pm 3 DS)
Türkçe	kısaltılmış ad	ACE Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri \pm 3SD)

ACE	
English	Angiotensin Converting Enzyme
Русский	ангиотензинпревращающего фермента
Português	Enzima de conversão da angiotensina
Español	Enzima convertidora de angiotensina
Italiano	Enzima convertitore dell'angiotensina
Türkçe	Anjiotensin Dönüştürücü Enzim