

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов измерительной системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадает во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1. **BS-120**: BS-120, BS-130;

7. **BS-380**: BS-380, BS-390;

2. **BS-200**: BS-200, BS-220;

8. **BS-400**: BS-400, BS-420;

3. **BS-200E**: BS-200E, BS-220E;

9. **BS-430**: BS-430, BS-450, BS-460;

4. **BS-230**: BS-230, BS-240;

10. **BS-480**: BS-480, BS-490;

5. **BS-240E**: BS-240E, BS-240Pro;

11. **BS-600**: BS-600, BS-620;

6. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S, **BS-330E(V35.00)**

12. **BS-600M**: BS-600M, BS-620M;

(Serial Number starts with "W8-" and software version starts

with "35.00"), **BS-350E(V35.00)**(Serial Number starts with

"W9-" and software version starts with "35.00");

13. **BS-800**: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

14. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

15. **BS-2800M**: BS-2600M, BS-2800M.

LOT : 057723006

EXP : 2024-09-04

Abbreviated name	ACE Control (L)						ACE Control (H)							
	Model	Unit	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value±3SD)		1 SD		
ACE	BS-120 ¹	U/L	60.4	45.4	—	75.4	5.0	BS-120 ¹	U/L	151	112	—	190	13
		µkat/L	1.01	0.76	—	1.26	0.08		µkat/L	2.52	1.87	—	3.17	0.22
	BS-200 ²	U/L	60.6	45.6	—	75.6	5.0	BS-200 ²	U/L	154	115	—	193	13
		µkat/L	1.01	0.76	—	1.26	0.08		µkat/L	2.57	1.92	—	3.22	0.22
	BS-200E ³	U/L	52.7	39.5	—	65.9	4.4	BS-200E ³	U/L	133	100	—	166	11
		µkat/L	0.880	0.660	—	1.101	0.073		µkat/L	2.22	1.67	—	2.77	0.18
	BS-230 ⁴	U/L	60.4	45.4	—	75.4	5.0	BS-230 ⁴	U/L	151	112	—	190	13
		µkat/L	1.01	0.76	—	1.26	0.08		µkat/L	2.52	1.87	—	3.17	0.22
	BS-240E ⁵	U/L	49.4	37.1	—	61.7	4.1	BS-240E ⁵	U/L	134	101	—	167	11
		µkat/L	0.825	0.620	—	1.030	0.068		µkat/L	2.24	1.69	—	2.79	0.18
	BS-360E ⁶	U/L	49.4	37.1	—	61.7	4.1	BS-360E ⁶	U/L	134	101	—	167	11
		µkat/L	0.825	0.620	—	1.030	0.068		µkat/L	2.24	1.69	—	2.79	0.18
	BS-380 ⁷	U/L	52.1	39.2	—	65.0	4.3	BS-380 ⁷	U/L	129	96	—	162	11
		µkat/L	0.870	0.655	—	1.086	0.072		µkat/L	2.15	1.60	—	2.71	0.18
	BS-400 ⁸	U/L	51.0	38.4	—	63.6	4.2	BS-400 ⁸	U/L	129	96	—	162	11
		µkat/L	0.852	0.641	—	1.062	0.070		µkat/L	2.15	1.60	—	2.71	0.18

Abbreviated name	ACE Control (L)						ACE Control (H)							
	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD		
ACE	BS-430 ⁹	U/L	49.4	37.1	—	61.7	4.1	BS-430 ⁹	U/L	134	101	—	167	11
		μ kat/L	0.825	0.620	—	1.030	0.068		μ kat/L	2.24	1.69	—	2.79	0.18
	BS-480 ¹⁰	U/L	49.4	37.1	—	61.7	4.1	BS-480 ¹⁰	U/L	134	101	—	167	11
		μ kat/L	0.825	0.620	—	1.030	0.068		μ kat/L	2.24	1.69	—	2.79	0.18
	BS-600 ¹¹	U/L	49.4	37.1	—	61.7	4.1	BS-600 ¹¹	U/L	134	101	—	167	11
		μ kat/L	0.825	0.620	—	1.030	0.068		μ kat/L	2.24	1.69	—	2.79	0.18
	BS-600M ¹²	U/L	52.1	39.2	—	65.0	4.3	BS-600M ¹²	U/L	134	101	—	167	11
		μ kat/L	0.870	0.655	—	1.086	0.072		μ kat/L	2.24	1.69	—	2.79	0.18
	BS-800 ¹³	U/L	49.4	37.1	—	61.7	4.1	BS-800 ¹³	U/L	134	101	—	167	11
		μ kat/L	0.825	0.620	—	1.030	0.068		μ kat/L	2.24	1.69	—	2.79	0.18
	BS-2000 ¹⁴	U/L	51.1	38.5	—	63.7	4.2	BS-2000 ¹⁴	U/L	132	99	—	165	11
		μ kat/L	0.853	0.643	—	1.064	0.070		μ kat/L	2.20	1.65	—	2.76	0.18
	BS-2800M ¹⁵	U/L	54.0	40.5	—	67.5	4.5	BS-2800M ¹⁵	U/L	137	104	—	170	11
		μ kat/L	0.902	0.676	—	1.127	0.075		μ kat/L	2.29	1.74	—	2.84	0.18

English	Abbreviated name	ACE Control	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)
Русский	сокращенное наименование	Контроль ACE	модель	Прибор	Результат анализа	Диапазон(результат анализа \pm 3CO)
Português	Nome abreviado	Controle ACE	Modelo	Unidade	Valores da análise	Faixa (Valores da análise \pm 3SD)
Español	nombre abreviado	Control ACE	modelo	Unidad	Valor de ensayo	Rango (Valor de ensayo \pm 3SD)
Italiano	abbreviazione	Controllo ACE	modelli	Unità	Valori di dosaggio	Intervallo(valore di concentrazione \pm 3 DS)
Türkçe	kısaltılmış ad	ACE Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri \pm 3SD)

ACE	
English	Angiotensin Converting Enzyme
Русский	ангиотензинпревращающего фермента
Português	Enzima de conversão da angiotensina
Español	Enzima convertidora de angiotensina
Italiano	Enzima convertitore dell'angiotensina
Türkçe	Anjiotensin Dönüştürücü Enzim