

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

- 1.BS-120: BS-120, BS-130, BS-180, BS-190;
- 2.BS-200: BS-200, BS-220, BS-330, BS-350;
- 3.BS-200E: BS-200E, BS-220E,;
- 4.BS-240: BS-230, BS-240;
- 5.BS-240E: BS240E, BS240Pro;
- 6.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00");
- 7.BS-380: BS-380, BS-390;
- 8.BS-400: BS-400, BS-420;

- 9.BS-430: BS-430, BS-450,BS-460;
- 10.BS-480: BS-480, BS-490;
- 11.BS-600: BS-600, BS-620;
- 12.BS-600M: BS-600M, BS-620M;
- 13.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;
- 14.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;
- 15.BS-2800M:BS-2600M, BS-2800M;
- 16.Applicable models of the chemistry shall be subject to the parameter list and instructions.

LOT : 057723010

🕒: 2024-11-13

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

ACE

English	Angiotensin Converting Enzyme
Русский	Ангиотензин-превращающего фермента
Português	Enzima de conversão da angiotensina
Español	Enzima convertidora de angiotensina
Italiano	Enzima convertitore dell'angiotensina
Türkçe	Anjiotensin Dönüştürücü Enzim

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
ACE	U/L	BS-120 ¹	55.2	41.4	—	69.0	4.6	BS-120 ¹	158	119	—	197	13
		BS-200 ²	52.6	39.4	—	65.8	4.4	BS-200 ²	158	119	—	197	13
		BS-200E ³	50.1	37.5	—	62.7	4.2	BS-200E ³	145	109	—	181	12
		BS-240 ⁴	55.2	41.4	—	69.0	4.6	BS-240 ⁴	158	119	—	197	13
		BS-240E ⁵	44.5	33.4	—	55.6	3.7	BS-240E ⁵	134	101	—	167	11
		BS-360E ⁶	44.5	33.4	—	55.6	3.7	BS-360E ⁶	134	101	—	167	11
		BS-380 ⁷	45.5	34.1	—	56.9	3.8	BS-380 ⁷	134	101	—	167	11
		BS-400 ⁸	46.3	34.9	—	57.7	3.8	BS-400 ⁸	133	100	—	166	11
		BS-430 ⁹	44.5	33.4	—	55.6	3.7	BS-430 ⁹	134	101	—	167	11
		BS-480 ¹⁰	44.5	33.4	—	55.6	3.7	BS-480 ¹⁰	134	101	—	167	11
		BS-600 ¹¹	44.5	33.4	—	55.6	3.7	BS-600 ¹¹	134	101	—	167	11
		BS-600M ¹²	46.7	35.0	—	58.4	3.9	BS-600M ¹²	141	105	—	177	12
		BS-800 ¹³	44.5	33.4	—	55.6	3.7	BS-800 ¹³	134	101	—	167	11
		BS-2000 ¹⁴	45.7	34.3	—	57.1	3.8	BS-2000 ¹⁴	136	103	—	169	11
		BS-2800M ¹⁵	48.2	36.2	—	60.2	4.0	BS-2800M ¹⁵	142	106	—	178	12
		BS-120 ¹	0.922	0.691	—	1.152	0.077	BS-120 ¹	2.64	1.99	—	3.29	0.22
		BS-200 ²	0.878	0.658	—	1.099	0.073	BS-200 ²	2.64	1.99	—	3.29	0.22

	BS-200E³	0.837	0.626	—	1.047	0.070	BS-200E³	2.42	1.82	—	3.02	0.20
	BS-240⁴	0.922	0.691	—	1.152	0.077	BS-240⁴	2.64	1.99	—	3.29	0.22
	BS-240E⁵	0.743	0.558	—	0.929	0.062	BS-240E⁵	2.24	1.69	—	2.79	0.18
	BS-360E⁶	0.743	0.558	—	0.929	0.062	BS-360E⁶	2.24	1.69	—	2.79	0.18
	BS-380⁷	0.760	0.569	—	0.950	0.063	BS-380⁷	2.24	1.69	—	2.79	0.18
	BS-400⁸	0.773	0.583	—	0.964	0.063	BS-400⁸	2.22	1.67	—	2.77	0.18
$\mu\text{kat/L}$	BS-430⁹	0.743	0.558	—	0.929	0.062	BS-430⁹	2.24	1.69	—	2.79	0.18
	BS-480¹⁰	0.743	0.558	—	0.929	0.062	BS-480¹⁰	2.24	1.69	—	2.79	0.18
	BS-600¹¹	0.743	0.558	—	0.929	0.062	BS-600¹¹	2.24	1.69	—	2.79	0.18
	BS-600M¹²	0.780	0.585	—	0.975	0.065	BS-600M¹²	2.35	1.75	—	2.96	0.20
	BS-800¹³	0.743	0.558	—	0.929	0.062	BS-800¹³	2.24	1.69	—	2.79	0.18
	BS-2000¹⁴	0.763	0.573	—	0.954	0.063	BS-2000¹⁴	2.27	1.72	—	2.82	0.18
	BS-2800M¹⁵	0.805	0.605	—	1.005	0.067	BS-2800M¹⁵	2.37	1.77	—	2.97	0.20

mindray

ACE Control

For use on: BS-2000, BS-2800M,
BS-600M, BS-800

LOT 057723010

2024-11-13

