

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.BS-240: BS-230, BS-240;

8.BS-600: BS-600, BS-620;

2.BS-240E: BS240E, BS240Pro;

9.BS-600M: BS-600M, BS-620M;

3.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-

10.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

330E(V35.00)(Serial Number starts with "W8-" and software version

starts with "35.00"), BS-350E(V35.00) (Serial Number starts with "W9-"

and software version starts with "35.00");

11.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

4.BS-380: BS-380, BS-390;

12.BS-2800M:BS-2600M, BS-2800M;

5.BS-400: BS-400, BS-420;

13.Applicable models of the chemistry shall be subject to the

6.BS-430: BS-430, BS-450,BS-460;

parameter list and instructions.

7.BS-480: BS-480, BS-490;

LOT : 057823010

: 2025-03-15

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа \pm 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise \pm 3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo \pm 3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione \pm 3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri \pm 3SD)

D-Dimer

English	D-Dimer
Русский	D-димера
Português	D-Dimer
Español	Dímero D
Italiano	D-Dimero
Türkçe	D-Dimer

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value \pm 3SD)		1 SD	Model	Assay Value	Range (Assay Value \pm 3SD)		1 SD		
D-Dimer	μ g/mL	BS-240 ¹	1.48	1.03	—	1.93	0.15	BS-240 ¹	15.3	12.9	—	17.7	0.8
		BS-240E ²	1.52	1.07	—	1.97	0.15	BS-240E ²	15.5	13.1	—	17.9	0.8
		BS-360E ³	1.52	1.07	—	1.97	0.15	BS-360E ³	15.5	13.1	—	17.9	0.8
		BS-380 ⁴	1.50	1.05	—	1.95	0.15	BS-380 ⁴	15.2	12.8	—	17.6	0.8
		BS-400 ⁵	1.50	1.05	—	1.95	0.15	BS-400 ⁵	15.2	12.8	—	17.6	0.8
		BS-430 ⁶	1.44	1.02	—	1.86	0.14	BS-430 ⁶	15.1	12.7	—	17.5	0.8
		BS-480 ⁷	1.52	1.07	—	1.97	0.15	BS-480 ⁷	15.5	13.1	—	17.9	0.8
		BS-600 ⁸	1.52	1.07	—	1.97	0.15	BS-600 ⁸	15.5	13.1	—	17.9	0.8
		BS-600M ⁹	1.40	0.98	—	1.82	0.14	BS-600M ⁹	15.2	12.8	—	17.6	0.8
		BS-800 ¹⁰	1.52	1.07	—	1.97	0.15	BS-800 ¹⁰	15.5	13.1	—	17.9	0.8
		BS-2000 ¹¹	1.43	1.01	—	1.85	0.14	BS-2000 ¹¹	15.1	12.7	—	17.5	0.8
		BS-2800M ¹²	1.40	0.98	—	1.82	0.14	BS-2800M ¹²	15.2	12.8	—	17.6	0.8

mindray

D-Dimer Control

For use on: BS-2000, BS-2800M,
BS-600M, BS-800

LOT 057823010

2025-03-15

