

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.BS-240: BS-230, BS-240;

8.BS-600: BS-600, BS-620;

2.BS-240E: BS240E, BS240Pro;

9.BS-600M: BS-600M, BS-620M;

3.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-

10.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

330E(V35.00)(Serial Number starts with "W8-" and software version

starts with "35.00"), BS-350E(V35.00) (Serial Number starts with "W9-"

and software version starts with "35.00");

11.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

4.BS-380: BS-380, BS-390;

12.BS-2800M:BS-2600M, BS-2800M;

5.BS-400: BS-400, BS-420;

13.For applicable analyzers of the analyte, please subject to the

6.BS-430: BS-430, BS-450,BS-460;

parameter sheet and instrument.

7.BS-480: BS-480, BS-490;

LOT : 057824003

: 2025-06-04

| English | Abbreviated name | Model | Unit | Assay Value | Range (Assay Value \pm 3SD) |
|-----------|--------------------------|---------|---------|--------------------|---|
| Русский | сокращенное наименование | модель | Прибор | Результат анализа | Диапазон(результат анализа \pm 3CO) |
| Português | Nome abreviado | Modelo | Unidade | Valores da análise | Faixa(Valores da análise \pm 3SD) |
| Español | nombre abreviado | modelo | Unidad | Valor de ensayo | Rango(Valor de ensayo \pm 3SD) |
| Italiano | abbreviazione | modelli | Unità | Valori di dosaggio | Intervallo (valore di concentrazione \pm 3 DS) |
| Türkçe | kısaltılmış ad | model | Ünite | Tayin Değeri | Aralık (Tayin Değeri \pm 3SD) |

D-Dimer

| | |
|-----------|----------|
| English | D-Dimer |
| Русский | D-димера |
| Português | D-Dimer |
| Español | Dímero D |
| Italiano | D-Dimero |
| Türkçe | D-Dimer |

| Abbreviated name | Unit | Level: L/N | | | | | Level: H/P | | | | | | |
|------------------|------------|------------------------|-------------|-------------------------------|---|------|------------|------------------------|-------------------------------|------|------|------|-----|
| | | Model | Assay Value | Range (Assay Value \pm 3SD) | | 1 SD | Model | Assay Value | Range (Assay Value \pm 3SD) | | 1 SD | | |
| D-Dimer | μ g/mL | BS-240 ¹ | 1.48 | 1.03 | — | 1.93 | 0.15 | BS-240 ¹ | 15.3 | 12.9 | — | 17.7 | 0.8 |
| | | BS-240E ² | 1.52 | 1.07 | — | 1.97 | 0.15 | BS-240E ² | 15.5 | 13.1 | — | 17.9 | 0.8 |
| | | BS-360E ³ | 1.52 | 1.07 | — | 1.97 | 0.15 | BS-360E ³ | 15.5 | 13.1 | — | 17.9 | 0.8 |
| | | BS-380 ⁴ | 1.50 | 1.05 | — | 1.95 | 0.15 | BS-380 ⁴ | 15.2 | 12.8 | — | 17.6 | 0.8 |
| | | BS-400 ⁵ | 1.50 | 1.05 | — | 1.95 | 0.15 | BS-400 ⁵ | 15.2 | 12.8 | — | 17.6 | 0.8 |
| | | BS-430 ⁶ | 1.44 | 1.02 | — | 1.86 | 0.14 | BS-430 ⁶ | 15.1 | 12.7 | — | 17.5 | 0.8 |
| | | BS-480 ⁷ | 1.52 | 1.07 | — | 1.97 | 0.15 | BS-480 ⁷ | 15.5 | 13.1 | — | 17.9 | 0.8 |
| | | BS-600 ⁸ | 1.52 | 1.07 | — | 1.97 | 0.15 | BS-600 ⁸ | 15.5 | 13.1 | — | 17.9 | 0.8 |
| | | BS-600M ⁹ | 1.40 | 0.98 | — | 1.82 | 0.14 | BS-600M ⁹ | 15.2 | 12.8 | — | 17.6 | 0.8 |
| | | BS-800 ¹⁰ | 1.52 | 1.07 | — | 1.97 | 0.15 | BS-800 ¹⁰ | 15.5 | 13.1 | — | 17.9 | 0.8 |
| | | BS-2000 ¹¹ | 1.43 | 1.01 | — | 1.85 | 0.14 | BS-2000 ¹¹ | 15.1 | 12.7 | — | 17.5 | 0.8 |
| | | BS-2800M ¹² | 1.40 | 0.98 | — | 1.82 | 0.14 | BS-2800M ¹² | 15.2 | 12.8 | — | 17.6 | 0.8 |

mindray
D-Dimer Control

For use on: BS-2000, BS-2800M,
BS-600M, BS-800

LOT 057824003

2025-06-04

