

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов измерительной системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадает во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

1. **BS-120**: BS-120, BS-130;

2. **BS-200**: BS-200, BS-220;

3. **BS-200E**: BS-200E, BS-220E;

4. **BS-230**: BS-230, BS-240;

5. **BS-240E**: BS-240E, BS-240Pro;

6. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

BS-330E (V35.00) (Serial Number starts with "W8-" and

software version starts with "35.00"), **BS-350E(V35.00)**

(Serial Number starts with "W9-" and software version

starts with "35.00");

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

7. **BS-380**: BS-380, BS-390;

8. **BS-400**: BS-400, BS-420;

9. **BS-430**: BS-430, BS-450, BS-460;

10. **BS-480**: BS-480, BS-490;

11. **BS-600**: BS-600, BS-620;

12. **BS-600M**: BS-600M, BS-620M;

13. **BS-800**: BS-800, BS-820, BS-800M, BS-820M,

BS-1800, BS-1800plus;

14. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

15. **BS-2800M**: BS-2600M, BS-2800M.

LOT : 057922007

: 2023-11-30

Abbreviated name	TRF Control (L)						TRF Control (H)							
	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD		
TRF	BS-120 ¹	g/L	1.01	0.86	—	1.16	0.05	BS-120 ¹	g/L	2.18	1.85	—	2.51	0.11
		μ mol/L	12.7	10.8	—	14.6	0.6		μ mol/L	27.5	23.3	—	31.6	1.4
	BS-200 ²	g/L	1.01	0.86	—	1.16	0.05	BS-200 ²	g/L	2.24	1.90	—	2.58	0.11
		μ mol/L	12.7	10.8	—	14.6	0.6		μ mol/L	28.2	23.9	—	32.5	1.4
	BS-200E ³	g/L	1.00	0.85	—	1.15	0.05	BS-200E ³	g/L	2.21	1.88	—	2.54	0.11
		μ mol/L	12.6	10.7	—	14.5	0.6		μ mol/L	27.8	23.7	—	32.0	1.4
	BS-230 ⁴	g/L	1.02	0.87	—	1.17	0.05	BS-230 ⁴	g/L	2.22	1.89	—	2.55	0.11
		μ mol/L	12.9	11.0	—	14.7	0.6		μ mol/L	28.0	23.8	—	32.1	1.4
	BS-240E ⁵	g/L	1.01	0.86	—	1.16	0.05	BS-240E ⁵	g/L	2.24	1.90	—	2.58	0.11
		μ mol/L	12.7	10.8	—	14.6	0.6		μ mol/L	28.2	23.9	—	32.5	1.4
	BS-360E ⁶	g/L	1.01	0.86	—	1.16	0.05	BS-360E ⁶	g/L	2.24	1.90	—	2.58	0.11
		μ mol/L	12.7	10.8	—	14.6	0.6		μ mol/L	28.2	23.9	—	32.5	1.4
	BS-380 ⁷	g/L	1.00	0.85	—	1.15	0.05	BS-380 ⁷	g/L	2.21	1.88	—	2.54	0.11
		μ mol/L	12.6	10.7	—	14.5	0.6		μ mol/L	27.8	23.7	—	32.0	1.4
	BS-400 ⁸	g/L	1.00	0.85	—	1.15	0.05	BS-400 ⁸	g/L	2.20	1.87	—	2.53	0.11
		μ mol/L	12.6	10.7	—	14.5	0.6		μ mol/L	27.7	23.6	—	31.9	1.4
	BS-430 ⁹	g/L	1.00	0.85	—	1.15	0.05	BS-430 ⁹	g/L	2.23	1.90	—	2.56	0.11
		μ mol/L	12.6	10.7	—	14.5	0.6		μ mol/L	28.1	23.9	—	32.3	1.4
	BS-480 ¹⁰	g/L	1.01	0.86	—	1.16	0.05	BS-480 ¹⁰	g/L	2.24	1.90	—	2.58	0.11
		μ mol/L	12.7	10.8	—	14.6	0.6		μ mol/L	28.2	23.9	—	32.5	1.4
BS-600 ¹¹	g/L	1.01	0.86	—	1.16	0.05	BS-600 ¹¹	g/L	2.24	1.90	—	2.58	0.11	
	μ mol/L	12.7	10.8	—	14.6	0.6		μ mol/L	28.2	23.9	—	32.5	1.4	

Abbreviated name	TRF Control (L)							TRF Control (H)						
	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)			1 SD	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)			1 SD
TRF	BS-600M ¹²	g/L	1.06	0.90	—	1.22	0.05	BS-600M ¹²	g/L	2.25	1.91	—	2.59	0.11
		μ mol/L	13.4	11.3	—	15.4	0.6		μ mol/L	28.4	24.1	—	32.6	1.4
	BS-800 ¹³	g/L	1.01	0.86	—	1.16	0.05	BS-800 ¹³	g/L	2.24	1.90	—	2.58	0.11
		μ mol/L	12.7	10.8	—	14.6	0.6		μ mol/L	28.2	23.9	—	32.5	1.4
	BS-2000 ¹⁴	g/L	1.01	0.86	—	1.16	0.05	BS-2000 ¹⁴	g/L	2.23	1.90	—	2.56	0.11
		μ mol/L	12.7	10.8	—	14.6	0.6		μ mol/L	28.1	23.9	—	32.3	1.4
	BS-2800M ¹⁵	g/L	1.06	0.90	—	1.22	0.05	BS-2800M ¹⁵	g/L	2.25	1.91	—	2.59	0.11
		μ mol/L	13.4	11.3	—	15.4	0.6		μ mol/L	28.4	24.1	—	32.6	1.4

English	Abbreviated name	TRF Control	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)
Русский	сокращенное наименование	Контроль трансферрина	модель	Прибор	Результат анализа	Диапазон (результат анализа \pm 3CO)
Português	Nome abreviado	Controle TRF	Modelo	Unidade	Valores da análise	Faixa (Valores da análise \pm 3SD)
Español	nombre abreviado	Control TRF	modelo	Unidad	Valor de ensayo	Rango (Valor de ensayo \pm 3SD)
Italiano	abbreviazione	Controllo TRF	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione \pm 3 DS)
Türkçe	kısaltılmış ad	TRF Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri \pm 3SD)

TRF	
English	Transferrin
Русский	трансферрина
Português	Transferrina
Español	Transferrina
Italiano	Transferrina
Türkçe	Transferin