

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.BS-120: BS-120, BS-130, BS-180, BS-190;

9.BS-430: BS-430, BS-450,BS-460;

2.BS-200: BS-200, BS-220, BS-330, BS-350;

10.BS-480: BS-480, BS-490;

3.BS-200E: BS-200E, BS-220E,;

11.BS-600: BS-600, BS-620;

4.BS-240: BS-230, BS-240;

12.BS-600M: BS-600M, BS-620M;

5.BS-240E: BS240E, BS240Pro;

13.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

6.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-

14.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00");

15.BS-2800M:BS-2600M, BS-2800M;

7.BS-380: BS-380, BS-390;

16.For applicable analyzers of the analyte, please subject to the parameter sheet and instrument.

8.BS-400: BS-400, BS-420;

LOT : 057924002

: 2026-03-14

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

TRF

English	Transferrin
Русский	трансферрина
Português	Transferrina
Español	Transferrina
Italiano	Transferrina
Türkçe	Transferin

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
TRF	g/L	BS-120 ¹	1.11	0.93	—	1.29	0.06	BS-120 ¹	2.43	2.07	—	2.79	0.12
		BS-200 ²	1.11	0.93	—	1.29	0.06	BS-200 ²	2.43	2.07	—	2.79	0.12
		BS-200E ³	1.14	0.96	—	1.32	0.06	BS-200E ³	2.42	2.06	—	2.78	0.12
		BS-240 ⁴	1.16	0.98	—	1.34	0.06	BS-240 ⁴	2.47	2.11	—	2.83	0.12
		BS-240E ⁵	1.13	0.95	—	1.31	0.06	BS-240E ⁵	2.42	2.06	—	2.78	0.12
		BS-360E ⁶	1.13	0.95	—	1.31	0.06	BS-360E ⁶	2.42	2.06	—	2.78	0.12
		BS-380 ⁷	1.14	0.96	—	1.32	0.06	BS-380 ⁷	2.42	2.06	—	2.78	0.12
		BS-400 ⁸	1.14	0.96	—	1.32	0.06	BS-400 ⁸	2.42	2.06	—	2.78	0.12
		BS-430 ⁹	1.13	0.95	—	1.31	0.06	BS-430 ⁹	2.42	2.06	—	2.78	0.12
		BS-480 ¹⁰	1.13	0.95	—	1.31	0.06	BS-480 ¹⁰	2.42	2.06	—	2.78	0.12
		BS-600 ¹¹	1.14	0.96	—	1.32	0.06	BS-600 ¹¹	2.42	2.06	—	2.78	0.12
		BS-600M ¹²	1.15	0.97	—	1.33	0.06	BS-600M ¹²	2.44	2.08	—	2.80	0.12
		BS-800 ¹³	1.13	0.95	—	1.31	0.06	BS-800 ¹³	2.42	2.06	—	2.78	0.12
		BS-2000 ¹⁴	1.19	1.01	—	1.37	0.06	BS-2000 ¹⁴	2.49	2.13	—	2.85	0.12
		BS-2800M ¹⁵	1.15	0.97	—	1.33	0.06	BS-2800M ¹⁵	2.44	2.08	—	2.80	0.12
		BS-120 ¹	14.0	11.7	—	16.3	0.8	BS-120 ¹	30.6	26.1	—	35.2	1.5
		BS-200 ²	14.0	11.7	—	16.3	0.8	BS-200 ²	30.6	26.1	—	35.2	1.5

	BS-200E³	14.4	12.1	—	16.6	0.8	BS-200E³	30.5	26.0	—	35.0	1.5
	BS-240⁴	14.6	12.3	—	16.9	0.8	BS-240⁴	31.1	26.6	—	35.7	1.5
	BS-240E⁵	14.2	12.0	—	16.5	0.8	BS-240E⁵	30.5	26.0	—	35.0	1.5
	BS-360E⁶	14.2	12.0	—	16.5	0.8	BS-360E⁶	30.5	26.0	—	35.0	1.5
	BS-380⁷	14.4	12.1	—	16.6	0.8	BS-380⁷	30.5	26.0	—	35.0	1.5
	BS-400⁸	14.4	12.1	—	16.6	0.8	BS-400⁸	30.5	26.0	—	35.0	1.5
$\mu\text{mol/L}$	BS-430⁹	14.2	12.0	—	16.5	0.8	BS-430⁹	30.5	26.0	—	35.0	1.5
	BS-480¹⁰	14.2	12.0	—	16.5	0.8	BS-480¹⁰	30.5	26.0	—	35.0	1.5
	BS-600¹¹	14.4	12.1	—	16.6	0.8	BS-600¹¹	30.5	26.0	—	35.0	1.5
	BS-600M¹²	14.5	12.2	—	16.8	0.8	BS-600M¹²	30.7	26.2	—	35.3	1.5
	BS-800¹³	14.2	12.0	—	16.5	0.8	BS-800¹³	30.5	26.0	—	35.0	1.5
	BS-2000¹⁴	15.0	12.7	—	17.3	0.8	BS-2000¹⁴	31.4	26.8	—	35.9	1.5
	BS-2800M¹⁵	14.5	12.2	—	16.8	0.8	BS-2800M¹⁵	30.7	26.2	—	35.3	1.5

mindray

TRF Control

For use on: BS-2000, BS-2800M,
BS-600M, BS-800

LOT 057924002

🕒 2026-03-14

