

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.BS-200E: BS-200E, BS-220E,;

9.BS-600: BS-600, BS-620;

2.BS-240: BS-230, BS-240;

10.BS-600M: BS-600M, BS-620M;

3.BS-240E: BS240E, BS240Pro;

11.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

4.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-

330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00");

12.BS-1000M: BS-1000M, BS-1100M;

5.BS-380: BS-380, BS-390;

13.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

6.BS-400: BS-400, BS-420;

14.BS-2800M:BS-2600M, BS-2800M;

7.BS-430: BS-430, BS-450,BS-460;

15.For applicable models of the analyte, please refer to the reagent parameter sheet and instruction.

8.BS-480: BS-480, BS-490;

LOT : 058024003

🕒: 2025-06-21

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

FER			IgE		
English		Ferritin			Immunoglobulin E
Русский		ферритин			Иммуноглобулин E
Português		Ferritina			Imunoglobulina E
Español		Ferritina			Inmunoglobulina E
Italiano		Ferritina			Immunoglobulina E
Türkçe		Ferritin			İmmünglobulin E

MYO		
English		Myoglobin
Русский		Миоглобин
Português		Mioglobina
Español		Mioglobina
Italiano		Mioglobina
Türkçe		Myoglobin

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
		BS-200E ¹	124	106	—	142	6	BS-200E ¹	468	399	—	537	23
		BS-240 ²	124	106	—	142	6	BS-240 ²	475	403	—	547	24
		BS-240E ³	124	106	—	142	6	BS-240E ³	468	399	—	537	23
		BS-360E ⁴	124	106	—	142	6	BS-360E ⁴	468	399	—	537	23
		BS-380 ⁵	123	105	—	141	6	BS-380 ⁵	457	388	—	526	23
	ng/mL	BS-400 ⁶	123	105	—	141	6	BS-400 ⁶	457	388	—	526	23
		BS-430 ⁷	124	106	—	142	6	BS-430 ⁷	468	399	—	537	23
		BS-480 ⁸	123	105	—	141	6	BS-480 ⁸	457	388	—	526	23
		BS-600 ⁹	123	105	—	141	6	BS-600 ⁹	457	388	—	526	23
		BS-600M ¹⁰	123	105	—	141	6	BS-600M ¹⁰	457	388	—	526	23
		BS-800 ¹¹	124	106	—	142	6	BS-800 ¹¹	468	399	—	537	23

FER	pmol/L	BS-1000M¹²	123	105	—	141	6	BS-1000M¹²	457	388	—	526	23
		BS-2000¹³	121	103	—	139	6	BS-2000¹³	451	382	—	520	23
		BS-2800M¹⁴	123	105	—	141	6	BS-2800M¹⁴	457	388	—	526	23
		BS-200E¹	279	238	—	319	13	BS-200E¹	1052	897	—	1207	52
		BS-240²	279	238	—	319	13	BS-240²	1067	906	—	1229	54
		BS-240E³	279	238	—	319	13	BS-240E³	1052	897	—	1207	52
		BS-360E⁴	279	238	—	319	13	BS-360E⁴	1052	897	—	1207	52
		BS-380⁵	276	236	—	317	13	BS-380⁵	1027	872	—	1182	52
		BS-400⁶	276	236	—	317	13	BS-400⁶	1027	872	—	1182	52
		BS-430⁷	279	238	—	319	13	BS-430⁷	1052	897	—	1207	52
		BS-480⁸	276	236	—	317	13	BS-480⁸	1027	872	—	1182	52
		BS-600⁹	276	236	—	317	13	BS-600⁹	1027	872	—	1182	52
		BS-600M¹⁰	276	236	—	317	13	BS-600M¹⁰	1027	872	—	1182	52
		BS-800¹¹	279	238	—	319	13	BS-800¹¹	1052	897	—	1207	52
		BS-1000M¹²	276	236	—	317	13	BS-1000M¹²	1027	872	—	1182	52
BS-2000¹³	272	231	—	312	13	BS-2000¹³	1013	858	—	1168	52		
BS-2800M¹⁴	276	236	—	317	13	BS-2800M¹⁴	1027	872	—	1182	52		

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
IgE	IU/mL	BS-240²	79.2	55.5	—	102.9	7.9	BS-240²	529	370	—	688	53
		BS-240E³	81.8	57.2	—	106.4	8.2	BS-240E³	521	365	—	677	52
		BS-360E⁴	82.6	57.7	—	107.5	8.3	BS-360E⁴	507	354	—	660	51
		BS-380⁵	81.6	57.0	—	106.2	8.2	BS-380⁵	523	367	—	679	52
		BS-400⁶	81.6	57.0	—	106.2	8.2	BS-400⁶	523	367	—	679	52
		BS-480⁸	81.6	57.0	—	106.2	8.2	BS-480⁸	523	367	—	679	52
		BS-600⁹	81.6	57.0	—	106.2	8.2	BS-600⁹	523	367	—	679	52
		BS-600M¹⁰	82.5	57.6	—	107.4	8.3	BS-600M¹⁰	523	367	—	679	52
		BS-800¹¹	81.6	57.0	—	106.2	8.2	BS-800¹¹	523	367	—	679	52
		BS-1000M¹²	82.5	57.6	—	107.4	8.3	BS-1000M¹²	523	367	—	679	52
		BS-2000¹³	78.5	54.8	—	102.2	7.9	BS-2000¹³	517	361	—	673	52
		BS-2800M¹⁴	82.5	57.6	—	107.4	8.3	BS-2800M¹⁴	523	367	—	679	52

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
MYO	ng/mL	BS-240²	105	90	—	120	5	BS-240²	403	343	—	463	20
		BS-240E³	106	91	—	121	5	BS-240E³	409	349	—	469	20
		BS-360E⁴	106	91	—	121	5	BS-360E⁴	409	349	—	469	20
		BS-380⁵	106	91	—	121	5	BS-380⁵	407	347	—	467	20
		BS-400⁶	106	91	—	121	5	BS-400⁶	407	347	—	467	20
		BS-430⁷	106	91	—	121	5	BS-430⁷	399	339	—	459	20
		BS-480⁸	106	91	—	121	5	BS-480⁸	407	347	—	467	20
		BS-600⁹	106	91	—	121	5	BS-600⁹	407	347	—	467	20
		BS-600M¹⁰	106	91	—	121	5	BS-600M¹⁰	407	347	—	467	20
		BS-800¹¹	106	91	—	121	5	BS-800¹¹	409	349	—	469	20
		BS-1000M¹²	106	91	—	121	5	BS-1000M¹²	407	347	—	467	20
		BS-2000¹³	108	93	—	123	5	BS-2000¹³	420	357	—	483	21
		BS-2800M¹⁴	106	91	—	121	5	BS-2800M¹⁴	407	347	—	467	20

mindray

Multimmun Control

For use on: BS-1000M, BS-2000,
BS-2800M, BS-600M,
BS-800

LOT 058024003

🕒 2025-06-21

