

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português : Planilha de valores do Sistema de Medição BS da Mindray

Español : Hoja de valores del sistema de medición Mindray BS

Italiano : Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe : Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português : A dados de cada grupo é a mesma.

Español : la datos de cada grupo es la misma.

Italiano : la dati di ogni gruppo è la stessa.

Türkçe : her grubun veri aynıdır.

1.**BS-200E**: BS-200E, BS-220E;

5.**BS-380**: BS-380, BS-390;

2.**BS-230**: BS-230, BS-240;

6.**BS-400**: BS-400, BS-420;

3.**BS-240E**: BS240E, BS240Pro;

7.**BS-430**: BS-430, BS-450, BS-460;

8.**BS-480**: BS-480, BS-490;

4.**BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

9.**BS-600**: BS-600, BS-620;

BS-330E(V35.00)(Serial Number starts with "W8-"

10.**BS-600M**: BS-600M, BS-620M;

and software version starts with "35.00"),

11.**BS-800**: BS-800, BS-820, BS-800M, BS-820M,

BS-1800, BS-1800plus;

BS-350E(V35.00) (Serial Number starts with"W9-"

12.**BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M

and software version starts with "35.00");

13.**BS-2800M**:BS-2600M, BS-2800M.



059022001



: 2023-07-25

Abbreviated name	Cystatin C Control (L)					Cystatin C Control (H)				
	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)	1 SD	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)	
CysC II	BS-200E ¹	mg/L	/	/ - /	/	BS-200E ¹	mg/L	/	/ - /	
	BS-230 ²	mg/L	/	/ - /	/	BS-230 ²	mg/L	/	/ - /	
	BS-240E ³	mg/L	/	/ - /	/	BS-240E ³	mg/L	/	/ - /	
	BS-360E ⁴	mg/L	/	/ - /	/	BS-360E ⁴	mg/L	/	/ - /	
	BS-380 ⁵	mg/L	/	/ - /	/	BS-380 ⁵	mg/L	/	/ - /	
	BS-400 ⁶	mg/L	/	/ - /	/	BS-400 ⁶	mg/L	/	/ - /	
	BS-430 ⁷	mg/L	/	/ - /	/	BS-430 ⁷	mg/L	/	/ - /	
	BS-480 ⁸	mg/L	/	/ - /	/	BS-480 ⁸	mg/L	/	/ - /	
	BS-600 ⁹	mg/L	/	/ - /	/	BS-600 ⁹	mg/L	/	/ - /	
	BS-600M ¹⁰	mg/L	0.697	0.592 - 0.802	0.035	BS-600M ¹⁰	mg/L	1.17	0.99 - 1.35	
	BS-800 ¹¹	mg/L	/	/ - /	/	BS-800 ¹¹	mg/L	/	/ - /	
	BS-2000 ¹²	mg/L	/	/ - /	/	BS-2000 ¹²	mg/L	/	/ - /	
	BS-2800M ¹³	mg/L	/	/ - /	/	BS-2800M ¹³	mg/L	/	/ - /	

English	Abbreviated name	Cystatin C Control	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)
Русский	сокращенное наименование	Контроль цистатина C	модель	Прибор	Результат анализа	Диапазон(результат анализа \pm 3С
Português	Nome abreviado	Controlo Cistatina C	Modelo	Unidade	Valores da análise	Faixa(Valores análise \pm 3SD
Español	nombre abreviado	Control cistatina C	modelo	Unidad	Valor de ensayo	Rango(Valor ensayo \pm 3SD
Italiano	abbreviazione	Controllo Cistatina C	modelli	Unità	Valori di dosaggio	Intervallo(valor concentrazione \pm 3SD
Türkçe	kısaltılmış ad	Cystatin C Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri \pm 3SD

CysC

English	Cystatin C
Русский	цистатина C
Português	Cistatina C
Español	cistatina C
Italiano	cistatina C
Türkçe	Cystatin C

M;

1 SD

/

/

/

/

/

/

/

/

/

0.06

/

/

/

3SD)

ьтат
O)

da
))

de
))

е di
3 DS)

1
)
