

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.**BS-200E**: BS-200E, BS-220E;

5.**BS-380**: BS-380, BS-390;

2.**BS-230**: BS-230, BS-240;

6.**BS-400**: BS-400, BS-420;

3.**BS-240E**: BS240E, BS240Pro;

7.**BS-430**: BS-430, BS-450, BS-460;

4.**BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

8.**BS-480**: BS-480, BS-490;

BS-330E(V35.00)(Serial Number starts with "W8-"

9.**BS-600**: BS-600, BS-620;

and software version starts with "35.00"),

10.**BS-600M**: BS-600M, BS-620M;

BS-350E(V35.00) (Serial Number starts with "W9-"

11.**BS-800**: BS-800, BS-820, BS-800M, BS-820M,

and software version starts with "35.00");

BS-1800, BS-1800plus;

12.**BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

13.**BS-2800M**:BS-2600M, BS-2800M.

LOT: 059022004

EXP: 2023-10-07

Abbreviated name	Cystatin C Control (L)					Cystatin C Control (H)				
	Model	Unit	Assay Value	Range (Assay Value±3SD)	1 SD	Model	Unit	Assay Value	Range (Assay Value±3SD)	1 SD
CysC II	BS-200E ¹	mg/L	0.744	0.632 — 0.856	0.037	BS-200E ¹	mg/L	1.24	1.05 — 1.43	0.06
	BS-230 ²	mg/L	0.732	0.622 — 0.842	0.037	BS-230 ²	mg/L	1.24	1.05 — 1.43	0.06
	BS-240E ³	mg/L	0.737	0.626 — 0.848	0.037	BS-240E ³	mg/L	1.24	1.05 — 1.43	0.06
	BS-360E ⁴	mg/L	0.741	0.630 — 0.852	0.037	BS-360E ⁴	mg/L	1.24	1.05 — 1.43	0.06
	BS-380 ⁵	mg/L	0.744	0.632 — 0.856	0.037	BS-380 ⁵	mg/L	1.24	1.05 — 1.43	0.06
	BS-400 ⁶	mg/L	0.723	0.615 — 0.831	0.036	BS-400 ⁶	mg/L	1.24	1.05 — 1.43	0.06
	BS-430 ⁷	mg/L	0.741	0.630 — 0.852	0.037	BS-430 ⁷	mg/L	1.24	1.05 — 1.43	0.06
	BS-480 ⁸	mg/L	0.740	0.629 — 0.851	0.037	BS-480 ⁸	mg/L	1.23	1.05 — 1.41	0.06
	BS-600 ⁹	mg/L	0.741	0.630 — 0.852	0.037	BS-600 ⁹	mg/L	1.24	1.05 — 1.43	0.06
	BS-600M ¹⁰	mg/L	0.763	0.649 — 0.877	0.038	BS-600M ¹⁰	mg/L	1.24	1.05 — 1.43	0.06
	BS-800 ¹¹	mg/L	0.741	0.630 — 0.852	0.037	BS-800 ¹¹	mg/L	1.24	1.05 — 1.43	0.06
	BS-2000 ¹²	mg/L	0.775	0.659 — 0.891	0.039	BS-2000 ¹²	mg/L	1.26	1.07 — 1.45	0.06
	BS-2800M ¹³	mg/L	0.744	0.632 — 0.856	0.037	BS-2800M ¹³	mg/L	1.23	1.05 — 1.41	0.06

English	Abbreviated name	Cystatin C Control	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	Контроль цистатина C	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Controlo Cistatina C	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	Control cistatina C	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	Controllo Cistatina C	modelli	Unità	Valori di dosaggio	Intervallo(valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	Cystatin C Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

CysC	
English	Cystatin C
Русский	цистатина C
Português	Cistatina C
Español	cistatina C
Italiano	cistatina C
Türkçe	Cystatin C