

Value sheet of Mindray BS Measurement System

Русский : Таблица результатов для системы BS компании Mindray

Português : Planilha de valores do Sistema de Medição BS da Mindray

Español : Hoja de valores del sistema de medición Mindray BS

Italiano : Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe : Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский : Данные совпадают во всех группах.

Português : A dados de cada grupo é a mesma.

Español : la datos de cada grupo es la misma.

Italiano : la dati di ogni gruppo è la stessa.

Türkçe : her grubun veri aynıdır.

1.**BS-200E**: BS-200E, BS-220E;

5.**BS-380**: BS-380, BS-390;

2.**BS-230**: BS-230, BS-240;

6.**BS-400**: BS-400, BS-420;

3.**BS-240E**: BS240E, BS240Pro;

7.**BS-430**: BS-430, BS-450, BS-460;

8.**BS-480**: BS-480, BS-490;

4.**BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

9.**BS-600**: BS-600, BS-620;

10.**BS-600M**: BS-600M, BS-620M;

11.**BS-800**: BS-800, BS-820, BS-800M, BS-820M,
BS-1800, BS-1800plus;

12.**BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

13.**BS-2800M**:BS-2600M, BS-2800M.

BS-330E(V35.00)(Serial Number starts with "W8-"

and software version starts with "35.00"),

BS-350E(V35.00) (Serial Number starts with"W9-"

and software version starts with "35.00");

LOT : 059022006

EXP : 2024-02-29

Abbreviated name	Cystatin C Control (L)							Cystatin C Control (H)						
	Model	Unit	Assay Value	Range (Assay Value±3SD)			1 SD	Model	Unit	Assay Value	Range (Assay Value±3SD)			1 SD
CysC II	BS-200E ¹	mg/L	0.683	0.581	—	0.785	0.034	BS-200E ¹	mg/L	1.15	0.98	—	1.32	0.06
	BS-230 ²	mg/L	0.677	0.575	—	0.779	0.034	BS-230 ²	mg/L	1.15	0.98	—	1.32	0.06
	BS-240E ³	mg/L	0.670	0.570	—	0.771	0.034	BS-240E ³	mg/L	1.14	0.97	—	1.31	0.06
	BS-360E ⁴	mg/L	0.684	0.581	—	0.787	0.034	BS-360E ⁴	mg/L	1.15	0.98	—	1.32	0.06
	BS-380 ⁵	mg/L	0.683	0.581	—	0.785	0.034	BS-380 ⁵	mg/L	1.15	0.98	—	1.32	0.06
	BS-400 ⁶	mg/L	0.676	0.575	—	0.777	0.034	BS-400 ⁶	mg/L	1.14	0.97	—	1.31	0.06
	BS-430 ⁷	mg/L	0.684	0.581	—	0.787	0.034	BS-430 ⁷	mg/L	1.15	0.98	—	1.32	0.06
	BS-480 ⁸	mg/L	0.675	0.574	—	0.776	0.034	BS-480 ⁸	mg/L	1.13	0.96	—	1.30	0.06
	BS-600 ⁹	mg/L	0.684	0.581	—	0.787	0.034	BS-600 ⁹	mg/L	1.15	0.98	—	1.32	0.06
	BS-600M ¹⁰	mg/L	0.701	0.596	—	0.806	0.035	BS-600M ¹⁰	mg/L	1.16	0.99	—	1.33	0.06
	BS-800 ¹¹	mg/L	0.684	0.581	—	0.787	0.034	BS-800 ¹¹	mg/L	1.15	0.98	—	1.32	0.06
	BS-2000 ¹²	mg/L	0.690	0.587	—	0.794	0.035	BS-2000 ¹²	mg/L	1.17	0.99	—	1.35	0.06
	BS-2800M ¹³	mg/L	0.699	0.594	—	0.804	0.035	BS-2800M ¹³	mg/L	1.15	0.98	—	1.32	0.06

English	Abbreviated name	Cystatin C Control	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	Контроль цистатина C	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Controlo Cistatina C	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	Control cistatina C	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	Controllo Cistatina C	modelli	Unità	Valori di dosaggio	Intervallo(valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	Cystatin C Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

CysC	
English	Cystatin C
Русский	цистатина C
Português	Cistatina C
Español	cistatina C
Italiano	cistatina C
Türkçe	Cystatin C