

## Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.**BS-200E**: BS-200E, BS-220E;

5.**BS-380**: BS-380, BS-390;

2.**BS-230**: BS-230, BS-240;

6.**BS-400**: BS-400, BS-420;

3.**BS-240E**: BS240E, BS240Pro;

7.**BS-430**: BS-430, BS-450, BS-460;

4.**BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

8.**BS-480**: BS-480, BS-490;

**BS-330E(V35.00)**(Serial Number starts with "W8-"

9.**BS-600**: BS-600, BS-620;

and software version starts with "35.00"),

10.**BS-600M**: BS-600M, BS-620M;

**BS-350E(V35.00)** (Serial Number starts with "W9-"

11.**BS-800**: BS-800, BS-820, BS-800M, BS-820M,

and software version starts with "35.00");

BS-1800, BS-1800plus;

12.**BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

13.**BS-2800M**:BS-2600M, BS-2800M.

**LOT**: 059023002

**EXP**: 2024-06-25

Abbreviated name	Cystatin C Control (L)					Cystatin C Control (H)				
	Model	Unit	Assay Value	Range (Assay Value $\pm$ 3SD)	1 SD	Model	Unit	Assay Value	Range (Assay Value $\pm$ 3SD)	1 SD
CysC II	<b>BS-200E</b> <sup>1</sup>	mg/L	0.806	0.685 — 0.927	0.040	<b>BS-200E</b> <sup>1</sup>	mg/L	1.33	1.13 — 1.53	0.07
	<b>BS-230</b> <sup>2</sup>	mg/L	0.815	0.693 — 0.937	0.041	<b>BS-230</b> <sup>2</sup>	mg/L	1.33	1.13 — 1.53	0.07
	<b>BS-240E</b> <sup>3</sup>	mg/L	0.795	0.676 — 0.914	0.040	<b>BS-240E</b> <sup>3</sup>	mg/L	1.31	1.11 — 1.51	0.07
	<b>BS-360E</b> <sup>4</sup>	mg/L	0.803	0.683 — 0.923	0.040	<b>BS-360E</b> <sup>4</sup>	mg/L	1.34	1.14 — 1.54	0.07
	<b>BS-380</b> <sup>5</sup>	mg/L	0.806	0.685 — 0.927	0.040	<b>BS-380</b> <sup>5</sup>	mg/L	1.33	1.13 — 1.53	0.07
	<b>BS-400</b> <sup>6</sup>	mg/L	0.804	0.683 — 0.925	0.040	<b>BS-400</b> <sup>6</sup>	mg/L	1.33	1.13 — 1.53	0.07
	<b>BS-430</b> <sup>7</sup>	mg/L	0.803	0.683 — 0.923	0.040	<b>BS-430</b> <sup>7</sup>	mg/L	1.34	1.14 — 1.54	0.07
	<b>BS-480</b> <sup>8</sup>	mg/L	0.802	0.682 — 0.922	0.040	<b>BS-480</b> <sup>8</sup>	mg/L	1.32	1.12 — 1.52	0.07
	<b>BS-600</b> <sup>9</sup>	mg/L	0.803	0.683 — 0.923	0.040	<b>BS-600</b> <sup>9</sup>	mg/L	1.34	1.14 — 1.54	0.07
	<b>BS-600M</b> <sup>10</sup>	mg/L	0.804	0.683 — 0.925	0.040	<b>BS-600M</b> <sup>10</sup>	mg/L	1.32	1.12 — 1.52	0.07
	<b>BS-800</b> <sup>11</sup>	mg/L	0.803	0.683 — 0.923	0.040	<b>BS-800</b> <sup>11</sup>	mg/L	1.34	1.14 — 1.54	0.07
	<b>BS-2000</b> <sup>12</sup>	mg/L	0.806	0.685 — 0.927	0.040	<b>BS-2000</b> <sup>12</sup>	mg/L	1.33	1.13 — 1.53	0.07
	<b>BS-2800M</b> <sup>13</sup>	mg/L	0.811	0.689 — 0.933	0.041	<b>BS-2800M</b> <sup>13</sup>	mg/L	1.33	1.13 — 1.53	0.07

English	Abbreviated name	Cystatin C Control	Model	Unit	Assay Value	Range (Assay Value $\pm$ 3SD)
Русский	сокращенное наименование	Контроль цистатина C	модель	Прибор	Результат анализа	Диапазон(результат анализа $\pm$ 3CO)
Português	Nome abreviado	Controlo Cistatina C	Modelo	Unidade	Valores da análise	Faixa(Valores da análise $\pm$ 3SD)
Español	nombre abreviado	Control cistatina C	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo $\pm$ 3SD)
Italiano	abbreviazione	Controllo Cistatina C	modelli	Unità	Valori di dosaggio	Intervallo(valore di concentrazione $\pm$ 3 DS)
Türkçe	kısaltılmış ad	Cystatin C Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri $\pm$ 3SD)

CysC	
English	Cystatin C
Русский	цистатина C
Português	Cistatina C
Español	cistatina C
Italiano	cistatina C
Türkçe	Cystatin C