

**Value sheet of Mindray BS Measurement System**

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.BS-200E: BS-200E, BS-220E,;

8.BS-480: BS-480, BS-490;

2.BS-240: BS-230, BS-240;

9.BS-600: BS-600, BS-620;

3.BS-240E: BS240E, BS240Pro;

10.BS-600M: BS-600M, BS-620M;

4.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-

11.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with "W9-"

and software version starts with "35.00");

12.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

5.BS-380: BS-380, BS-390;

13.BS-2800M:BS-2600M, BS-2800M;

6.BS-400: BS-400, BS-420;

14.Applicable models of the chemistry shall be subject to the parameter list and instructions.

7.BS-430: BS-430, BS-450,BS-460;

LOT : 059024001

⌚: 2025-03-12

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

**CysC**

English	Cystatin C
Русский	цистатин C
Português	Cistatina C
Español	cistatina C
Italiano	cistatina C
Türkçe	Cystatin C

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
CysC II	mg/L	BS-200E <sup>1</sup>	0.721	0.613	—	0.829	0.036	BS-200E <sup>1</sup>	1.23	1.05	—	1.41	0.06
		BS-240 <sup>2</sup>	0.691	0.586	—	0.796	0.035	BS-240 <sup>2</sup>	1.24	1.06	—	1.42	0.06
		BS-240E <sup>3</sup>	0.761	0.647	—	0.875	0.038	BS-240E <sup>3</sup>	1.28	1.10	—	1.46	0.06
		BS-360E <sup>4</sup>	0.726	0.618	—	0.834	0.036	BS-360E <sup>4</sup>	1.25	1.07	—	1.43	0.06
		BS-380 <sup>5</sup>	0.721	0.613	—	0.829	0.036	BS-380 <sup>5</sup>	1.23	1.05	—	1.41	0.06
		BS-400 <sup>6</sup>	0.712	0.604	—	0.820	0.036	BS-400 <sup>6</sup>	1.24	1.06	—	1.42	0.06
		BS-430 <sup>7</sup>	0.726	0.618	—	0.834	0.036	BS-430 <sup>7</sup>	1.25	1.07	—	1.43	0.06
		BS-480 <sup>8</sup>	0.724	0.616	—	0.832	0.036	BS-480 <sup>8</sup>	1.24	1.06	—	1.42	0.06
		BS-600 <sup>9</sup>	0.726	0.618	—	0.834	0.036	BS-600 <sup>9</sup>	1.25	1.07	—	1.43	0.06
		BS-600M <sup>10</sup>	0.755	0.641	—	0.869	0.038	BS-600M <sup>10</sup>	1.24	1.06	—	1.42	0.06
		BS-800 <sup>11</sup>	0.726	0.618	—	0.834	0.036	BS-800 <sup>11</sup>	1.25	1.07	—	1.43	0.06
		BS-2000 <sup>12</sup>	0.730	0.619	—	0.841	0.037	BS-2000 <sup>12</sup>	1.26	1.08	—	1.44	0.06
		BS-2800M <sup>13</sup>	0.750	0.636	—	0.864	0.038	BS-2800M <sup>13</sup>	1.24	1.06	—	1.42	0.06

**mindray**  
Cystatin C Control

For use on: BS-2000, BS-2800M,  
BS-600M, BS-800

LOT 059024001

2025-03-12

