

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1. **BS-200E**: BS-200E, BS-220E;

5. **BS-380**: BS-380, BS-390;

2. **BS-230**: BS-230, BS-240;

6. **BS-400**: BS-400, BS-420;

3. **BS-240E**: BS240E, BS240Pro;

7. **BS-430**: BS-430, BS-450, BS-460;

4. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

8. **BS-480**: BS-480, BS-490;

BS-330E(V35.00)(Serial Number starts with "W8-"

9. **BS-600**: BS-600, BS-620;

and software version starts with "35.00"),

10. **BS-800**: BS-800, BS-820, BS-800M, BS-820M,
BS-1800, BS-1800plus;

BS-350E(V35.00) (Serial Number starts with "W9-"

11. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

and software version starts with "35.00");

12. **BS-2800M**: BS-2600M, BS-2800M.

LOT: 059122002

EXP: 2023-08-17

Abbreviated name	β2-Microglobulin Control (L)						β2-Microglobulin Control (H)					
	Model	Unit	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value±3SD)		1 SD
β2-MGII	BS-200E ¹	mg/L	1.02	0.87	— 1.17	0.05	BS-200E ¹	mg/L	3.02	2.57	— 3.47	0.15
		nmol/L	86.4	73.7	— 99.1	4.2		nmol/L	256	218	— 294	13
	BS-230 ²	mg/L	1.04	0.88	— 1.20	0.05	BS-230 ²	mg/L	3.04	2.58	— 3.50	0.15
		nmol/L	88.1	74.5	— 101.6	4.2		nmol/L	257	219	— 296	13
	BS-240E ³	mg/L	1.04	0.88	— 1.20	0.05	BS-240E ³	mg/L	3.04	2.58	— 3.50	0.15
		nmol/L	88.1	74.5	— 101.6	4.2		nmol/L	257	219	— 296	13
	BS-360E ⁴	mg/L	1.03	0.88	— 1.18	0.05	BS-360E ⁴	mg/L	2.97	2.52	— 3.42	0.15
		nmol/L	87.2	74.5	— 99.9	4.2		nmol/L	252	213	— 290	13
	BS-380 ⁵	mg/L	1.02	0.87	— 1.17	0.05	BS-380 ⁵	mg/L	3.02	2.57	— 3.47	0.15
		nmol/L	86.4	73.7	— 99.1	4.2		nmol/L	256	218	— 294	13
	BS-400 ⁶	mg/L	1.02	0.87	— 1.17	0.05	BS-400 ⁶	mg/L	3.02	2.57	— 3.47	0.15
		nmol/L	86.4	73.7	— 99.1	4.2		nmol/L	256	218	— 294	13
	BS-430 ⁷	mg/L	1.05	0.89	— 1.21	0.05	BS-430 ⁷	mg/L	3.02	2.57	— 3.47	0.15
		nmol/L	88.9	75.4	— 102.5	4.2		nmol/L	256	218	— 294	13
	BS-480 ⁸	mg/L	1.03	0.88	— 1.18	0.05	BS-480 ⁸	mg/L	3.02	2.57	— 3.47	0.15
		nmol/L	87.2	74.5	— 99.9	4.2		nmol/L	256	218	— 294	13
	BS-600 ⁹	mg/L	1.03	0.88	— 1.18	0.05	BS-600 ⁹	mg/L	3.01	2.56	— 3.46	0.15
		nmol/L	87.2	74.5	— 99.9	4.2		nmol/L	255	217	— 293	13
	BS-800 ¹⁰	mg/L	1.03	0.88	— 1.18	0.05	BS-800 ¹⁰	mg/L	3.04	2.58	— 3.50	0.15
		nmol/L	87.2	74.5	— 99.9	4.2		nmol/L	257	219	— 296	13
	BS-2000 ¹¹	mg/L	1.04	0.88	— 1.20	0.05	BS-2000 ¹¹	mg/L	3.12	2.65	— 3.59	0.16
		nmol/L	88.1	74.5	— 101.6	4.2		nmol/L	264	224	— 304	14
	BS-2800M ¹²	mg/L	1.06	0.90	— 1.22	0.05	BS-2800M ¹²	mg/L	3.04	2.58	— 3.50	0.15
		nmol/L	89.8	76.2	— 103.3	4.2		nmol/L	257	219	— 296	13

β2-Microglobulin Control

English	Abbreviated name	β2-Microglobulin Control	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	Контроль β2-микрोगлобулина	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Controlo β2-Microglobulina	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	Control β2-microglobulina	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	Controllo β2 Microglobulina	modelli	Unità	Valori di dosaggio	Intervallo(valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	β2-Mikroglobulin Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

β2-MG	
English	β2-Microglobulin
Русский	β2-микрोगлобулина
Português	β2-Microglobulina
Español	β2-microglobulina
Italiano	β2-Microglobulina
Türkçe	β2-Mikroglobulin