

Value sheet of Mindray BS Measurement System

Русский : Таблица результатов для системы BS компании Mindray

Português : Planilha de valores do Sistema de Medição BS da Mindray

Español : Hoja de valores del sistema de medición Mindray BS

Italiano : Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe : Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский : Данные совпадают во всех группах.

Português : A dados de cada grupo é a mesma.

Español : la datos de cada grupo es la misma.

Italiano : la dati di ogni gruppo è la stessa.

Türkçe : her grubun veri aynıdır.

1.**BS-200E**: BS-200E, BS-220E;

5.**BS-380**: BS-380, BS-390;

2.**BS-230**: BS-230, BS-240;

6.**BS-400**: BS-400, BS-420;

3.**BS-240E**: BS240E, BS240Pro;

7.**BS-430**: BS-430, BS-450, BS-460;

4.**BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

8.**BS-480**: BS-480, BS-490;

BS-330E(V35.00)(Serial Number starts with "W8-"

9.**BS-600**: BS-600, BS-620;

and software version starts with "35.00"),

10.**BS-600M**: BS-600M, BS-620M;

BS-350E(V35.00) (Serial Number starts with "W9-"

11.**BS-800**: BS-800, BS-820, BS-800M, BS-820M,

BS-1800, BS-1800plus;

and software version starts with "35.00");

12.**BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

13.**BS-2800M**: BS-2600M, BS-2800M.

LOT : 059122005

EXP : 2024-05-09

Abbreviated name	β2-Microglobulin Control (L)						β2-Microglobulin Control (H)							
	Model	Unit	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value±3SD)		1 SD		
β2-MGII	BS-200E ¹	mg/L	1.09	0.93	—	1.25	0.05	BS-200E ¹	mg/L	2.95	2.51	—	3.39	0.15
		nmol/L	92.3	78.8	—	105.9	4.2		nmol/L	250	213	—	287	13
	BS-230 ²	mg/L	1.09	0.93	—	1.25	0.05	BS-230 ²	mg/L	2.99	2.54	—	3.44	0.15
		nmol/L	92.3	78.8	—	105.9	4.2		nmol/L	253	215	—	291	13
	BS-240E ³	mg/L	1.09	0.93	—	1.25	0.05	BS-240E ³	mg/L	2.93	2.49	—	3.37	0.15
		nmol/L	92.3	78.8	—	105.9	4.2		nmol/L	248	211	—	285	13
	BS-360E ⁴	mg/L	1.10	0.94	—	1.27	0.06	BS-360E ⁴	mg/L	2.99	2.54	—	3.44	0.15
		nmol/L	93.2	79.6	—	107.6	5.1		nmol/L	253	215	—	291	13
	BS-380 ⁵	mg/L	1.09	0.93	—	1.25	0.05	BS-380 ⁵	mg/L	2.95	2.51	—	3.39	0.15
		nmol/L	92.3	78.8	—	105.9	4.2		nmol/L	250	213	—	287	13
	BS-400 ⁶	mg/L	1.09	0.93	—	1.25	0.05	BS-400 ⁶	mg/L	2.95	2.51	—	3.39	0.15
		nmol/L	92.3	78.8	—	105.9	4.2		nmol/L	250	213	—	287	13
	BS-430 ⁷	mg/L	1.08	0.92	—	1.24	0.05	BS-430 ⁷	mg/L	2.98	2.53	—	3.43	0.15
		nmol/L	91.5	77.9	—	105.0	4.2		nmol/L	252	214	—	291	13
	BS-480 ⁸	mg/L	1.09	0.93	—	1.25	0.05	BS-480 ⁸	mg/L	3.01	2.56	—	3.46	0.15
		nmol/L	92.3	78.8	—	105.9	4.2		nmol/L	255	217	—	293	13
	BS-600 ⁹	mg/L	1.09	0.93	—	1.25	0.05	BS-600 ⁹	mg/L	2.98	2.53	—	3.43	0.15
		nmol/L	92.3	78.8	—	105.9	4.2		nmol/L	252	214	—	291	13
	BS-600M ¹⁰	mg/L	1.08	0.92	—	1.24	0.05	BS-600M ¹⁰	mg/L	2.96	2.52	—	3.40	0.15
		nmol/L	91.5	77.9	—	105.0	4.2		nmol/L	251	213	—	288	13
BS-800 ¹¹	mg/L	1.10	0.94	—	1.27	0.06	BS-800 ¹¹	mg/L	2.97	2.52	—	3.42	0.15	
	nmol/L	93.2	79.6	—	107.6	5.1		nmol/L	252	213	—	290	13	
BS-2000 ¹²	mg/L	1.08	0.92	—	1.24	0.05	BS-2000 ¹²	mg/L	3.03	2.58	—	3.48	0.15	
	nmol/L	91.5	77.9	—	105.0	4.2		nmol/L	257	219	—	295	13	
BS-2800M ¹³	mg/L	1.08	0.92	—	1.24	0.05	BS-2800M ¹³	mg/L	2.96	2.52	—	3.40	0.15	
	nmol/L	91.5	77.9	—	105.0	4.2		nmol/L	251	213	—	288	13	

β2-Microglobulin Control

English	Abbreviated name	β2-Microglobulin Control	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	Контроль β2-микροглобулина	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Controlo β2-Microglobulina	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	Control β2-microglobulina	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	Controllo β2 Microglobulina	modelli	Unità	Valori di dosaggio	Intervallo(valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	β2-Mikroglobulin Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

β2-MG

English	β2-Microglobulin
Русский	β2-микροглобулина
Português	β2-Microglobulina
Español	β2-microglobulina
Italiano	β2-Microglobulina
Türkçe	β2-Mikroglobulin