

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.**BS-200E**: BS-200E, BS-220E;

5.**BS-380**: BS-380, BS-390;

2.**BS-230**: BS-230, BS-240;

6.**BS-400**: BS-400, BS-420;

3.**BS-240E**: BS240E, BS240Pro;

7.**BS-430**: BS-430, BS-450, BS-460;

4.**BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

8.**BS-480**: BS-480, BS-490;

9.**BS-600**: BS-600, BS-620;

BS-330E(V35.00)(Serial Number starts with "W8-"

10.**BS-600M**: BS-600M, BS-620M;

and software version starts with "35.00"),

11.**BS-800**: BS-800, BS-820, BS-800M, BS-820M,

BS-1800, BS-1800plus;

BS-350E(V35.00) (Serial Number starts with "W9-"

12.**BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

and software version starts with "35.00");

13.**BS-2800M**: BS-2600M, BS-2800M.

LOT: 059123003

EXP: 2024-11-10

Abbreviated name	β2-Microglobulin Control (L)							β2-Microglobulin Control (H)						
	Model	Unit	Assay Value	Range (Assay Value±3SD)		1 SD		Model	Unit	Assay Value	Range (Assay Value±3SD)		1 SD	
β2-MGII	BS-200E ¹	mg/L	1.01	0.86	—	1.16	0.05	BS-200E ¹	mg/L	2.90	2.47	—	3.34	0.15
		nmol/L	85.5	72.8	—	98.3	4.2		nmol/L	246	209	—	283	13
	BS-230 ²	mg/L	1.03	0.88	—	1.18	0.05	BS-230 ²	mg/L	2.90	2.47	—	3.34	0.15
		nmol/L	87.2	74.5	—	99.9	4.2		nmol/L	246	209	—	283	13
	BS-240E ³	mg/L	1.03	0.88	—	1.18	0.05	BS-240E ³	mg/L	2.89	2.46	—	3.32	0.14
		nmol/L	87.2	74.5	—	99.9	4.2		nmol/L	245	208	—	281	12
	BS-360E ⁴	mg/L	1.04	0.88	—	1.20	0.05	BS-360E ⁴	mg/L	2.86	2.43	—	3.29	0.14
		nmol/L	88.1	74.5	—	101.6	4.2		nmol/L	242	206	—	279	12
	BS-380 ⁵	mg/L	1.01	0.86	—	1.16	0.05	BS-380 ⁵	mg/L	2.90	2.47	—	3.34	0.15
		nmol/L	85.5	72.8	—	98.3	4.2		nmol/L	246	209	—	283	13
	BS-400 ⁶	mg/L	1.01	0.86	—	1.16	0.05	BS-400 ⁶	mg/L	2.90	2.47	—	3.34	0.15
		nmol/L	85.5	72.8	—	98.3	4.2		nmol/L	246	209	—	283	13
	BS-430 ⁷	mg/L	1.03	0.88	—	1.18	0.05	BS-430 ⁷	mg/L	2.88	2.45	—	3.31	0.14
		nmol/L	87.2	74.5	—	99.9	4.2		nmol/L	244	208	—	280	12
	BS-480 ⁸	mg/L	1.05	0.89	—	1.21	0.05	BS-480 ⁸	mg/L	2.88	2.45	—	3.31	0.14
		nmol/L	88.9	75.4	—	102.5	4.2		nmol/L	244	208	—	280	12
	BS-600 ⁹	mg/L	1.04	0.88	—	1.20	0.05	BS-600 ⁹	mg/L	2.87	2.44	—	3.30	0.14
		nmol/L	88.1	74.5	—	101.6	4.2		nmol/L	243	207	—	280	12
	BS-600M ¹⁰	mg/L	1.06	0.90	—	1.22	0.05	BS-600M ¹⁰	mg/L	2.86	2.43	—	3.29	0.14
		nmol/L	89.8	76.2	—	103.3	4.2		nmol/L	242	206	—	279	12
	BS-800 ¹¹	mg/L	1.05	0.89	—	1.21	0.05	BS-800 ¹¹	mg/L	2.87	2.44	—	3.30	0.14
		nmol/L	88.9	75.4	—	102.5	4.2		nmol/L	243	207	—	280	12
	BS-2000 ¹²	mg/L	1.05	0.89	—	1.21	0.05	BS-2000 ¹²	mg/L	2.97	2.52	—	3.42	0.15
		nmol/L	88.9	75.4	—	102.5	4.2		nmol/L	252	213	—	290	13
	BS-2800M ¹³	mg/L	1.06	0.90	—	1.22	0.05	BS-2800M ¹³	mg/L	2.86	2.43	—	3.29	0.14
		nmol/L	89.8	76.2	—	103.3	4.2		nmol/L	242	206	—	279	12

β2-Microglobulin Control



English	Abbreviated name	β2-Microglobulin Control	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	Контроль β2-микροглобулина	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Controlo β2-Microglobulina	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	Control β2-microglobulina	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	Controllo β2 Microglobulina	modelli	Unità	Valori di dosaggio	Intervallo(valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	β2-Mikroglobulin Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

β2-MG	
English	β2-Microglobulin
Русский	β2-микροглобулина
Português	β2-Microglobulina
Español	β2-microglobulina
Italiano	β2-Microglobulina
Türkçe	β2-Mikroglobulin