

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

- 1.BS-200E: BS-200E, BS-220E,;
- 2.BS-240: BS-230, BS-240;
- 3.BS-240E: BS240E, BS240Pro;
- 4.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00");
- 5.BS-380: BS-380, BS-390;
- 6.BS-400: BS-400, BS-420;
- 7.BS-430: BS-430, BS-450,BS-460;

- 8.BS-480: BS-480, BS-490;
- 9.BS-600: BS-600, BS-620;
- 10.BS-600M: BS-600M, BS-620M;
- 11.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;
- 12.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;
- 13.BS-2800M:BS-2600M, BS-2800M;
- 14.For applicable analyzers of the analyte, please subject to the parameter sheet and instrument.

LOT : 059124002

Expires: 2025-07-14

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

β2-MG

English	β2-Microglobulin
Русский	β2-микрोगлобулина
Português	β2-Microglobulina
Español	β2-microglobulina
Italiano	β2-Microglobulina
Türkçe	β2-Mikroglobulin

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
β2-MG II	mg/L	BS-200E ¹	0.962	0.818	—	1.106	0.048	BS-200E ¹	2.67	2.28	—	3.06	0.13
		BS-240 ²	0.974	0.827	—	1.121	0.049	BS-240 ²	2.66	2.27	—	3.05	0.13
		BS-240E ³	0.969	0.825	—	1.113	0.048	BS-240E ³	2.64	2.25	—	3.03	0.13
		BS-360E ⁴	0.961	0.817	—	1.105	0.048	BS-360E ⁴	2.63	2.24	—	3.02	0.13
		BS-380 ⁵	0.962	0.818	—	1.106	0.048	BS-380 ⁵	2.67	2.28	—	3.06	0.13
		BS-400 ⁶	0.962	0.818	—	1.106	0.048	BS-400 ⁶	2.67	2.28	—	3.06	0.13
		BS-430 ⁷	0.966	0.822	—	1.110	0.048	BS-430 ⁷	2.67	2.28	—	3.06	0.13
		BS-480 ⁸	0.969	0.825	—	1.113	0.048	BS-480 ⁸	2.70	2.28	—	3.12	0.14
		BS-600 ⁹	0.952	0.808	—	1.096	0.048	BS-600 ⁹	2.62	2.23	—	3.01	0.13
		BS-600M ¹⁰	0.961	0.817	—	1.105	0.048	BS-600M ¹⁰	2.65	2.26	—	3.04	0.13
		BS-800 ¹¹	0.944	0.803	—	1.085	0.047	BS-800 ¹¹	2.70	2.28	—	3.12	0.14
		BS-2000 ¹²	0.961	0.817	—	1.105	0.048	BS-2000 ¹²	2.75	2.33	—	3.17	0.14
		BS-2800M ¹³	0.961	0.817	—	1.105	0.048	BS-2800M ¹³	2.65	2.26	—	3.04	0.13
BS-200E ¹	81.5	69.3	—	93.7	4.1	BS-200E ¹	226	193	—	259	11		
BS-240 ²	82.5	70.0	—	94.9	4.2	BS-240 ²	225	192	—	258	11		
BS-240E ³	82.1	69.9	—	94.3	4.1	BS-240E ³	224	191	—	257	11		
BS-360E ⁴	81.4	69.2	—	93.6	4.1	BS-360E ⁴	223	190	—	256	11		
BS-380 ⁵	81.5	69.3	—	93.7	4.1	BS-380 ⁵	226	193	—	259	11		

β2-Microglobulin Control



nmol/L	BS-400⁶	81.5	69.3	—	93.7	4.1	BS-400⁶	226	193	—	259	11
	BS-430⁷	81.8	69.6	—	94.0	4.1	BS-430⁷	226	193	—	259	11
	BS-480⁸	82.1	69.9	—	94.3	4.1	BS-480⁸	229	193	—	264	12
	BS-600⁹	80.6	68.4	—	92.8	4.1	BS-600⁹	222	189	—	255	11
	BS-600M¹⁰	81.4	69.2	—	93.6	4.1	BS-600M¹⁰	224	191	—	257	11
	BS-800¹¹	80.0	68.0	—	91.9	4.0	BS-800¹¹	229	193	—	264	12
	BS-2000¹²	81.4	69.2	—	93.6	4.1	BS-2000¹²	233	197	—	268	12
	BS-2800M¹³	81.4	69.2	—	93.6	4.1	BS-2800M¹³	224	191	—	257	11

mindray

β 2-Microglobulin Control

For use on: BS-2000, BS-2800M,
BS-600M, BS-800

LOT 059124002

🕒 2025-07-14

