

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

1. **BS-200E**: BS-200E, BS-220E;

2. **BS-230**: BS-230, BS-240;

3. **BS-240E**: BS-240E, BS-240Pro;

4. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

BS-330E(V35.00)(Serial Number starts with "W8-"

and software version starts with "35.00"),

BS-350E(V35.00) (Serial Number starts with "W9-"

and software version starts with "35.00");

5. **BS-380**: BS-380, BS-390;

Русский: Данные совпадают во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

6. **BS-400**: BS-400, BS-420;

7. **BS-430**: BS-430, BS-450, BS-460;

8. **BS-480**: BS-480, BS-490;

9. **BS-600**: BS-600, BS-620;

10. **BS-600M**: BS-600M, BS-620M;

11. **BS-800**: BS-800, BS-820, BS-800M, BS-820M,
BS-1800, BS-1800plus;

12. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

13. **BS-2800M**: BS-2600M, BS-2800M.

14. **Control target value of RF reagent kit corresponds to RF III, and control target value of RF II reagent kit corresponds to RF II.**

LOT: 059222001

EXP: 2023-06-10

Abbreviated name	Unit	Model	Assay Value	Range (Assay Value±3SD)	1 SD	Model	Assay Value	Range (Assay Value±3SD)	1 SD
ASO II	IU/mL	BS-200E ¹	/	/ - /	/	BS-480 ⁸	/	/ - /	/
		BS-230 ²	/	/ - /	/	BS-600 ⁹	/	/ - /	/
		BS-240E ³	/	/ - /	/	BS-600M ¹⁰	104	67 - 141	12
		BS-360E ⁴	/	/ - /	/	BS-800 ¹¹	/	/ - /	/
		BS-380 ⁵	/	/ - /	/	BS-2000 ¹²	/	/ - /	/
		BS-400 ⁶	/	/ - /	/	BS-2800M ¹³	/	/ - /	/
		BS-430 ⁷	/	/ - /	/				
RF II	IU/mL	BS-200E ¹	/	/ - /	/	BS-480 ⁸	/	/ - /	/
		BS-230 ²	/	/ - /	/	BS-600 ⁹	/	/ - /	/
		BS-240E ³	/	/ - /	/	BS-600M ¹⁰	26.2	17.0 - 35.4	3.1
		BS-360E ⁴	/	/ - /	/	BS-800 ¹¹	/	/ - /	/
		BS-380 ⁵	/	/ - /	/	BS-2000 ¹²	/	/ - /	/
		BS-400 ⁶	/	/ - /	/	BS-2800M ¹³	/	/ - /	/
		BS-430 ⁷	/	/ - /	/				
RF III	IU/mL	BS-800 ¹¹	/	/ - /	/	BS-2800M ¹³	/	/ - /	/
		BS-2000 ¹²	/	/ - /	/				

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

ASO			RF		
English	Antistreptolysin "O"		Rheumatoid Factor		
Русский	антистрептолизина O		ревматоидные факторы		
Português	Antiestreptolisina "O"		Fator Reumatoide		
Español	antiestreptolisina "O"		factor reumatoide		
Italiano	Anti-Streptolisina "O"		fattore reumatoide		
Türkçe	Antistreptolisin "O"		Romatoid Faktör		

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray
 Português: Planilha de valores do Sistema de Medição BS da Mindray
 Español: Hoja de valores del sistema de medición Mindray BS
 Italiano: Scheda dei valori del sistema di misurazione BS di Mindray
 Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

1. **BS-200E**: BS-200E, BS-220E;

2. **BS-230**: BS-230, BS-240;

3. **BS-240E**: BS-240E, BS-240Pro;

4. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

BS-330E(V35.00)(Serial Number starts with "W8-"
and software version starts with "35.00");

BS-350E(V35.00) (Serial Number starts with "W9-"
and software version starts with "35.00");

5. **BS-380**: BS-380, BS-390;

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

Español: la datos de cada grupo es la misma.

6. **BS-400**: BS-400, BS-420;

7. **BS-430**: BS-430, BS-450, BS-460;

8. **BS-480**: BS-480, BS-490;

9. **BS-600**: BS-600, BS-620;

10. **BS-600M**: BS-600M, BS-620M;

11. **BS-800**: BS-800, BS-820, BS-800M, BS-820M,
BS-1800, BS-1800plus;

12. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

13. **BS-2800M**: BS-2600M, BS-2800M.

14. **Control target value of RF reagent kit corresponds to RF III, and control target value of RF II reagent kit corresponds to RF II.**

LOT: 059222001

EXP: 2023-06-10

Abbreviated name	Unit	Model	Assay Value	Range (Assay Value \pm 3SD)	1 SD	Model	Assay Value	Range (Assay Value \pm 3SD)	1 SD
ASO II	IU/mL	BS-200E ¹	/	/ - /	/	BS-480 ⁸	/	/ - /	/
		BS-230 ²	/	/ - /	/	BS-600 ⁹	/	/ - /	/
		BS-240E ³	/	/ - /	/	BS-600M ¹⁰	340	221 - 459	40
		BS-360E ⁴	/	/ - /	/	BS-800 ¹¹	/	/ - /	/
		BS-380 ⁵	/	/ - /	/	BS-2000 ¹²	/	/ - /	/
		BS-400 ⁶	/	/ - /	/	BS-2800M ¹³	/	/ - /	/
		BS-430 ⁷	/	/ - /	/				
RF II	IU/mL	BS-200E ¹	/	/ - /	/	BS-480 ⁸	/	/ - /	/
		BS-230 ²	/	/ - /	/	BS-600 ⁹	/	/ - /	/
		BS-240E ³	/	/ - /	/	BS-600M ¹⁰	72.1	46.8 - 97.4	8.4
		BS-360E ⁴	/	/ - /	/	BS-800 ¹¹	/	/ - /	/
		BS-380 ⁵	/	/ - /	/	BS-2000 ¹²	/	/ - /	/
		BS-400 ⁶	/	/ - /	/	BS-2800M ¹³	/	/ - /	/
		BS-430 ⁷	/	/ - /	/				
RF III	IU/mL	BS-800 ¹¹	/	/ - /	/	BS-2800M ¹³	/	/ - /	/
		BS-2000 ¹²	/	/ - /	/				

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа \pm 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise \pm 3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo \pm 3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione \pm 3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri \pm 3SD)

ASO			RF		
English	Antistreptolysin "O"		Rheumatoid Factor		
Русский	антistreптолизина O		ревматоидные факторы		
Português	Antistreptolisina "O"		Fator Reumatoide		
Español	antistreptolisina "O"		factor reumatoide		
Italiano	Anti-Streptolisina "O"		fattore reumatoide		
Türkçe	Antistreptolisin "O"		Romatoid Faktör		