

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

1. **BS-200E**: BS-200E, BS-220E;

2. **BS-230**: BS-230, BS-240;

3. **BS-240E**: BS-240E, BS-240Pro;

4. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

BS-330E(V35.00)(Serial Number starts with "W8-"

and software version starts with "35.00"),

BS-350E(V35.00) (Serial Number starts with "W9-"

and software version starts with "35.00");

5. **BS-380**: BS-380, BS-390;

Русский: Данные совпадают во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

6. **BS-400**: BS-400, BS-420;

7. **BS-430**: BS-430, BS-450, BS-460;

8. **BS-480**: BS-480, BS-490;

9. **BS-600**: BS-600, BS-620;

10. **BS-600M**: BS-600M, BS-620M;

11. **BS-800**: BS-800, BS-820, BS-800M, BS-820M,
BS-1800, BS-1800plus;

12. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

13. **BS-2800M**: BS-2600M, BS-2800M.

14. **Control target value of RF reagent kit corresponds to RF III, and control target value of RF II reagent kit corresponds to RF II.**

LOT: 059222003

Exp. Date: 2023-08-21

Abbreviated name	Unit	Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
ASO II	IU/mL	BS-200E ¹	114	74	—	154	13	BS-480 ⁸	112	73	—	151	13
		BS-230 ²	113	73	—	153	13	BS-600 ⁹	112	73	—	151	13
		BS-240E ³	112	73	—	151	13	BS-600M ¹⁰	110	71	—	149	13
		BS-360E ⁴	112	73	—	151	13	BS-800 ¹¹	112	73	—	151	13
		BS-380 ⁵	114	74	—	154	13	BS-2000 ¹²	110	71	—	149	13
		BS-400 ⁶	114	74	—	154	13	BS-2800M ¹³	110	71	—	149	13
		BS-430 ⁷	112	73	—	151	13						
RF II	IU/mL	BS-200E ¹	27.5	17.8	—	37.2	3.2	BS-480 ⁸	27.6	17.9	—	37.3	3.2
		BS-230 ²	27.2	17.7	—	36.7	3.2	BS-600 ⁹	28.2	18.3	—	38.1	3.3
		BS-240E ³	27.9	18.1	—	37.7	3.3	BS-600M ¹⁰	26.5	17.2	—	35.8	3.1
		BS-360E ⁴	27.6	17.9	—	37.3	3.2	BS-800 ¹¹	27.6	17.9	—	37.3	3.2
		BS-380 ⁵	27.2	17.7	—	36.7	3.2	BS-2000 ¹²	28.9	18.8	—	39.0	3.4
		BS-400 ⁶	27.2	17.7	—	36.7	3.2	BS-2800M ¹³	/	/	—	/	/
		BS-430 ⁷	27.6	17.9	—	37.3	3.2						
RF III	IU/mL	BS-800 ¹¹	39.6	25.7	—	53.5	4.6	BS-2800M ¹³	37.7	24.5	—	50.9	4.4
		BS-2000 ¹²	39.6	25.7	—	53.5	4.6						

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

ASO		RF	
English	Antistreptolysin "O"	Rheumatoid Factor	
Русский	антистрептолизина O	ревматоидные факторы	
Português	Antiestreptolisina "O"	Fator Reumatoide	
Español	antiestreptolisina "O"	factor reumatoide	
Italiano	Anti-Streptolisina "O"	fattore reumatoide	
Türkçe	Antistreptolisin "O"	Romatoid Faktör	

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray
 Português: Planilha de valores do Sistema de Medição BS da Mindray
 Español: Hoja de valores del sistema de medición Mindray BS
 Italiano: Scheda dei valori del sistema di misurazione BS di Mindray
 Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

1. **BS-200E**: BS-200E, BS-220E;

2. **BS-230**: BS-230, BS-240;

3. **BS-240E**: BS-240E, BS-240Pro;

4. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

BS-330E(V35.00)(Serial Number starts with "W8-"
and software version starts with "35.00"),

BS-350E(V35.00) (Serial Number starts with "W9-"
and software version starts with "35.00");

5. **BS-380**: BS-380, BS-390;

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

Español: la datos de cada grupo es la misma.

6. **BS-400**: BS-400, BS-420;

7. **BS-430**: BS-430, BS-450, BS-460;

8. **BS-480**: BS-480, BS-490;

9. **BS-600**: BS-600, BS-620;

10. **BS-600M**: BS-600M, BS-620M;

11. **BS-800**: BS-800, BS-820, BS-800M, BS-820M,
BS-1800, BS-1800plus;

12. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

13. **BS-2800M**: BS-2600M, BS-2800M.

14. **Control target value of RF reagent kit corresponds to RF III, and control target value of RF II reagent kit corresponds to RF II.**

LOT: 059222003

: 2023-08-21

Abbreviated name	Unit	Model	Assay Value	Range (Assay Value±3SD)			1 SD	Model	Assay Value	Range (Assay Value±3SD)			1 SD
ASO II	IU/mL	BS-200E ¹	337	219	—	455	39	BS-480 ⁸	335	217	—	453	39
		BS-230 ²	332	215	—	449	39	BS-600 ⁹	335	217	—	453	39
		BS-240E ³	335	217	—	453	39	BS-600M ¹⁰	331	215	—	447	39
		BS-360E ⁴	335	217	—	453	39	BS-800 ¹¹	335	217	—	453	39
		BS-380 ⁵	337	219	—	455	39	BS-2000 ¹²	333	216	—	450	39
		BS-400 ⁶	337	219	—	455	39	BS-2800M ¹³	331	215	—	447	39
		BS-430 ⁷	335	217	—	453	39						
RF II	IU/mL	BS-200E ¹	72.2	46.9	—	97.5	8.4	BS-480 ⁸	71.1	46.1	—	96.1	8.3
		BS-230 ²	71.2	46.2	—	96.2	8.3	BS-600 ⁹	71.4	46.3	—	96.5	8.4
		BS-240E ³	74.8	48.5	—	101.1	8.8	BS-600M ¹⁰	67.2	43.6	—	90.8	7.9
		BS-360E ⁴	71.1	46.1	—	96.1	8.3	BS-800 ¹¹	71.1	46.1	—	96.1	8.3
		BS-380 ⁵	70.0	45.4	—	94.6	8.2	BS-2000 ¹²	71.9	46.7	—	97.1	8.4
		BS-400 ⁶	70.0	45.4	—	94.6	8.2	BS-2800M ¹³	/	/	—	/	/
		BS-430 ⁷	71.1	46.1	—	96.1	8.3						
RF III	IU/mL	BS-800 ¹¹	95.3	61.8	—	128.8	11.2	BS-2800M ¹³	94.9	61.6	—	128.2	11.1
		BS-2000 ¹²	96.8	62.8	—	130.8	11.3						

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

ASO			RF		
English	Antistreptolysin "O"		Rheumatoid Factor		
Русский	антistreптолизина O		ревматоидные факторы		
Português	Antistreptolisina "O"		Fator Reumatoide		
Español	antistreptolisina "O"		factor reumatoide		
Italiano	Anti-Streptolisina "O"		fattore reumatoide		
Türkçe	Antistreptolisın "O"		Romatoid Faktör		