

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

1. **BS-200E**: BS-200E, BS-220E;

2. **BS-230**: BS-230, BS-240;

3. **BS-240E**: BS-240E, BS-240Pro;

4. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

BS-330E(V35.00)(Serial Number starts with "W8-"

and software version starts with "35.00"),

BS-350E(V35.00) (Serial Number starts with "W9-"

and software version starts with "35.00");

5. **BS-380**: BS-380, BS-390;

Русский: Данные совпадают во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

6. **BS-400**: BS-400, BS-420;

7. **BS-430**: BS-430, BS-450, BS-460;

8. **BS-480**: BS-480, BS-490;

9. **BS-600**: BS-600, BS-620;

10. **BS-600M**: BS-600M, BS-620M;

11. **BS-800**: BS-800, BS-820, BS-800M, BS-820M,

BS-1800, BS-1800plus;

12. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

13. **BS-2800M**: BS-2600M, BS-2800M.

14. **Control target value of RF reagent kit corresponds to RF III, and control target value of RF II reagent kit corresponds to RF II.**

LOT: 059222005

Exp. Date: 2024-02-18

Abbreviated name	Unit	Model	Assay Value	Range (Assay Value±3SD)			1 SD	Model	Assay Value	Range (Assay Value±3SD)			1 SD
ASO II	IU/mL	BS-200E ¹	142	92	—	192	17	BS-480 ⁸	142	92	—	192	17
		BS-230 ²	143	93	—	193	17	BS-600 ⁹	142	92	—	192	17
		BS-240E ³	142	92	—	192	17	BS-600M ¹⁰	135	88	—	182	16
		BS-360E ⁴	142	92	—	192	17	BS-800 ¹¹	142	92	—	192	17
		BS-380 ⁵	142	92	—	192	17	BS-2000 ¹²	136	88	—	184	16
		BS-400 ⁶	142	92	—	192	17	BS-2800M ¹³	135	88	—	182	16
		BS-430 ⁷	142	92	—	192	17						
RF II	IU/mL	BS-200E ¹	34.9	22.7	—	47.1	4.1	BS-480 ⁸	37.4	24.3	—	50.5	4.4
		BS-230 ²	38.5	25.0	—	52.0	4.5	BS-600 ⁹	37.8	24.5	—	51.1	4.4
		BS-240E ³	38.4	24.9	—	51.9	4.5	BS-600M ¹⁰	36.6	23.8	—	49.4	4.3
		BS-360E ⁴	37.4	24.3	—	50.5	4.4	BS-800 ¹¹	37.4	24.3	—	50.5	4.4
		BS-380 ⁵	38.3	24.9	—	51.7	4.5	BS-2000 ¹²	37.5	24.3	—	50.7	4.4
		BS-400 ⁶	38.3	24.9	—	51.7	4.5	BS-2800M ¹³	/	/	—	/	/
		BS-430 ⁷	37.4	24.3	—	50.5	4.4						
RF III	IU/mL	BS-800 ¹¹	52.0	33.7	—	70.3	6.1	BS-2800M ¹³	51.3	33.3	—	69.3	6.0
		BS-2000 ¹²	52.5	34.1	—	70.9	6.1						

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

ASO			RF		
English	Antistreptolysin "O"		Rheumatoid Factor		
Русский	антistreптолизина O		ревматоидные факторы		
Português	AntiStreptolisina "O"		Fator Reumatoide		
Español	antiStreptolisina "O"		factor reumatoide		
Italiano	Anti-Streptolisina "O"		fattore reumatoide		
Türkçe	Antistreptolisin "O"		Romatoid Faktör		

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

1. **BS-200E**: BS-200E, BS-220E;

2. **BS-230**: BS-230, BS-240;

3. **BS-240E**: BS-240E, BS-240Pro;

4. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

BS-330E(V35.00)(Serial Number starts with "W8-"

and software version starts with "35.00"),

BS-350E(V35.00) (Serial Number starts with "W9-"

and software version starts with "35.00");

5. **BS-380**: BS-380, BS-390;

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

Español: la datos de cada grupo es la misma.

6. **BS-400**: BS-400, BS-420;

7. **BS-430**: BS-430, BS-450, BS-460;

8. **BS-480**: BS-480, BS-490;

9. **BS-600**: BS-600, BS-620;

10. **BS-600M**: BS-600M, BS-620M;

11. **BS-800**: BS-800, BS-820, BS-800M, BS-820M,
BS-1800, BS-1800plus;

12. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

13. **BS-2800M**: BS-2600M, BS-2800M.

14. Control target value of RF reagent kit corresponds to RF III, and control target value of RF II reagent kit corresponds to RF II.

LOT: 059222005

EXP: 2024-02-18

Abbreviated name	Unit	Model	Assay Value	Range (Assay Value±3SD)			1 SD	Model	Assay Value	Range (Assay Value±3SD)			1 SD
ASO II	IU/mL	BS-200E ¹	356	231	—	481	42	BS-480 ⁸	354	230	—	478	41
		BS-230 ²	345	224	—	466	40	BS-600 ⁹	354	230	—	478	41
		BS-240E ³	354	230	—	478	41	BS-600M ¹⁰	346	225	—	467	40
		BS-360E ⁴	354	230	—	478	41	BS-800 ¹¹	354	230	—	478	41
		BS-380 ⁵	356	231	—	481	42	BS-2000 ¹²	346	225	—	467	40
		BS-400 ⁶	356	231	—	481	42	BS-2800M ¹³	346	225	—	467	40
		BS-430 ⁷	354	230	—	478	41						
RF II	IU/mL	BS-200E ¹	105	68	—	142	12	BS-480 ⁸	93.6	60.7	—	126.5	11.0
		BS-230 ²	93.7	60.8	—	126.6	11.0	BS-600 ⁹	93.2	60.5	—	125.9	10.9
		BS-240E ³	93.4	60.6	—	126.2	10.9	BS-600M ¹⁰	92.5	60.0	—	125.0	10.8
		BS-360E ⁴	93.6	60.7	—	126.5	11.0	BS-800 ¹¹	93.6	60.7	—	126.5	11.0
		BS-380 ⁵	94.5	61.3	—	127.7	11.1	BS-2000 ¹²	94.6	61.4	—	127.8	11.1
		BS-400 ⁶	94.5	61.3	—	127.7	11.1	BS-2800M ¹³	/	/	—	/	/
		BS-430 ⁷	93.6	60.7	—	126.5	11.0						
RF III	IU/mL	BS-800 ¹¹	110	71	—	149	13	BS-2800M ¹³	106	69	—	143	12
		BS-2000 ¹²	113	73	—	153	13						

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

ASO			RF		
English	Antistreptolysin "O"		Rheumatoid Factor		
Русский	антистрептолизина O		ревматоидные факторы		
Português	Antiestreptolisina "O"		Fator Reumatoide		
Español	antiestreptolisina "O"		factor reumatoide		
Italiano	Anti-Streptolisina "O"		fattore reumatoide		
Türkçe	Antistreptolisin "O"		Romatoid Faktör		