

Value sheet of Mindray BS Measurement System

Русский : Таблица результатов для системы BS компании Mindray
 Português : Planilha de valores do Sistema de Medição BS da Mindray
 Español : Hoja de valores del sistema de medición Mindray BS
 Italiano : Scheda dei valori del sistema di misurazione BS di Mindray
 Türkçe : Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Português : A dados de cada grupo é a mesma.

Italiano : la dati di ogni gruppo è la stessa.

1.**BS-200E**: BS-200E, BS-220E;

2.**BS-230**: BS-230, BS-240;

3.**BS-240E**: BS-240E, BS-240Pro;

4.**BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

BS-330E(V35.00)(Serial Number starts with "W8-"

and software version starts with "35.00"),

BS-350E(V35.00) (Serial Number starts with"W9-"

and software version starts with "35.00");

5.**BS-380**: BS-380, BS-390;

Русский: Данные совпадают во всех группах.

Español : la datos de cada grupo es la misma.

Türkçe : her grubun veri aynıdır.

6.**BS-400**: BS-400, BS-420;

7.**BS-430**: BS-430, BS-450, BS-460;

8.**BS-480**: BS-480, BS-490;

9.**BS-600**: BS-600, BS-620;

10.**BS-600M**: BS-600M, BS-620M;

11.**BS-800**: BS-800, BS-820, BS-800M, BS-820M,
BS-1800, BS-1800plus;

12.**BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

13.**BS-2800M**:BS-2600M, BS-2800M.

14.**Control target value of RF reagent kit corresponds to RF III, and control target value of RF II reagent kit corresponds to RF II.**

LOT : 059223001

EXP : 2024-08-13

Abbreviated name	Unit	Model	Assay Value	Range (Assay Value±3SD)			1 SD	Model	Assay Value	Range (Assay Value±3SD)			1 SD
ASO II	IU/mL	BS-200E ¹	128	83	—	173	15	BS-480 ⁸	123	80	—	166	14
		BS-230 ²	126	82	—	170	15	BS-600 ⁹	123	80	—	166	14
		BS-240E ³	123	80	—	166	14	BS-600M ¹⁰	125	81	—	169	15
		BS-360E ⁴	123	80	—	166	14	BS-800 ¹¹	123	80	—	166	14
		BS-380 ⁵	128	83	—	173	15	BS-2000 ¹²	121	79	—	163	14
		BS-400 ⁶	128	83	—	173	15	BS-2800M ¹³	125	81	—	169	15
		BS-430 ⁷	123	80	—	166	14						
RF II	IU/mL	BS-200E ¹	35.6	23.1	—	48.1	4.2	BS-480 ⁸	34.1	22.1	—	46.1	4.0
		BS-230 ²	33.6	21.8	—	45.4	3.9	BS-600 ⁹	34.8	22.6	—	47.0	4.1
		BS-240E ³	34.4	22.3	—	46.5	4.0	BS-600M ¹⁰	33.5	21.7	—	45.3	3.9
		BS-360E ⁴	34.1	22.1	—	46.1	4.0	BS-800 ¹¹	34.1	22.1	—	46.1	4.0
		BS-380 ⁵	34.4	22.3	—	46.5	4.0	BS-2000 ¹²	34.2	22.2	—	46.2	4.0
		BS-400 ⁶	34.4	22.3	—	46.5	4.0	BS-2800M ¹³	/	/	—	/	/
		BS-430 ⁷	34.1	22.1	—	46.1	4.0						
RF III	IU/mL	BS-800 ¹¹	44.7	29.0	—	60.4	5.2	BS-2800M ¹³	44.3	28.8	—	59.8	5.2
		BS-2000 ¹²	45.3	29.4	—	61.2	5.3						

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

ASO			RF		
English	Antistreptolysin "O"		Rheumatoid Factor		
Русский	антistreптолизина O		ревматоидные факторы		
Português	Antiestreptolisina "O"		Fator Reumatoide		
Español	antiestreptolisina "O"		factor reumatoide		
Italiano	Anti-Streptolisina "O"		fattore reumatoide		
Türkçe	Antistreptolisin "O"		Romatoid Faktör		

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray
 Português: Planilha de valores do Sistema de Medição BS da Mindray
 Español: Hoja de valores del sistema de medición Mindray BS
 Italiano: Scheda dei valori del sistema di misurazione BS di Mindray
 Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

1. **BS-200E**: BS-200E, BS-220E;

2. **BS-230**: BS-230, BS-240;

3. **BS-240E**: BS-240E, BS-240Pro;

4. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

BS-330E(V35.00)(Serial Number starts with "W8-"
and software version starts with "35.00"),

BS-350E(V35.00) (Serial Number starts with "W9-"
and software version starts with "35.00");

5. **BS-380**: BS-380, BS-390;

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

Español: la datos de cada grupo es la misma.

6. **BS-400**: BS-400, BS-420;

7. **BS-430**: BS-430, BS-450, BS-460;

8. **BS-480**: BS-480, BS-490;

9. **BS-600**: BS-600, BS-620;

10. **BS-600M**: BS-600M, BS-620M;

11. **BS-800**: BS-800, BS-820, BS-800M, BS-820M,
BS-1800, BS-1800plus;

12. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

13. **BS-2800M**: BS-2600M, BS-2800M.

14. **Control target value of RF reagent kit corresponds to RF III, and control target value of RF II reagent kit corresponds to RF II.**

LOT: 059223001

2024-08-13

Abbreviated name	Unit	Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
ASO II	IU/mL	BS-200E ¹	332	215	—	449	39	BS-480 ⁸	328	213	—	443	38
		BS-230 ²	335	217	—	453	39	BS-600 ⁹	328	213	—	443	38
		BS-240E ³	328	213	—	443	38	BS-600M ¹⁰	329	214	—	444	38
		BS-360E ⁴	328	213	—	443	38	BS-800 ¹¹	328	213	—	443	38
		BS-380 ⁵	332	215	—	449	39	BS-2000 ¹²	320	208	—	432	37
		BS-400 ⁶	332	215	—	449	39	BS-2800M ¹³	329	214	—	444	38
		BS-430 ⁷	328	213	—	443	38						
RF II	IU/mL	BS-200E ¹	85.8	55.7	—	115.9	10.0	BS-480 ⁸	80.3	52.1	—	108.5	9.4
		BS-230 ²	79.4	51.5	—	107.3	9.3	BS-600 ⁹	83.3	54.1	—	112.5	9.7
		BS-240E ³	84.9	55.1	—	114.7	9.9	BS-600M ¹⁰	78.7	51.1	—	106.3	9.2
		BS-360E ⁴	80.3	52.1	—	108.5	9.4	BS-800 ¹¹	80.3	52.1	—	108.5	9.4
		BS-380 ⁵	80.2	52.0	—	108.4	9.4	BS-2000 ¹²	80.0	51.9	—	108.1	9.4
		BS-400 ⁶	80.2	52.0	—	108.4	9.4	BS-2800M ¹³	/	/	—	/	/
		BS-430 ⁷	80.3	52.1	—	108.5	9.4						
RF III	IU/mL	BS-800 ¹¹	104	67	—	141	12	BS-2800M ¹³	103	67	—	139	12
		BS-2000 ¹²	123	80	—	166	14						

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

	ASO	RF
English	Antistreptolysin "O"	Rheumatoid Factor
Русский	антистрептолизина O	ревматоидные факторы
Português	AntiStreptolisina "O"	Fator Reumatoide
Español	antiestreptolisina "O"	factor reumatoide
Italiano	Anti-Streptolisina "O"	fattore reumatoide
Türkçe	Antistreptolisin "O"	Romatoid Faktör