

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

1. **BS-200E**: BS-200E, BS-220E;2. **BS-230**: BS-230, BS-240;3. **BS-240E**: BS-240E, BS-240Pro;4. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,**BS-330E(V35.00)**(Serial Number starts with "W8-"

and software version starts with "35.00"),

BS-350E(V35.00) (Serial Number starts with "W9-"

and software version starts with "35.00");

5. **BS-380**: BS-380, BS-390;

Русский: Данные совпадают во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

6. **BS-400**: BS-400, BS-420;7. **BS-430**: BS-430, BS-450, BS-460;8. **BS-480**: BS-480, BS-490;9. **BS-600**: BS-600, BS-620;10. **BS-600M**: BS-600M, BS-620M;11. **BS-800**: BS-800, BS-820, BS-800M, BS-820M,
BS-1800, BS-1800plus;12. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;13. **BS-2800M**: BS-2600M, BS-2800M.

14. Control target value of RF reagent kit corresponds to RF III, and control target value of RF II reagent kit corresponds to RF II.

LOT : 059223002

: 2024-08-13

Abbreviated name	Unit	Model	Assay Value	Range (Assay Value \pm 3SD)			1 SD	Model	Assay Value	Range (Assay Value \pm 3SD)			1 SD
ASO II	IU/mL	BS-200E ¹	128	83	—	173	15	BS-480 ⁸	123	80	—	166	14
		BS-230 ²	126	82	—	170	15	BS-600 ⁹	123	80	—	166	14
		BS-240E ³	123	80	—	166	14	BS-600M ¹⁰	125	81	—	169	15
		BS-360E ⁴	123	80	—	166	14	BS-800 ¹¹	123	80	—	166	14
		BS-380 ⁵	128	83	—	173	15	BS-2000 ¹²	121	79	—	163	14
		BS-400 ⁶	128	83	—	173	15	BS-2800M ¹³	125	81	—	169	15
		BS-430 ⁷	123	80	—	166	14						
RF II	IU/mL	BS-200E ¹	35.6	23.1	—	48.1	4.2	BS-480 ⁸	34.1	22.1	—	46.1	4.0
		BS-230 ²	33.6	21.8	—	45.4	3.9	BS-600 ⁹	34.8	22.6	—	47.0	4.1
		BS-240E ³	34.4	22.3	—	46.5	4.0	BS-600M ¹⁰	33.5	21.7	—	45.3	3.9
		BS-360E ⁴	34.1	22.1	—	46.1	4.0	BS-800 ¹¹	34.1	22.1	—	46.1	4.0
		BS-380 ⁵	34.4	22.3	—	46.5	4.0	BS-2000 ¹²	34.2	22.2	—	46.2	4.0
		BS-400 ⁶	34.4	22.3	—	46.5	4.0	BS-2800M ¹³	/	/	—	/	/
		BS-430 ⁷	34.1	22.1	—	46.1	4.0						
RF III	IU/mL	BS-800 ¹¹	44.7	29.0	—	60.4	5.2	BS-2800M ¹³	44.3	28.8	—	59.8	5.2
		BS-2000 ¹²	45.3	29.4	—	61.2	5.3						

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа \pm 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise \pm 3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo \pm 3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione \pm 3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri \pm 3SD)

ASO			RF		
English	Antistreptolysin "O"		Rheumatoid Factor		
Русский	антistreптолизина O		ревматоидные факторы		
Português	Antiestreptolisina "O"		Fator Reumatoide		
Español	antiestreptolisina "O"		factor reumatoide		
Italiano	Anti-Streptolisina "O"		fattore reumatoide		
Türkçe	Antistreptolisin "O"		Romatoid Faktör		

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray
 Português: Planilha de valores do Sistema de Medição BS da Mindray
 Español: Hoja de valores del sistema de medición Mindray BS
 Italiano: Scheda dei valori del sistema di misurazione BS di Mindray
 Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

1. **BS-200E**: BS-200E, BS-220E;

2. **BS-230**: BS-230, BS-240;

3. **BS-240E**: BS-240E, BS-240Pro;

4. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

BS-330E(V35.00)(Serial Number starts with "W8-"
and software version starts with "35.00"),

BS-350E(V35.00) (Serial Number starts with "W9-"
and software version starts with "35.00");

5. **BS-380**: BS-380, BS-390;

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

Español: la datos de cada grupo es la misma.

6. **BS-400**: BS-400, BS-420;

7. **BS-430**: BS-430, BS-450, BS-460;

8. **BS-480**: BS-480, BS-490;

9. **BS-600**: BS-600, BS-620;

10. **BS-600M**: BS-600M, BS-620M;

11. **BS-800**: BS-800, BS-820, BS-800M, BS-820M,
BS-1800, BS-1800plus;

12. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

13. **BS-2800M**: BS-2600M, BS-2800M.

14. Control target value of RF reagent kit corresponds to RF III, and control target value of RF II reagent kit corresponds to RF II.

LOT: 059223002

EXP: 2024-08-13

Abbreviated name	Unit	Model	Assay Value	Range (Assay Value \pm 3SD)		1 SD	Model	Assay Value	Range (Assay Value \pm 3SD)		1 SD		
ASO II	IU/mL	BS-200E ¹	332	215	—	449	39	BS-480 ⁸	328	213	—	443	38
		BS-230 ²	335	217	—	453	39	BS-600 ⁹	328	213	—	443	38
		BS-240E ³	328	213	—	443	38	BS-600M ¹⁰	329	214	—	444	38
		BS-360E ⁴	328	213	—	443	38	BS-800 ¹¹	328	213	—	443	38
		BS-380 ⁵	332	215	—	449	39	BS-2000 ¹²	320	208	—	432	37
		BS-400 ⁶	332	215	—	449	39	BS-2800M ¹³	329	214	—	444	38
		BS-430 ⁷	328	213	—	443	38						
RF II	IU/mL	BS-200E ¹	85.8	55.7	—	115.9	10.0	BS-480 ⁸	80.3	52.1	—	108.5	9.4
		BS-230 ²	79.4	51.5	—	107.3	9.3	BS-600 ⁹	83.3	54.1	—	112.5	9.7
		BS-240E ³	84.9	55.1	—	114.7	9.9	BS-600M ¹⁰	78.7	51.1	—	106.3	9.2
		BS-360E ⁴	80.3	52.1	—	108.5	9.4	BS-800 ¹¹	80.3	52.1	—	108.5	9.4
		BS-380 ⁵	80.2	52.0	—	108.4	9.4	BS-2000 ¹²	80.0	51.9	—	108.1	9.4
		BS-400 ⁶	80.2	52.0	—	108.4	9.4	BS-2800M ¹³	/	/	—	/	/
		BS-430 ⁷	80.3	52.1	—	108.5	9.4						
RF III	IU/mL	BS-800 ¹¹	104	67	—	141	12	BS-2800M ¹³	103	67	—	139	12
		BS-2000 ¹²	123	80	—	166	14						

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа \pm 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise \pm 3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo \pm 3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione \pm 3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri \pm 3SD)

	ASO	RF
English	Antistreptolysin "O"	Rheumatoid Factor
Русский	антистрептолизина O	ревматоидные факторы
Português	AntiStreptolisina "O"	Fator Reumatoide
Español	antiestreptolisina "O"	factor reumatoide
Italiano	Anti-Streptolisina "O"	fattore reumatoide
Türkçe	Antistreptolisin "O"	Romatoid Faktör