

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

1. **BS-200E**: BS-200E, BS-220E;

2. **BS-230**: BS-230, BS-240;

3. **BS-240E**: BS-240E, BS-240Pro;

4. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

BS-330E(V35.00)(Serial Number starts with "W8-"

and software version starts with "35.00"),

BS-350E(V35.00) (Serial Number starts with "W9-"

and software version starts with "35.00");

5. **BS-380**: BS-380, BS-390;

Русский: Данные совпадают во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

6. **BS-400**: BS-400, BS-420;

7. **BS-430**: BS-430, BS-450, BS-460;

8. **BS-480**: BS-480, BS-490;

9. **BS-600**: BS-600, BS-620;

10. **BS-600M**: BS-600M, BS-620M;

11. **BS-800**: BS-800, BS-820, BS-800M, BS-820M,
BS-1800, BS-1800plus;

12. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

13. **BS-2800M**: BS-2600M, BS-2800M.

14. **Control target value of RF reagent kit corresponds to RF III, and control target value of RF II reagent kit corresponds to RF II.**

LOT : 059223004

EXP : 2024-11-09

Abbreviated name	Unit	Model	Assay Value	Range (Assay Value±3SD)			1 SD	Model	Assay Value	Range (Assay Value±3SD)			1 SD
ASO II	IU/mL	BS-200E ¹	107	69	—	145	13	BS-480 ⁸	108	70	—	146	13
		BS-230 ²	108	70	—	146	13	BS-600 ⁹	108	70	—	146	13
		BS-240E ³	108	70	—	146	13	BS-600M ¹⁰	104	67	—	141	12
		BS-360E ⁴	108	70	—	146	13	BS-800 ¹¹	108	70	—	146	13
		BS-380 ⁵	107	69	—	145	13	BS-2000 ¹²	104	67	—	141	12
		BS-400 ⁶	107	69	—	145	13	BS-2800M ¹³	104	67	—	141	12
		BS-430 ⁷	108	70	—	146	13						
RF II	IU/mL	BS-200E ¹	34.4	22.3	—	46.5	4.0	BS-480 ⁸	32.9	21.4	—	44.4	3.8
		BS-230 ²	33.2	21.5	—	44.9	3.9	BS-600 ⁹	33.9	22.0	—	45.8	4.0
		BS-240E ³	32.8	21.3	—	44.3	3.8	BS-600M ¹⁰	32.2	20.9	—	43.5	3.8
		BS-360E ⁴	32.9	21.4	—	44.4	3.8	BS-800 ¹¹	32.9	21.4	—	44.4	3.8
		BS-380 ⁵	33.9	22.0	—	45.8	4.0	BS-2000 ¹²	32.7	21.2	—	44.2	3.8
		BS-400 ⁶	33.9	22.0	—	45.8	4.0	BS-2800M ¹³	/	/	—	/	/
		BS-430 ⁷	32.9	21.4	—	44.4	3.8						
RF III	IU/mL	BS-800 ¹¹	41.3	26.8	—	55.8	4.8	BS-2800M ¹³	40.0	26.0	—	54.0	4.7
		BS-2000 ¹²	40.4	26.2	—	54.6	4.7						

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

ASO			RF		
English	Antistreptolysin "O"		Rheumatoid Factor		
Русский	антистрептолизина O		ревматоидные факторы		
Português	Antiestreptolisina "O"		Fator Reumatoide		
Español	antiestreptolisina "O"		factor reumatoide		
Italiano	Anti-Streptolisina "O"		fattore reumatoide		
Türkçe	Antistreptolisin "O"		Romatoid Faktör		

Value sheet of Mindray BS Measurement System

Русский : Таблица результатов для системы BS компании Mindray
 Português : Planilha de valores do Sistema de Medição BS da Mindray
 Español : Hoja de valores del sistema de medición Mindray BS
 Italiano : Scheda dei valori del sistema di misurazione BS di Mindray
 Türkçe : Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский : Данные совпадают во всех группах.

Español : la datos de cada grupo es la misma.

Türkçe : her grubun veri aynıdır.

1.**BS-200E**: BS-200E, BS-220E;

2.**BS-230**: BS-230, BS-240;

3.**BS-240E**: BS-240E, BS-240Pro;

4.**BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

BS-330E(V35.00)(Serial Number starts with "W8-"
and software version starts with "35.00"),

BS-350E(V35.00) (Serial Number starts with "W9-"
and software version starts with "35.00");

5.**BS-380**: BS-380, BS-390;

Português : A dados de cada grupo é a mesma.

Italiano : la dati di ogni gruppo è la stessa.

Español : la datos de cada grupo es la misma.

6.**BS-400**: BS-400, BS-420;

7.**BS-430**: BS-430, BS-450, BS-460;

8.**BS-480**: BS-480, BS-490;

9.**BS-600**: BS-600, BS-620;

10.**BS-600M**: BS-600M, BS-620M;

11.**BS-800**: BS-800, BS-820, BS-800M, BS-820M,
BS-1800, BS-1800plus;

12.**BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

13.**BS-2800M**:BS-2600M, BS-2800M.

14.**Control target value of RF reagent kit corresponds to RF III, and control target value of RF II reagent kit corresponds to RF II.**

LOT : 059223004

EXP : 2024-11-09

Abbreviated name	Unit	Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
ASO II	IU/mL	BS-200E ¹	324	210	—	438	38	BS-480 ⁸	325	211	—	439	38
		BS-230 ²	321	208	—	434	38	BS-600 ⁹	325	211	—	439	38
		BS-240E ³	325	211	—	439	38	BS-600M ¹⁰	319	207	—	431	37
		BS-360E ⁴	325	211	—	439	38	BS-800 ¹¹	325	211	—	439	38
		BS-380 ⁵	324	210	—	438	38	BS-2000 ¹²	318	206	—	430	37
		BS-400 ⁶	324	210	—	438	38	BS-2800M ¹³	319	207	—	431	37
		BS-430 ⁷	325	211	—	439	38						
RF II	IU/mL	BS-200E ¹	84.4	54.8	—	114.0	9.9	BS-480 ⁸	81.7	53.0	—	110.4	9.6
		BS-230 ²	80.4	52.2	—	108.6	9.4	BS-600 ⁹	80.8	52.4	—	109.2	9.5
		BS-240E ³	83.8	54.4	—	113.2	9.8	BS-600M ¹⁰	80.4	52.2	—	108.6	9.4
		BS-360E ⁴	81.7	53.0	—	110.4	9.6	BS-800 ¹¹	81.7	53.0	—	110.4	9.6
		BS-380 ⁵	81.1	52.6	—	109.6	9.5	BS-2000 ¹²	78.6	51.0	—	106.2	9.2
		BS-400 ⁶	81.1	52.6	—	109.6	9.5	BS-2800M ¹³	/	/	—	/	/
		BS-430 ⁷	81.7	53.0	—	110.4	9.6						
RF III	IU/mL	BS-800 ¹¹	105	68	—	142	12	BS-2800M ¹³	104	67	—	141	12
		BS-2000 ¹²	104	67	—	141	12						

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

ASO			RF		
English	Antistreptolysin "O"		Rheumatoid Factor		
Русский	антistreптолизина O		ревматоидные факторы		
Português	Antistreptolisina "O"		Fator Reumatoide		
Español	antistreptolisina "O"		factor reumatoide		
Italiano	Anti-Streptolisina "O"		fattore reumatoide		
Türkçe	Antistreptolisin "O"		Romatoid Faktör		