

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray
 Português: Planilha de valores do Sistema de Medição BS da Mindray
 Español: Hoja de valores del sistema de medición Mindray BS
 Italiano: Scheda dei valori del sistema di misurazione BS di Mindray
 Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

1. **BS-200E**: BS-200E, BS-220E;

2. **BS-230**: BS-230, BS-240;

3. **BS-240E**: BS-240E, BS-240Pro;

4. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

BS-330E(V35.00)(Serial Number starts with "W8-"

and software version starts with "35.00"),

BS-350E(V35.00) (Serial Number starts with "W9-"

and software version starts with "35.00");

5. **BS-380**: BS-380, BS-390;

Русский: Данные совпадают во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

6. **BS-400**: BS-400, BS-420;

7. **BS-430**: BS-430, BS-450, BS-460;

8. **BS-480**: BS-480, BS-490;

9. **BS-600**: BS-600, BS-620;

10. **BS-600M**: BS-600M, BS-620M;

11. **BS-800**: BS-800, BS-820, BS-800M, BS-820M,
BS-1800, BS-1800plus;

12. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

13. **BS-2800M**: BS-2600M, BS-2800M.

14. **Control target value of RF reagent kit corresponds to RF III, and control target value of RF II reagent kit corresponds to RF II.**

LOT: 059223009

EXP: 2025-03-24

Abbreviated name	Unit	Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
ASO II	IU/mL	BS-200E ¹	123	81	—	165	14	BS-480 ⁸	123	81	—	165	14
		BS-230 ²	126	81	—	171	15	BS-600 ⁹	123	81	—	165	14
		BS-240E ³	123	81	—	165	14	BS-600M ¹⁰	119	77	—	161	14
		BS-360E ⁴	123	81	—	165	14	BS-800 ¹¹	123	81	—	165	14
		BS-380 ⁵	123	81	—	165	14	BS-2000 ¹²	119	77	—	161	14
		BS-400 ⁶	123	81	—	165	14	BS-2800M ¹³	119	77	—	161	14
		BS-430 ⁷	123	81	—	165	14						
RF II	IU/mL	BS-200E ¹	35.8	23.2	—	48.4	4.2	BS-480 ⁸	36.3	23.7	—	48.9	4.2
		BS-230 ²	35.0	22.7	—	47.3	4.1	BS-600 ⁹	35.1	22.8	—	47.4	4.1
		BS-240E ³	35.8	23.2	—	48.4	4.2	BS-600M ¹⁰	36.5	23.6	—	49.4	4.3
		BS-360E ⁴	36.3	23.7	—	48.9	4.2	BS-800 ¹¹	36.3	23.7	—	48.9	4.2
		BS-380 ⁵	37.0	24.1	—	49.9	4.3	BS-2000 ¹²	34.2	22.2	—	46.2	4.0
		BS-400 ⁶	37.0	24.1	—	49.9	4.3	BS-2800M ¹³	/	/	—	/	/
		BS-430 ⁷	36.3	23.7	—	48.9	4.2						
RF III	IU/mL	BS-800 ¹¹	44.8	29.2	—	60.4	5.2	BS-2800M ¹³	43.0	28.0	—	58.0	5.0
		BS-2000 ¹²	44.4	28.8	—	60.0	5.2						

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

ASO		RF	
English	Antistreptolysin "O"	Rheumatoid Factor	
Русский	антистрептолизина O	ревматоидные факторы	
Português	Antiestreptolisina "O"	Fator Reumatoide	
Español	antiestreptolisina "O"	factor reumatoide	
Italiano	Anti-Streptolisina "O"	fattore reumatoide	
Türkçe	Antistreptolisin "O"	Romatoid Faktör	

Value sheet of Mindray BS Measurement System

Русский : Таблица результатов для системы BS компании Mindray
 Português : Planilha de valores do Sistema de Medição BS da Mindray
 Español : Hoja de valores del sistema de medición Mindray BS
 Italiano : Scheda dei valori del sistema di misurazione BS di Mindray
 Türkçe : Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский : Данные совпадают во всех группах.

Español : la datos de cada grupo es la misma.

Türkçe : her grubun veri aynıdır.

1.**BS-200E**: BS-200E, BS-220E;

2.**BS-230**: BS-230, BS-240;

3.**BS-240E**: BS-240E, BS-240Pro;

4.**BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00");

BS-350E(V35.00) (Serial Number starts with "W9-" and software version starts with "35.00");

5.**BS-380**: BS-380, BS-390;

Português : A dados de cada grupo é a mesma.

Italiano : la dati di ogni gruppo è la stessa.

Español : la datos de cada grupo es la misma.

6.**BS-400**: BS-400, BS-420;

7.**BS-430**: BS-430, BS-450, BS-460;

8.**BS-480**: BS-480, BS-490;

9.**BS-600**: BS-600, BS-620;

10.**BS-600M**: BS-600M, BS-620M;

11.**BS-800**: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

12.**BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

13.**BS-2800M**:BS-2600M, BS-2800M.

14.**Control target value of RF reagent kit corresponds to RF III, and control target value of RF II reagent kit corresponds to RF II.**

LOT : 059223009



: 2025-03-24

Abbreviated name	Unit	Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
ASO II	IU/mL	BS-200E ¹	328	214	—	442	38	BS-480 ⁸	328	214	—	442	38
		BS-230 ²	331	214	—	448	39	BS-600 ⁹	328	214	—	442	38
		BS-240E ³	328	214	—	442	38	BS-600M ¹⁰	321	207	—	435	38
		BS-360E ⁴	328	214	—	442	38	BS-800 ¹¹	328	214	—	442	38
		BS-380 ⁵	328	214	—	442	38	BS-2000 ¹²	321	207	—	435	38
		BS-400 ⁶	328	214	—	442	38	BS-2800M ¹³	321	207	—	435	38
		BS-430 ⁷	328	214	—	442	38						
RF II	IU/mL	BS-200E ¹	83.0	53.9	—	112.1	9.7	BS-480 ⁸	84.1	54.7	—	113.5	9.8
		BS-230 ²	85.0	55.3	—	114.7	9.9	BS-600 ⁹	79.8	51.9	—	107.7	9.3
		BS-240E ³	84.1	54.7	—	113.5	9.8	BS-600M ¹⁰	83.6	54.2	—	113.0	9.8
		BS-360E ⁴	84.1	54.7	—	113.5	9.8	BS-800 ¹¹	84.1	54.7	—	113.5	9.8
		BS-380 ⁵	85.1	55.1	—	115.1	10.0	BS-2000 ¹²	79.1	51.2	—	107.0	9.3
		BS-400 ⁶	85.1	55.1	—	115.1	10.0	BS-2800M ¹³	/	/	—	/	/
		BS-430 ⁷	84.1	54.7	—	113.5	9.8						
RF III	IU/mL	BS-800 ¹¹	103	67	—	139	12	BS-2800M ¹³	101	65	—	137	12
		BS-2000 ¹²	107	68	—	146	13						

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3СО)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

ASO			RF		
English	Antistreptolysin "O"		Rheumatoid Factor		
Русский	антистрептолизина O		ревматоидные факторы		
Português	Antiestreptolisina "O"		Fator Reumatoide		
Español	antiestreptolisina "O"		factor reumatoide		
Italiano	Anti-Streptolisina "O"		fattore reumatoide		
Türkçe	Antistreptolisin "O"		Romatoid Faktör		