

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.BS-200E: BS-200E, BS-220E,;

8.BS-480: BS-480, BS-490;

2.BS-240: BS-230, BS-240;

9.BS-600: BS-600, BS-620;

3.BS-240E: BS240E, BS240Pro;

10.BS-600M: BS-600M, BS-620M;

4.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-

11.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with "W9-"

12.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

and software version starts with "35.00");

13.BS-2800M:BS-2600M, BS-2800M;

5.BS-380: BS-380, BS-390;

14.Applicable models of the chemistry shall be subject to the parameter list and instructions.

6.BS-400: BS-400, BS-420;

7.BS-430: BS-430, BS-450,BS-460;

LOT : 059224001

: 2025-03-24

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

ASO

RF

English	Antistreptolysin "O"	Rheumatoid Factor
Русский	антистрептолизина O	ревматоидные факторы
Português	Antiestreptolisina "O"	Fator Reumatoide
Español	antiestreptolisina "O"	factor reumatoide
Italiano	Anti-Streptolisina "O"	fattore reumatoide
Türkçe	Antistreptolisin "O"	Romatoid Faktör

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
ASO II	IU/mL	BS-200E ¹	123	81	—	165	14	BS-200E ¹	328	214	—	442	38
		BS-240 ²	126	81	—	171	15	BS-240 ²	331	214	—	448	39
		BS-240E ³	123	81	—	165	14	BS-240E ³	328	214	—	442	38
		BS-360E ⁴	123	81	—	165	14	BS-360E ⁴	328	214	—	442	38
		BS-380 ⁵	123	81	—	165	14	BS-380 ⁵	328	214	—	442	38
		BS-400 ⁶	123	81	—	165	14	BS-400 ⁶	328	214	—	442	38
		BS-430 ⁷	123	81	—	165	14	BS-430 ⁷	328	214	—	442	38
		BS-480 ⁸	123	81	—	165	14	BS-480 ⁸	328	214	—	442	38
		BS-600 ⁹	123	81	—	165	14	BS-600 ⁹	328	214	—	442	38
		BS-600M ¹⁰	119	77	—	161	14	BS-600M ¹⁰	321	207	—	435	38
		BS-800 ¹¹	123	81	—	165	14	BS-800 ¹¹	328	214	—	442	38
		BS-2000 ¹²	119	77	—	161	14	BS-2000 ¹²	321	207	—	435	38
		BS-2800M ¹³	119	77	—	161	14	BS-2800M ¹³	321	207	—	435	38

Abbreviated name	Unit	Level: L/N					Level: H/P					
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD	
BS-200E ¹		35.8	23.2	—	48.4	4.2	BS-200E ¹	83.0	53.9	—	112.1	9.7
BS-240 ²		35.0	22.7	—	47.3	4.1	BS-240 ²	85.0	55.3	—	114.7	9.9

RF II	IU/mL	BS-240E³	35.8	23.2	—	48.4	4.2	BS-240E³	84.1	54.7	—	113.5	9.8
		BS-360E⁴	36.3	23.7	—	48.9	4.2	BS-360E⁴	84.1	54.7	—	113.5	9.8
		BS-380⁵	37.0	24.1	—	49.9	4.3	BS-380⁵	85.1	55.1	—	115.1	10.0
		BS-400⁶	37.0	24.1	—	49.9	4.3	BS-400⁶	85.1	55.1	—	115.1	10.0
		BS-430⁷	36.3	23.7	—	48.9	4.2	BS-430⁷	84.1	54.7	—	113.5	9.8
		BS-480⁸	36.3	23.7	—	48.9	4.2	BS-480⁸	84.1	54.7	—	113.5	9.8
		BS-600⁹	35.1	22.8	—	47.4	4.1	BS-600⁹	79.8	51.9	—	107.7	9.3
		BS-600M¹⁰	36.5	23.6	—	49.4	4.3	BS-600M¹⁰	83.6	54.2	—	113.0	9.8
		BS-800¹¹	36.3	23.7	—	48.9	4.2	BS-800¹¹	84.1	54.7	—	113.5	9.8
		BS-2000¹²	34.2	22.2	—	46.2	4.0	BS-2000¹²	79.1	51.2	—	107.0	9.3
		Abbreviated name	Unit	Level: L/N					Level: H/P				
				Model	Assay Value	Range (Assay Value \pm 3SD)		1 SD	Model	Assay Value	Range (Assay Value \pm 3SD)		1 SD
RF III	IU/mL	BS-800¹¹	44.8	29.2	—	60.4	5.2	BS-800¹¹	103	67	—	139	12
		BS-2000¹²	44.4	28.8	—	60.0	5.2	BS-2000¹²	107	68	—	146	13
		BS-2800M¹³	43.0	28.0	—	58.0	5.0	BS-2800M¹³	101	65	—	137	12

mindray

ASO/CRP/RF Triple Control

For use on: BS-2000, BS-2800M,
BS-600M, BS-800

LOT 059224001

2025-03-24

