

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.BS-200E: BS-200E, BS-220E,;

8.BS-480: BS-480, BS-490;

2.BS-240: BS-230, BS-240;

9.BS-600: BS-600, BS-620;

3.BS-240E: BS240E, BS240Pro;

10.BS-600M: BS-600M, BS-620M;

4.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-

11.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with "W9-" and software version starts with "35.00");

12.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

5.BS-380: BS-380, BS-390;

13.BS-2800M:BS-2600M, BS-2800M;

6.BS-400: BS-400, BS-420;

14.For applicable analyzers of the analyte, please subject to the parameter sheet and instrument.

7.BS-430: BS-430, BS-450,BS-460;

LOT : 059224002

: 2025-07-25

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kisaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)
ASO			CRP		
English		Antistreptolysin "O"			C- Reactive protein
Русский		антистрептолизина O			C-реактивный белок
Português		Antiestreptolisina "O"			proteína C-reactiva
Español		antiestreptolisina "O"			proteína C reactiva
Italiano		Anti-Streptolisina "O"			proteina C-reattiva
Türkçe		Antistreptolisin "O"			C-Reaktif proteini
RF					
English		Rheumatoid Factor			
Русский		ревматоидные факторы			
Português		Fator Reumatoide			
Español		factor reumatoide			
Italiano		fattore reumatoide			
Türkçe		Romatoid Faktör			

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
ASO II	IU/mL	BS-200E ¹	124	79	—	169	15	BS-200E ¹	331	214	—	448	39
		BS-240 ²	123	81	—	165	14	BS-240 ²	324	210	—	438	38
		BS-240E ³	123	81	—	165	14	BS-240E ³	329	215	—	443	38
		BS-360E ⁴	123	81	—	165	14	BS-360E ⁴	329	215	—	443	38
		BS-380 ⁵	124	79	—	169	15	BS-380 ⁵	331	214	—	448	39
		BS-400 ⁶	124	79	—	169	15	BS-400 ⁶	331	214	—	448	39
		BS-430 ⁷	123	81	—	165	14	BS-430 ⁷	329	215	—	443	38
		BS-480 ⁸	123	81	—	165	14	BS-480 ⁸	329	215	—	443	38
		BS-600 ⁹	123	81	—	165	14	BS-600 ⁹	329	215	—	443	38
		BS-600M ¹⁰	118	76	—	160	14	BS-600M ¹⁰	319	208	—	430	37
		BS-800 ¹¹	123	81	—	165	14	BS-800 ¹¹	329	215	—	443	38

BS-2000¹²	119	77	—	161	14	BS-2000¹²	319	208	—	430	37
BS-2800M¹³	118	76	—	160	14	BS-2800M¹³	319	208	—	430	37

Abbreviated name	Unit	Level: L/N				Level: H/P			
		Model	Assay Value	Range (Assay Value±3SD)	1 SD	Model	Assay Value	Range (Assay Value±3SD)	1 SD
CRP II	mg/L	BS-240²	16.4	11.6 — 21.2	1.6	BS-240²	41.1	28.8 — 53.4	4.1
		BS-240E³	17.7	12.3 — 23.1	1.8	BS-240E³	43.4	30.5 — 56.3	4.3
		BS-360E⁴	17.5	12.1 — 22.9	1.8	BS-360E⁴	42.8	29.9 — 55.7	4.3
		BS-380⁵	17.3	12.2 — 22.4	1.7	BS-380⁵	42.9	30.0 — 55.8	4.3
		BS-400⁶	16.3	11.5 — 21.1	1.6	BS-400⁶	41.8	29.2 — 54.4	4.2
		BS-430⁷	16.2	11.4 — 21.0	1.6	BS-430⁷	41.8	29.2 — 54.4	4.2
		BS-480⁸	16.8	11.7 — 21.9	1.7	BS-480⁸	42.1	29.5 — 54.7	4.2
		BS-600⁹	18.2	12.8 — 23.6	1.8	BS-600⁹	43.8	30.6 — 57.0	4.4
		BS-600M¹⁰	17.5	12.1 — 22.9	1.8	BS-600M¹⁰	42.8	29.9 — 55.7	4.3
		BS-800¹¹	17.5	12.1 — 22.9	1.8	BS-800¹¹	42.8	29.9 — 55.7	4.3
	nmol/L	BS-2000¹²	17.7	12.3 — 23.1	1.8	BS-2000¹²	43.2	30.3 — 56.1	4.3
		BS-240²	156	110 — 202	15	BS-240²	391	274 — 508	39
		BS-240E³	169	117 — 220	17	BS-240E³	413	290 — 536	41
		BS-360E⁴	167	115 — 218	17	BS-360E⁴	407	285 — 530	41
		BS-380⁵	165	116 — 213	16	BS-380⁵	408	286 — 531	41
		BS-400⁶	155	109 — 201	15	BS-400⁶	398	278 — 518	40
		BS-430⁷	154	109 — 200	15	BS-430⁷	398	278 — 518	40
		BS-480⁸	160	111 — 208	16	BS-480⁸	401	281 — 521	40
		BS-600⁹	173	122 — 225	17	BS-600⁹	417	291 — 543	42
		BS-600M¹⁰	167	115 — 218	17	BS-600M¹⁰	407	285 — 530	41
BS-800¹¹	167	115 — 218	17	BS-800¹¹	407	285 — 530	41		
BS-2000¹²	169	117 — 220	17	BS-2000¹²	411	288 — 534	41		

Abbreviated name	Unit	Level: L/N				Level: H/P			
		Model	Assay Value	Range (Assay Value±3SD)	1 SD	Model	Assay Value	Range (Assay Value±3SD)	1 SD
RF II	IU/mL	BS-200E¹	33.7	22.0 — 45.4	3.9	BS-200E¹	84.4	54.7 — 114.1	9.9
		BS-240²	32.5	21.1 — 43.9	3.8	BS-240²	80.0	51.8 — 108.2	9.4
		BS-240E³	32.7	21.3 — 44.1	3.8	BS-240E³	81.0	52.5 — 109.5	9.5
		BS-360E⁴	33.5	21.8 — 45.2	3.9	BS-360E⁴	81.9	53.1 — 110.7	9.6
		BS-380⁵	33.1	21.4 — 44.8	3.9	BS-380⁵	81.9	53.1 — 110.7	9.6
		BS-400⁶	33.1	21.4 — 44.8	3.9	BS-400⁶	81.9	53.1 — 110.7	9.6
		BS-430⁷	33.5	21.8 — 45.2	3.9	BS-430⁷	81.9	53.1 — 110.7	9.6
		BS-480⁸	33.5	21.8 — 45.2	3.9	BS-480⁸	81.9	53.1 — 110.7	9.6
		BS-600⁹	33.3	21.6 — 45.0	3.9	BS-600⁹	81.9	53.1 — 110.7	9.6
		BS-600M¹⁰	34.1	22.1 — 46.1	4.0	BS-600M¹⁰	80.7	52.5 — 108.9	9.4
		BS-800¹¹	33.5	21.8 — 45.2	3.9	BS-800¹¹	81.9	53.1 — 110.7	9.6
		BS-2000¹²	32.4	21.0 — 43.8	3.8	BS-2000¹²	78.8	51.2 — 106.4	9.2

Abbreviated name	Unit	Level: L/N				Level: H/P			
		Model	Assay Value	Range (Assay Value±3SD)	1 SD	Model	Assay Value	Range (Assay Value±3SD)	1 SD
RF III	IU/mL	BS-800¹¹	44.0	28.7 — 59.3	5.1	BS-800¹¹	112	73 — 151	13
		BS-2000¹²	43.1	28.1 — 58.1	5.0	BS-2000¹²	107	68 — 146	13
		BS-2800M¹³	42.9	27.9 — 57.9	5.0	BS-2800M¹³	109	70 — 148	13

mindray

ASO/CRP/RF Triple Control

For use on: BS-2000, BS-2800M,
BS-600M, BS-800

LOT 059224002

🕒 2025-07-25

