

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

1. **BS-120**: BS-120, BS-130;

2. **BS-200**: BS-200, BS-220;

3. **BS-200E**: BS-200E, BS-220E;

4. **BS-230**: BS-230, BS-240;

5. **BS-240E**: BS240E, BS240Pro;

6. **BS-330**: BS-330, BS-350;

7. **BS-330E**: **BS-330E**(Serial Number starts with "XQ-"),

BS-350E(Serial Number starts with "XS-");

LOT: **060522001**

🕒: **2023-05-23**

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

8. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S, **BS-330E(V35.00)**(Serial Number starts with "W8-" and software version starts with "35.00"), **BS-350E(V35.00)** (Serial Number starts with "W9-" and software version starts with "35.00");

9. **BS-380**: BS-380, BS-390;

10. **BS-400**: BS-400, BS-420;

11. **BS-430**: BS-430, BS-450, BS-460;

12. **BS-480**: BS-480, BS-490;

13. **BS-600**: BS-600, BS-620;

14. **BS-800**: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

15. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

16. **BS-2800M**: BS-2800M, BS-2800M

Abbreviated name	TPUC Control (L)					TPUC Control (H)								
	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD		
TPUC	BS-120 ¹	mg/L	132	112	—	152	7	BS-120 ¹	mg/L	483	411	—	555	24
	BS-200 ²	mg/L	133	113	—	153	7	BS-200 ²	mg/L	480	408	—	552	24
	BS-200E ³	mg/L	131	111	—	151	7	BS-200E ³	mg/L	495	421	—	569	25
	BS-230 ⁴	mg/L	133	113	—	153	7	BS-230 ⁴	mg/L	480	408	—	552	24
	BS-240E ⁵	mg/L	131	111	—	151	7	BS-240E ⁵	mg/L	490	417	—	564	25
	BS-330 ⁶	mg/L	133	113	—	153	7	BS-330 ⁶	mg/L	480	408	—	552	24
	BS-330E ⁷	mg/L	131	111	—	151	7	BS-330E ⁷	mg/L	495	421	—	569	25
	BS-360E ⁸	mg/L	131	111	—	151	7	BS-360E ⁸	mg/L	490	417	—	564	25
	BS-380 ⁹	mg/L	128	109	—	147	6	BS-380 ⁹	mg/L	485	412	—	558	24
	BS-400 ¹⁰	mg/L	130	111	—	150	7	BS-400 ¹⁰	mg/L	492	418	—	566	25
	BS-430 ¹¹	mg/L	131	111	—	151	7	BS-430 ¹¹	mg/L	490	417	—	564	25
	BS-480 ¹²	mg/L	131	111	—	151	7	BS-480 ¹²	mg/L	490	417	—	564	25
	BS-600 ¹³	mg/L	131	111	—	151	7	BS-600 ¹³	mg/L	492	418	—	566	25
	BS-800 ¹⁴	mg/L	131	111	—	151	7	BS-800 ¹⁴	mg/L	490	417	—	564	25
	BS-2000 ¹⁵	mg/L	136	116	—	156	7	BS-2000 ¹⁵	mg/L	507	431	—	583	25
	BS-2800M ¹⁶	mg/L	134	114	—	154	7	BS-2800M ¹⁶	mg/L	493	419	—	567	25

English	Abbreviated name	TPUC Control	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)
Русский	сокращенное наименование	Контроль общего белка в моче/СМЖ	модель	Прибор	Результат анализа	Диапазон(результат анализа \pm 3CO)
Português	Nome abreviado	Controle PTUL	Modelo	Unidade	Valores da análise	Faixa(Valores da análise \pm 3SD)
Español	nombre abreviado	Control TPUC	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo \pm 3SD)
Italiano	abbreviazione	Controllo TPUC	modelli	Unità	Valori di dosaggio	Intervallo(valore di concentrazione \pm 3 DS)
Türkçe	kısaltılmış ad	TPUC Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri \pm 3SD)

TPUC

English	Total Protein in Urine/CSF (TPUC)
Русский	Общий белок в моче/СМЖ (TPUC)
Português	Proteína total em urina e LCR (PTUL)
Español	proteína total en orina/LCR (TPUC)
Italiano	Proteine Totali Urinarie/CSF (TPUC)
Türkçe	BOS/İdrar Protein (TPUC)