

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

1. **BS-120**: BS-120, BS-130;

2. **BS-200**: BS-200, BS-220;

3. **BS-200E**: BS-200E, BS-220E;

4. **BS-230**: BS-230, BS-240;

5. **BS-240E**: BS240E, BS240Pro;

6. **BS-330**: BS-330, BS-350;

7. **BS-330E**: **BS-330E**(Serial Number starts with "XQ-"),

BS-350E(Serial Number starts with "XS-");

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

8. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S, **BS-330E(V35.00)**(Serial Number starts with "W8-" and software version starts with "35.00"), **BS-350E(V35.00)** (Serial Number starts with "W9-" and software version starts with "35.00");

9. **BS-380**: BS-380, BS-390;

10. **BS-400**: BS-400, BS-420;

11. **BS-430**: BS-430, BS-450, BS-460;

12. **BS-480**: BS-480, BS-490;

13. **BS-600**: BS-600, BS-620;

14. **BS-800**: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

15. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

16. **BS-2800M**: BS-2600M, BS-2800M.

LOT: 060522002

: 2023-09-12

Abbreviated name	TPUC Control (L)					TPUC Control (H)								
	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD		
TPUC	BS-120 ¹	mg/L	132	112	—	152	7	BS-120 ¹	mg/L	490	417	—	564	25
	BS-200 ²	mg/L	133	113	—	153	7	BS-200 ²	mg/L	488	415	—	561	24
	BS-200E ³	mg/L	131	111	—	151	7	BS-200E ³	mg/L	497	422	—	572	25
	BS-230 ⁴	mg/L	134	114	—	154	7	BS-230 ⁴	mg/L	485	412	—	558	24
	BS-240E ⁵	mg/L	135	115	—	155	7	BS-240E ⁵	mg/L	502	427	—	577	25
	BS-330 ⁶	mg/L	133	113	—	153	7	BS-330 ⁶	mg/L	488	415	—	561	24
	BS-330E ⁷	mg/L	131	111	—	151	7	BS-330E ⁷	mg/L	497	422	—	572	25
	BS-360E ⁸	mg/L	135	115	—	155	7	BS-360E ⁸	mg/L	502	427	—	577	25
	BS-380 ⁹	mg/L	129	110	—	148	6	BS-380 ⁹	mg/L	496	422	—	570	25
	BS-400 ¹⁰	mg/L	129	110	—	148	6	BS-400 ¹⁰	mg/L	498	423	—	573	25
	BS-430 ¹¹	mg/L	135	115	—	155	7	BS-430 ¹¹	mg/L	502	427	—	577	25
	BS-480 ¹²	mg/L	135	115	—	155	7	BS-480 ¹²	mg/L	502	427	—	577	25
	BS-600 ¹³	mg/L	133	113	—	153	7	BS-600 ¹³	mg/L	500	425	—	575	25
	BS-800 ¹⁴	mg/L	135	115	—	155	7	BS-800 ¹⁴	mg/L	502	427	—	577	25
	BS-2000 ¹⁵	mg/L	137	116	—	158	7	BS-2000 ¹⁵	mg/L	514	437	—	591	26
	BS-2800M ¹⁶	mg/L	133	113	—	153	7	BS-2800M ¹⁶	mg/L	495	421	—	569	25

English	Abbreviated name	TPUC Control	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)
Русский	сокращенное наименование	Контроль общего белка в моче/СМЖ	модель	Прибор	Результат анализа	Диапазон(результат анализа \pm 3CO)
Português	Nome abreviado	Controle PTUL	Modelo	Unidade	Valores da análise	Faixa(Valores da análise \pm 3SD)
Español	nombre abreviado	Control TPUC	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo \pm 3SD)
Italiano	abbreviazione	Controllo TPUC	modelli	Unità	Valori di dosaggio	Intervallo(valore di concentrazione \pm 3 DS)
Türkçe	kısaltılmış ad	TPUC Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri \pm 3SD)

TPUC

English	Total Protein in Urine/CSF (TPUC)
Русский	Общий белок в моче/СМЖ (TPUC)
Português	Proteína total em urina e LCR (PTUL)
Español	proteína total en orina/LCR (TPUC)
Italiano	Proteine Totali Urinarie/CSF (TPUC)
Türkçe	BOS/İdrar Protein (TPUC)