

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

1. **BS-120**: BS-120, BS-130;

2. **BS-200**: BS-200, BS-220;

3. **BS-200E**: BS-200E, BS-220E;

4. **BS-230**: BS-230, BS-240;

5. **BS-240E**: BS240E, BS240Pro;

6. **BS-330**: BS-330, BS-350;

7. **BS-330E**: **BS-330E**(Serial Number starts with "XQ-"),

BS-350E(Serial Number starts with "XS-");

8. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S, **BS-330E(V35.00)**(Serial Number starts with "W8-" and software version starts with "35.00"), **BS-350E(V35.00)** (Serial Number starts with "W9-" and software version starts with "35.00");

9. **BS-380**: BS-380, BS-390;

10. **BS-400**: BS-400, BS-420;

11. **BS-430**: BS-430, BS-450, BS-460;

12. **BS-480**: BS-480, BS-490;

13. **BS-600**: BS-600, BS-620;

14. **BS-600M**: BS-600M, BS-620M;

15. **BS-800**: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

16. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

17. **BS-2800M**: BS-2600M, BS-2800M.

LOT: 060522003

🕒: 2023-09-12

| Abbreviated name | TPUC Control (L) | | | | | TPUC Control (H) | | | | | | | | |
|------------------|-------------------------------|------|-------------|-------------------------------|---|------------------|-------|-------------------------------|-------------|-------------------------------|-----|------|-----|----|
| | Model | Unit | Assay Value | Range (Assay Value \pm 3SD) | | 1 SD | Model | Unit | Assay Value | Range (Assay Value \pm 3SD) | | 1 SD | | |
| TPUC | BS-120 ¹ | mg/L | 132 | 112 | — | 152 | 7 | BS-120 ¹ | mg/L | 490 | 417 | — | 564 | 25 |
| | BS-200 ² | mg/L | 133 | 113 | — | 153 | 7 | BS-200 ² | mg/L | 488 | 415 | — | 561 | 24 |
| | BS-200E ³ | mg/L | 131 | 111 | — | 151 | 7 | BS-200E ³ | mg/L | 497 | 422 | — | 572 | 25 |
| | BS-230 ⁴ | mg/L | 134 | 114 | — | 154 | 7 | BS-230 ⁴ | mg/L | 485 | 412 | — | 558 | 24 |
| | BS-240E ⁵ | mg/L | 135 | 115 | — | 155 | 7 | BS-240E ⁵ | mg/L | 502 | 427 | — | 577 | 25 |
| | BS-330 ⁶ | mg/L | 133 | 113 | — | 153 | 7 | BS-330 ⁶ | mg/L | 488 | 415 | — | 561 | 24 |
| | BS-330E ⁷ | mg/L | 131 | 111 | — | 151 | 7 | BS-330E ⁷ | mg/L | 497 | 422 | — | 572 | 25 |
| | BS-360E ⁸ | mg/L | 135 | 115 | — | 155 | 7 | BS-360E ⁸ | mg/L | 502 | 427 | — | 577 | 25 |
| | BS-380 ⁹ | mg/L | 129 | 110 | — | 148 | 6 | BS-380 ⁹ | mg/L | 496 | 422 | — | 570 | 25 |
| | BS-400 ¹⁰ | mg/L | 129 | 110 | — | 148 | 6 | BS-400 ¹⁰ | mg/L | 498 | 423 | — | 573 | 25 |
| | BS-430 ¹¹ | mg/L | 135 | 115 | — | 155 | 7 | BS-430 ¹¹ | mg/L | 502 | 427 | — | 577 | 25 |
| | BS-480 ¹² | mg/L | 135 | 115 | — | 155 | 7 | BS-480 ¹² | mg/L | 502 | 427 | — | 577 | 25 |
| | BS-600 ¹³ | mg/L | 133 | 113 | — | 153 | 7 | BS-600 ¹³ | mg/L | 500 | 425 | — | 575 | 25 |
| | BS-600M ¹⁴ | mg/L | 133 | 113 | — | 153 | 7 | BS-600M ¹⁴ | mg/L | 495 | 421 | — | 569 | 25 |
| | BS-800 ¹⁵ | mg/L | 135 | 115 | — | 155 | 7 | BS-800 ¹⁵ | mg/L | 502 | 427 | — | 577 | 25 |
| | BS-2000 ¹⁶ | mg/L | 137 | 116 | — | 158 | 7 | BS-2000 ¹⁶ | mg/L | 514 | 437 | — | 591 | 26 |
| | BS-2800M ¹⁷ | mg/L | 133 | 113 | — | 153 | 7 | BS-2800M ¹⁷ | mg/L | 495 | 421 | — | 569 | 25 |

| English | Abbreviated name | TPUC Control | Model | Unit | Assay Value | Range (Assay Value \pm 3SD) |
|-----------|--------------------------|----------------------------------|---------|---------|--------------------|---|
| Русский | сокращенное наименование | Контроль общего белка в моче/СМЖ | модель | Прибор | Результат анализа | Диапазон(результат анализа \pm 3CO) |
| Português | Nome abreviado | Controle PTUL | Modelo | Unidade | Valores da análise | Faixa(Valores da análise \pm 3SD) |
| Español | nombre abreviado | Control TPUC | modelo | Unidad | Valor de ensayo | Rango(Valor de ensayo \pm 3SD) |
| Italiano | abbreviazione | Controllo TPUC | modelli | Unità | Valori di dosaggio | Intervallo(valore di concentrazione \pm 3 DS) |
| Türkçe | kısaltılmış ad | TPUC Kontrolü | model | Ünite | Tayin Değeri | Aralık (Tayin Değeri \pm 3SD) |

TPUC

| | |
|------------------|--------------------------------------|
| English | Total Protein in Urine/CSF (TPUC) |
| Русский | Общий белок в моче/СМЖ (TPUC) |
| Português | Proteína total em urina e LCR (PTUL) |
| Español | proteína total en orina/LCR (TPUC) |
| Italiano | Proteine Totali Urinarie/CSF (TPUC) |
| Türkçe | BOS/İdrar Protein (TPUC) |