

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

1.**BS-120**: BS-120, BS-130;

2.**BS-200**: BS-200, BS-220;

3.**BS-200E**: BS-200E, BS-220E;

4.**BS-230**: BS-230, BS-240;

5.**BS-240E**: BS240E, BS240Pro;

6.**BS-330**: BS-330, BS-350;

7.**BS-330E**: **BS-330E**(Serial Number starts with "XQ-"),

BS-350E(Serial Number starts with "XS-");

8.**BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,
BS-330E(V35.00)(Serial Number starts with "W8-"
and software version starts with "35.00"),
BS-350E(V35.00) (Serial Number starts with "W9-"
and software version starts with "35.00");

9.**BS-380**: BS-380, BS-390;

10.**BS-400**: BS-400, BS-420;

11.**BS-430**: BS-430, BS-450, BS-460;

12.**BS-480**: BS-480, BS-490;

13.**BS-600**: BS-600, BS-620;

14.**BS-600M**: BS-600M, BS-620M;

15.**BS-800**: BS-800, BS-820, BS-800M, BS-820M,

BS-1800, BS-1800plus;

16.**BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

17.**BS-2800M**:BS-2600M, BS-2800M.

LOT: 060523004

: 2024-08-23

Abbreviated name	TPUC Control (L)					TPUC Control (H)								
	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD		
TPUC	BS-120 ¹	mg/L	123	105	—	141	6	BS-120 ¹	mg/L	475	404	—	546	24
	BS-200 ²	mg/L	124	105	—	143	6	BS-200 ²	mg/L	477	405	—	549	24
	BS-200E ³	mg/L	122	104	—	140	6	BS-200E ³	mg/L	480	408	—	552	24
	BS-230 ⁴	mg/L	126	107	—	145	6	BS-230 ⁴	mg/L	475	404	—	546	24
	BS-240E ⁵	mg/L	123	105	—	141	6	BS-240E ⁵	mg/L	487	414	—	560	24
	BS-330 ⁶	mg/L	124	105	—	143	6	BS-330 ⁶	mg/L	477	405	—	549	24
	BS-330E ⁷	mg/L	122	104	—	140	6	BS-330E ⁷	mg/L	480	408	—	552	24
	BS-360E ⁸	mg/L	123	105	—	141	6	BS-360E ⁸	mg/L	487	414	—	560	24
	BS-380 ⁹	mg/L	121	103	—	139	6	BS-380 ⁹	mg/L	484	411	—	557	24
	BS-400 ¹⁰	mg/L	122	104	—	140	6	BS-400 ¹⁰	mg/L	489	416	—	562	24
	BS-430 ¹¹	mg/L	123	105	—	141	6	BS-430 ¹¹	mg/L	487	414	—	560	24
	BS-480 ¹²	mg/L	123	105	—	141	6	BS-480 ¹²	mg/L	487	414	—	560	24
	BS-600 ¹³	mg/L	124	105	—	143	6	BS-600 ¹³	mg/L	489	416	—	562	24
	BS-600M ¹⁴	mg/L	126	107	—	145	6	BS-600M ¹⁴	mg/L	489	416	—	562	24
	BS-800 ¹⁵	mg/L	123	105	—	141	6	BS-800 ¹⁵	mg/L	487	414	—	560	24
	BS-2000 ¹⁶	mg/L	126	107	—	145	6	BS-2000 ¹⁶	mg/L	502	427	—	577	25
	BS-2800M ¹⁷	mg/L	126	107	—	145	6	BS-2800M ¹⁷	mg/L	489	416	—	562	24

English	Abbreviated name	TPUC Control	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)
Русский	сокращенное наименование	Контроль общего белка в моче/СМЖ	модель	Прибор	Результат анализа	Диапазон(результат анализа \pm 3CO)
Português	Nome abreviado	Controle PTUL	Modelo	Unidade	Valores da análise	Faixa(Valores da análise \pm 3SD)
Español	nombre abreviado	Control TPUC	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo \pm 3SD)
Italiano	abbreviazione	Controllo TPUC	modelli	Unità	Valori di dosaggio	Intervallo(valore di concentrazione \pm 3 DS)
Türkçe	kısaltılmış ad	TPUC Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri \pm 3SD)

TPUC

English	Total Protein in Urine/CSF (TPUC)
Русский	Общий белок в моче/СМЖ (TPUC)
Português	Proteína total em urina e LCR (PTUL)
Español	proteína total en orina/LCR (TPUC)
Italiano	Proteine Totali Urinarie/CSF (TPUC)
Türkçe	BOS/İdrar Protein (TPUC)