

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

- 1.BS-120: BS-120, BS-130, BS-180, BS-190;
- 2.BS-200: BS-200, BS-220, BS-330, BS-350;
- 3.BS-200E: BS-200E, BS-220E,;
- 4.BS-240: BS-230, BS-240;
- 5.BS-240E: BS240E, BS240Pro;
- 6.BS-330E:BS-330E(Serial Number starts with "XQ-"), BS-350E(Serial Number starts with "XS-")
- 7.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00");
- 8.BS-380: BS-380, BS-390;

- 9.BS-400: BS-400, BS-420;
- 10.BS-430: BS-430, BS-450,BS-460;
- 11.BS-480: BS-480, BS-490;
- 12.BS-600: BS-600, BS-620;
- 13.BS-600M: BS-600M, BS-620M;
- 14.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;
- 15.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;
- 16.BS-2800M:BS-2600M, BS-2800M;
- 17.Applicable models of the chemistry shall be subject to the parameter list and instructions.

LOT : 060523005

🕒: 2025-01-02

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

TPUC

English	Total Protein in Urine/CSF (TPUC)
Русский	Общий белок в моче/СМЖ (TPUC)
Português	Proteína total em urina e LCR (PTUL)
Español	proteína total en orina/LCR (TPUC)
Italiano	Proteine Totali Urinarie/CSF (TPUC)
Türkçe	BOS/İdrar Protein (TPUC)

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
TPUC	mg/L	BS-120 ¹	126	108	—	144	6	BS-120 ¹	461	392	—	530	23
		BS-200 ²	126	108	—	144	6	BS-200 ²	469	400	—	538	23
		BS-200E ³	122	104	—	140	6	BS-200E ³	461	392	—	530	23
		BS-240 ⁴	125	107	—	143	6	BS-240 ⁴	456	387	—	525	23
		BS-240E ⁵	126	108	—	144	6	BS-240E ⁵	470	398	—	542	24
		BS-330E ⁶	122	104	—	140	6	BS-330E ⁶	461	392	—	530	23
		BS-360E ⁷	126	108	—	144	6	BS-360E ⁷	470	398	—	542	24
		BS-380 ⁸	127	109	—	145	6	BS-380 ⁸	480	408	—	552	24
		BS-400 ⁹	126	108	—	144	6	BS-400 ⁹	480	408	—	552	24
		BS-430 ¹⁰	126	108	—	144	6	BS-430 ¹⁰	470	398	—	542	24
		BS-480 ¹¹	126	108	—	144	6	BS-480 ¹¹	470	398	—	542	24
		BS-600 ¹²	129	111	—	147	6	BS-600 ¹²	486	414	—	558	24
		BS-600M ¹³	132	111	—	153	7	BS-600M ¹³	488	416	—	560	24
		BS-800 ¹⁴	126	108	—	144	6	BS-800 ¹⁴	470	398	—	542	24
		BS-2000 ¹⁵	132	111	—	153	7	BS-2000 ¹⁵	493	418	—	568	25
		BS-2800M ¹⁶	132	111	—	153	7	BS-2800M ¹⁶	488	416	—	560	24

mindray
TPUC Control

For use on: BS-2000, BS-2800M,
BS-600M, BS-800

LOT 060523005

2025-01-02

