

### Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

- 1.BS-120: BS-120, BS-130, BS-180, BS-190;
- 2.BS-200: BS-200, BS-220, BS-330, BS-350;
- 3.BS-200E: BS-200E, BS-220E,;
- 4.BS-240: BS-230, BS-240;
- 5.BS-240E: BS240E, BS240Pro;
- 6.BS-330E:BS-330E(Serial Number starts with "XQ-"), BS-350E(Serial Number starts with "XS-")
- 7.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00");
- 8.BS-380: BS-380, BS-390;

- 9.BS-400: BS-400, BS-420;
- 10.BS-430: BS-430, BS-450,BS-460;
- 11.BS-480: BS-480, BS-490;
- 12.BS-600: BS-600, BS-620;
- 13.BS-600M: BS-600M, BS-620M;
- 14.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;
- 15.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;
- 16.BS-2800M:BS-2600M, BS-2800M;
- 17.Applicable models of the chemistry shall be subject to the parameter list and instructions.

LOT : 060523006

🕒: 2025-01-02

| English   | Abbreviated name         | Model   | Unit    | Assay Value        | Range (Assay Value±3SD)                    |
|-----------|--------------------------|---------|---------|--------------------|--|
| Русский   | сокращенное наименование | модель  | Прибор  | Результат анализа  | Диапазон(результат анализа ± 3CO)          |
| Português | Nome abreviado           | Modelo  | Unidade | Valores da análise | Faixa(Valores da análise ±3SD)             |
| Español   | nombre abreviado         | modelo  | Unidad  | Valor de ensayo    | Rango(Valor de ensayo ±3SD)                |
| Italiano  | abbreviazione            | modelli | Unità   | Valori di dosaggio | Intervallo (valore di concentrazione±3 DS) |
| Türkçe    | kısaltılmış ad           | model   | Ünite   | Tayin Değeri       | Aralık (Tayin Değeri±3SD)                  |

**TPUC**

|           |                                      |
|-----------|--------------------------------------|
| English   | Total Protein in Urine/CSF (TPUC)    |
| Русский   | Общий белок в моче/СМЖ (TPUC)        |
| Português | Proteína total em urina e LCR (PTUL) |
| Español   | proteína total en orina/LCR (TPUC)   |
| Italiano  | Proteine Totali Urinarie/CSF (TPUC)  |
| Türkçe    | BOS/İdrar Protein (TPUC)             |

| Abbreviated name | Unit | Level: L/N             |             |                         |   |      | Level: H/P |                        |                         |     |      |     |    |
|------------------|------|------------------------|-------------|-------------------------|---|------|------------|------------------------|-------------------------|-----|------|-----|----|
|                  |      | Model                  | Assay Value | Range (Assay Value±3SD) |   | 1 SD | Model      | Assay Value            | Range (Assay Value±3SD) |     | 1 SD |     |    |
| TPUC             | mg/L | BS-120 <sup>1</sup>    | 126         | 108                     | — | 144  | 6          | BS-120 <sup>1</sup>    | 461                     | 392 | —    | 530 | 23 |
|                  |      | BS-200 <sup>2</sup>    | 126         | 108                     | — | 144  | 6          | BS-200 <sup>2</sup>    | 469                     | 400 | —    | 538 | 23 |
|                  |      | BS-200E <sup>3</sup>   | 122         | 104                     | — | 140  | 6          | BS-200E <sup>3</sup>   | 461                     | 392 | —    | 530 | 23 |
|                  |      | BS-240 <sup>4</sup>    | 125         | 107                     | — | 143  | 6          | BS-240 <sup>4</sup>    | 456                     | 387 | —    | 525 | 23 |
|                  |      | BS-240E <sup>5</sup>   | 126         | 108                     | — | 144  | 6          | BS-240E <sup>5</sup>   | 470                     | 398 | —    | 542 | 24 |
|                  |      | BS-330E <sup>6</sup>   | 122         | 104                     | — | 140  | 6          | BS-330E <sup>6</sup>   | 461                     | 392 | —    | 530 | 23 |
|                  |      | BS-360E <sup>7</sup>   | 126         | 108                     | — | 144  | 6          | BS-360E <sup>7</sup>   | 470                     | 398 | —    | 542 | 24 |
|                  |      | BS-380 <sup>8</sup>    | 127         | 109                     | — | 145  | 6          | BS-380 <sup>8</sup>    | 480                     | 408 | —    | 552 | 24 |
|                  |      | BS-400 <sup>9</sup>    | 126         | 108                     | — | 144  | 6          | BS-400 <sup>9</sup>    | 480                     | 408 | —    | 552 | 24 |
|                  |      | BS-430 <sup>10</sup>   | 126         | 108                     | — | 144  | 6          | BS-430 <sup>10</sup>   | 470                     | 398 | —    | 542 | 24 |
|                  |      | BS-480 <sup>11</sup>   | 126         | 108                     | — | 144  | 6          | BS-480 <sup>11</sup>   | 470                     | 398 | —    | 542 | 24 |
|                  |      | BS-600 <sup>12</sup>   | 129         | 111                     | — | 147  | 6          | BS-600 <sup>12</sup>   | 486                     | 414 | —    | 558 | 24 |
|                  |      | BS-600M <sup>13</sup>  | 132         | 111                     | — | 153  | 7          | BS-600M <sup>13</sup>  | 488                     | 416 | —    | 560 | 24 |
|                  |      | BS-800 <sup>14</sup>   | 126         | 108                     | — | 144  | 6          | BS-800 <sup>14</sup>   | 470                     | 398 | —    | 542 | 24 |
|                  |      | BS-2000 <sup>15</sup>  | 132         | 111                     | — | 153  | 7          | BS-2000 <sup>15</sup>  | 493                     | 418 | —    | 568 | 25 |
|                  |      | BS-2800M <sup>16</sup> | 132         | 111                     | — | 153  | 7          | BS-2800M <sup>16</sup> | 488                     | 416 | —    | 560 | 24 |

**mindray**  
TPUC Control

For use on: BS-2000, BS-2800M,  
BS-600M, BS-800

LOT 060523006

2025-01-02

