

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

- 1.BS-120: BS-120, BS-130, BS-180, BS-190;
- 2.BS-200: BS-200, BS-220, BS-330, BS-350;
- 3.BS-240: BS-230, BS-240;
- 4.BS-240E: BS240E, BS240Pro;
- 5.BS-330E:BS-330E(Serial Number starts with "XQ-"), BS-350E(Serial Number starts with "XS-")
- 6.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00");
- 7.BS-380: BS-380, BS-390;
- 8.BS-400: BS-400, BS-420;

- 9.BS-430: BS-430, BS-450,BS-460;
- 10.BS-480: BS-480, BS-490;
- 11.BS-600: BS-600, BS-620;
- 12.BS-600M: BS-600M, BS-620M;
- 13.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;
- 14.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;
- 15.BS-2800M:BS-2600M, BS-2800M;
- 16.Applicable models of the chemistry shall be subject to the parameter list and instructions.

LOT : 060823008

🕒: 2024-11-29

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

HCY

English	Homocysteine (HCY)
Русский	Гомоцистеин (HCY)
Português	Homocisteína (HCY)
Español	Homocisteína (HCY)
Italiano	Omocisteina (HCY)
Türkçe	Homocistain (HCY)

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
HCY II	µmol/L	BS-120 ¹	14.4	11.4	—	17.4	1.0	BS-120 ¹	30.9	24.6	—	37.2	2.1
		BS-200 ²	15.1	12.1	—	18.1	1.0	BS-200 ²	31.2	24.9	—	37.5	2.1
		BS-240 ³	15.1	12.1	—	18.1	1.0	BS-240 ³	30.5	24.5	—	36.5	2.0
		BS-240E ⁴	15.4	12.4	—	18.4	1.0	BS-240E ⁴	31.0	24.7	—	37.3	2.1
		BS-330E ⁵	14.8	11.8	—	17.8	1.0	BS-330E ⁵	31.2	24.9	—	37.5	2.1
		BS-360E ⁶	15.4	12.4	—	18.4	1.0	BS-360E ⁶	31.0	24.7	—	37.3	2.1
		BS-380 ⁷	15.4	12.4	—	18.4	1.0	BS-380 ⁷	31.2	24.9	—	37.5	2.1
		BS-400 ⁸	15.4	12.4	—	18.4	1.0	BS-400 ⁸	31.2	24.9	—	37.5	2.1
		BS-430 ⁹	15.4	12.4	—	18.4	1.0	BS-430 ⁹	31.0	24.7	—	37.3	2.1
		BS-480 ¹⁰	15.2	12.2	—	18.2	1.0	BS-480 ¹⁰	30.8	24.5	—	37.1	2.1
		BS-600 ¹¹	15.4	12.4	—	18.4	1.0	BS-600 ¹¹	31.0	24.7	—	37.3	2.1
		BS-600M ¹²	15.5	12.5	—	18.5	1.0	BS-600M ¹²	30.9	24.6	—	37.2	2.1
		BS-800 ¹³	15.4	12.4	—	18.4	1.0	BS-800 ¹³	31.0	24.7	—	37.3	2.1
		BS-2000 ¹⁴	15.7	12.7	—	18.7	1.0	BS-2000 ¹⁴	32.5	25.9	—	39.1	2.2
		BS-2800M ¹⁵	15.7	12.7	—	18.7	1.0	BS-2800M ¹⁵	31.4	25.1	—	37.7	2.1

mindray
HCY Control

For use on: BS-2000, BS-2800M,
BS-600M, BS-800

LOT 060823008

2024-11-29

