

## Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.BS-120: BS-120, BS-130, BS-180, BS-190;

9.BS-430: BS-430, BS-450,BS-460;

2.BS-200: BS-200, BS-220, BS-330, BS-350;

10.BS-480: BS-480, BS-490;

3.BS-240: BS-230, BS-240;

11.BS-600: BS-600, BS-620;

4.BS-240E: BS240E, BS240Pro;

12.BS-600M: BS-600M, BS-620M;

5.BS-330E:BS-330E(Serial Number starts with "XQ-"), BS-350E(Serial Number starts with "XS-")

13.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

6.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-

14.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00");

15.BS-2800M:BS-2600M, BS-2800M;

7.BS-380: BS-380, BS-390;

16.For applicable analyzers of the analyte, please subject to the parameter sheet and instrument.

8.BS-400: BS-400, BS-420;

LOT : 060824003

🕒: 2025-05-24

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

### HCY

English	Homocysteine (HCY)
Русский	Гомоцистеин (HCY)
Português	Homocisteína (HCY)
Español	Homocisteína (HCY)
Italiano	Omocisteina (HCY)
Türkçe	Homosistain (HCY)

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
HCY II	µmol/L	BS-120 <sup>1</sup>	17.1	13.8	—	20.4	1.1	BS-120 <sup>1</sup>	35.0	28.1	—	41.9	2.3
		BS-200 <sup>2</sup>	16.9	13.6	—	20.2	1.1	BS-200 <sup>2</sup>	34.8	27.9	—	41.7	2.3
		BS-240 <sup>3</sup>	16.5	13.2	—	19.8	1.1	BS-240 <sup>3</sup>	34.3	27.4	—	41.2	2.3
		BS-240E <sup>4</sup>	16.8	13.5	—	20.1	1.1	BS-240E <sup>4</sup>	35.0	28.1	—	41.9	2.3
		BS-330E <sup>5</sup>	16.5	13.2	—	19.8	1.1	BS-330E <sup>5</sup>	34.6	27.7	—	41.5	2.3
		BS-360E <sup>6</sup>	16.8	13.5	—	20.1	1.1	BS-360E <sup>6</sup>	35.0	28.1	—	41.9	2.3
		BS-380 <sup>7</sup>	16.2	12.9	—	19.5	1.1	BS-380 <sup>7</sup>	34.1	27.2	—	41.0	2.3
		BS-400 <sup>8</sup>	16.2	12.9	—	19.5	1.1	BS-400 <sup>8</sup>	34.1	27.2	—	41.0	2.3
		BS-430 <sup>9</sup>	16.8	13.5	—	20.1	1.1	BS-430 <sup>9</sup>	35.0	28.1	—	41.9	2.3
		BS-480 <sup>10</sup>	16.8	13.5	—	20.1	1.1	BS-480 <sup>10</sup>	34.6	27.7	—	41.5	2.3
		BS-600 <sup>11</sup>	16.8	13.5	—	20.1	1.1	BS-600 <sup>11</sup>	35.0	28.1	—	41.9	2.3
		BS-600M <sup>12</sup>	16.8	13.5	—	20.1	1.1	BS-600M <sup>12</sup>	34.7	27.8	—	41.6	2.3
		BS-800 <sup>13</sup>	16.8	13.5	—	20.1	1.1	BS-800 <sup>13</sup>	35.0	28.1	—	41.9	2.3
		BS-2000 <sup>14</sup>	16.8	13.5	—	20.1	1.1	BS-2000 <sup>14</sup>	34.6	27.7	—	41.5	2.3
		BS-2800M <sup>15</sup>	16.8	13.5	—	20.1	1.1	BS-2800M <sup>15</sup>	34.4	27.5	—	41.3	2.3

**mindray**  
HCY Control

For use on: BS-2000, BS-2800M,  
BS-600M, BS-800

LOT 060824003

2025-05-24

