

# Lp(a) II Control

## Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

1.BS-120: BS-120, BS-130;

2.BS-180: BS-180, BS-190;

3.BS-200: BS-200, BS-220;

4.BS-200E: BS-200E, BS-220E;

5.BS-230: BS-230, BS-240;

6.BS-240E: BS240E, BS240Pro;

7.BS-300: BS-300, BS-320;

8.BS-330: BS-330, BS-350;

9.BS-330E:BS-330E(Serial Number starts with "XQ-"),

BS-350E(Serial Number starts with "XS-");

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

10.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S,  
BS-330E(V35.00)(Serial Number starts with "W8-"  
and software version starts with "35.00"),  
BS-350E(V35.00) (Serial Number starts with"W9-"  
and software version starts with "35.00");

11.BS-380: BS-380, BS-390;

12.BS-400: BS-400, BS-420;

13.BS-430: BS-430, BS-450, BS-460;

14.BS-480: BS-480, BS-490;

15.BS-600: BS-600, BS-620;

16.BS-800: BS-800, BS-820, BS-800M, BS-820M,  
BS-1800, BS-1800plus;

17.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M.

LOT: 061121005

: 2023-03-06

Abbreviated name	Lp(a) II Control (L)					Lp(a) II Control (H)								
	Model	Unit	Assay Value	Range (Assay Value $\pm$ 3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value $\pm$ 3SD)		1 SD		
Lp(a) II	BS-120 <sup>1</sup>	nmol/L	45.4	30.3	—	60.5	5.0	BS-120 <sup>1</sup>	nmol/L	117	78	—	156	13
	BS-180 <sup>2</sup>	nmol/L	45.4	30.3	—	60.5	5.0	BS-180 <sup>2</sup>	nmol/L	117	78	—	156	13
	BS-200 <sup>3</sup>	nmol/L	47.2	31.5	—	62.9	5.2	BS-200 <sup>3</sup>	nmol/L	117	78	—	156	13
	BS-200E <sup>4</sup>	nmol/L	45.5	30.3	—	60.7	5.1	BS-200E <sup>4</sup>	nmol/L	115	77	—	153	13
	BS-230 <sup>5</sup>	nmol/L	48.7	32.5	—	64.9	5.4	BS-230 <sup>5</sup>	nmol/L	115	77	—	153	13
	BS-240E <sup>6</sup>	nmol/L	45.6	30.4	—	60.8	5.1	BS-240E <sup>6</sup>	nmol/L	115	77	—	153	13
	BS-300 <sup>7</sup>	nmol/L	47.8	31.9	—	63.7	5.3	BS-300 <sup>7</sup>	nmol/L	119	79	—	159	13
	BS-330 <sup>8</sup>	nmol/L	47.2	31.5	—	62.9	5.2	BS-330 <sup>8</sup>	nmol/L	117	78	—	156	13
	BS-330E <sup>9</sup>	nmol/L	45.5	30.3	—	60.7	5.1	BS-330E <sup>9</sup>	nmol/L	115	77	—	153	13
	BS-360E <sup>10</sup>	nmol/L	43.9	29.3	—	58.5	4.9	BS-360E <sup>10</sup>	nmol/L	114	76	—	152	13
	BS-380 <sup>11</sup>	nmol/L	45.0	30.0	—	60.0	5.0	BS-380 <sup>11</sup>	nmol/L	116	77	—	155	13
	BS-400 <sup>12</sup>	nmol/L	46.8	31.2	—	62.4	5.2	BS-400 <sup>12</sup>	nmol/L	119	79	—	159	13
	BS-430 <sup>13</sup>	nmol/L	46.3	30.9	—	61.7	5.1	BS-430 <sup>13</sup>	nmol/L	116	77	—	155	13
	BS-480 <sup>14</sup>	nmol/L	47.5	31.7	—	63.3	5.3	BS-480 <sup>14</sup>	nmol/L	120	80	—	160	13
	BS-600 <sup>15</sup>	nmol/L	49.1	32.7	—	65.5	5.5	BS-600 <sup>15</sup>	nmol/L	121	81	—	161	13
	BS-800 <sup>16</sup>	nmol/L	47.2	31.5	—	62.9	5.2	BS-800 <sup>16</sup>	nmol/L	117	78	—	156	13
	BS-2000 <sup>17</sup>	nmol/L	45.4	30.3	—	60.5	5.0	BS-2000 <sup>17</sup>	nmol/L	115	77	—	153	13

English	Abbreviated name	Lp(a) II Control	Model	Unit	Assay Value	Range (Assay Value $\pm$ 3SD)
Русский	сокращенное наименование	Lp(a) II Контроль	модель	Прибор	Результат анализа	Диапазон(результат анализа $\pm$ 3CO)
Português	Nome abreviado	Controle da Lp(a) II	Modelo	Unidade	Valores da análise	Faixa(Valores da análise $\pm$ 3SD)
Español	nombre abreviado	Control Lp(a) II	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo $\pm$ 3SD)
Italiano	abbreviazione	Controllo Lp(a) II	modelli	Unità	Valori di dosaggio	Intervallo(valore di concentrazione $\pm$ 3SD)
Türkçe	kisaltılmış ad	Lp(a) II Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri $\pm$ 3SD)

Lp(a) II Control	
English	Lipoprotein (a) Control
Русский	Липопротеин (а) Контроль
Português	Controle da Lipoproteína (a)
Español	Control de Lipoproteína (a)
Italiano	Controllo della Lipoproteina (a)
Türkçe	Lipoprotein (a) Kontrolü