

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray
 Português: Planilha de valores do Sistema de Medição BS da Mindray
 Español: Hoja de valores del sistema de medición Mindray BS
 Italiano: Scheda dei valori del sistema di misurazione BS di Mindray
 Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

1.BS-120: BS-120, BS-130;

2.BS-180: BS-180, BS-190;

3.BS-200: BS-200, BS-220;

4.BS-200E: BS-200E, BS-220E;

5.BS-230: BS-230, BS-240;

6.BS-240E: BS240E, BS240Pro;

7.BS-300: BS-300, BS-320;

8.BS-330: BS-330, BS-350;

9.BS-330E:BS-330E(Serial Number starts with "XQ-"),

BS-350E(Serial Number starts with "XS-");

10.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S,
 BS-330E(V35.00)(Serial Number starts with "W8-"
 and software version starts with "35.00"),
 BS-350E(V35.00) (Serial Number starts with"W9-"
 and software version starts with "35.00");

11.BS-380: BS-380, BS-390;

12.BS-400: BS-400, BS-420;

13.BS-430: BS-430, BS-450, BS-460;

14.BS-480: BS-480, BS-490;

15.BS-600: BS-600, BS-620;

16.BS-800: BS-800, BS-820, BS-800M, BS-820M,
 BS-1800, BS-1800plus;

17.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

18.BS-2800M:BS-2600M, BS-2800M.

LOT: 061121006

: 2023-04-09

Abbreviate d name	Lp(a) II Control (L)					Lp(a) II Control (H)				
	Model	Unit	Assay Value	Range (Assay Value±3SD)	1 SD	Model	Unit	Assay Value	Range (Assay Value±3SD)	1 SD
Lp(a) II	BS-120 ¹	nmol/L	45.4	30.3 — 60.5	5.0	BS-120 ¹	nmol/L	117	78 — 156	13
	BS-180 ²	nmol/L	45.4	30.3 — 60.5	5.0	BS-180 ²	nmol/L	117	78 — 156	13
	BS-200 ³	nmol/L	47.2	31.5 — 62.9	5.2	BS-200 ³	nmol/L	117	78 — 156	13
	BS-200E ⁴	nmol/L	45.5	30.3 — 60.7	5.1	BS-200E ⁴	nmol/L	115	77 — 153	13
	BS-230 ⁵	nmol/L	48.7	32.5 — 64.9	5.4	BS-230 ⁵	nmol/L	115	77 — 153	13
	BS-240E ⁶	nmol/L	45.6	30.4 — 60.8	5.1	BS-240E ⁶	nmol/L	115	77 — 153	13
	BS-300 ⁷	nmol/L	47.8	31.9 — 63.7	5.3	BS-300 ⁷	nmol/L	119	79 — 159	13
	BS-330 ⁸	nmol/L	47.2	31.5 — 62.9	5.2	BS-330 ⁸	nmol/L	117	78 — 156	13
	BS-330E ⁹	nmol/L	45.5	30.3 — 60.7	5.1	BS-330E ⁹	nmol/L	115	77 — 153	13
	BS-360E ¹⁰	nmol/L	43.9	29.3 — 58.5	4.9	BS-360E ¹⁰	nmol/L	114	76 — 152	13
	BS-380 ¹¹	nmol/L	45.0	30.0 — 60.0	5.0	BS-380 ¹¹	nmol/L	116	77 — 155	13
	BS-400 ¹²	nmol/L	46.8	31.2 — 62.4	5.2	BS-400 ¹²	nmol/L	119	79 — 159	13
	BS-430 ¹³	nmol/L	46.3	30.9 — 61.7	5.1	BS-430 ¹³	nmol/L	116	77 — 155	13
	BS-480 ¹⁴	nmol/L	47.5	31.7 — 63.3	5.3	BS-480 ¹⁴	nmol/L	120	80 — 160	13
	BS-600 ¹⁵	nmol/L	49.1	32.7 — 65.5	5.5	BS-600 ¹⁵	nmol/L	121	81 — 161	13
	BS-800 ¹⁶	nmol/L	47.2	31.5 — 62.9	5.2	BS-800 ¹⁶	nmol/L	117	78 — 156	13
	BS-2000 ¹⁷	nmol/L	45.4	30.3 — 60.5	5.0	BS-2000 ¹⁷	nmol/L	115	77 — 153	13
	BS-2800M ¹⁸	nmol/L	44.3	29.5 — 59.1	4.9	BS-2800M ¹⁸	nmol/L	116	77 — 155	13

English	Abbreviated name	Lp(a) II Control	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	Lp(a) II Контроль	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Controle da Lp(a) II	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	Control Lp(a) II	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	Controllo Lp(a) II	modelli	Unità	Valori di dosaggio	Intervallo(valore di concentrazione±3DS)
Türkçe	kısaltılmış ad	Lp(a) II Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

Lp(a) II Control

English

Lipoprotein (a) Control

Lp(a) II Control

mindray

Русский	Липопротеин (а) Контроль
Português	Controle da Lipoproteína (a)
Español	Control de Lipoproteína (a)
Italiano	Controllo della Lipoproteina (a)
Türkçe	Lipoprotein (a) Kontrolü