

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

1.BS-120: BS-120, BS-130;

2.BS-180: BS-180, BS-190;

3.BS-200: BS-200, BS-220;

4.BS-200E: BS-200E, BS-220E;

5.BS-230: BS-230, BS-240;

6.BS-240E: BS240E, BS240Pro;

7.BS-300: BS-300, BS-320;

8.BS-330: BS-330, BS-350;

9.BS-330E:BS-330E(Serial Number starts with "XQ-"),

BS-350E(Serial Number starts with "XS-");

10.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00");

11.BS-380: BS-380, BS-390;

12.BS-400: BS-400, BS-420;

13.BS-430: BS-430, BS-450, BS-460;

14.BS-480: BS-480, BS-490;

15.BS-600: BS-600, BS-620;

16.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

17.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

18.BS-2800M:BS-2600M, BS-2800M.

LOT : 061122002

⌚ : 2023-08-25

Abbreviated name	Lp(a) II Control (L)					Lp(a) II Control (H)						
	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD
Lp(a) II	BS-120 ¹	nmol/L	41.6	27.7	— 55.5	4.6	BS-120 ¹	nmol/L	113	75	— 151	13
	BS-180 ²	nmol/L	41.6	27.7	— 55.5	4.6	BS-180 ²	nmol/L	113	75	— 151	13
	BS-200 ³	nmol/L	44.0	29.3	— 58.7	4.9	BS-200 ³	nmol/L	112	75	— 149	12
	BS-200E ⁴	nmol/L	42.4	28.3	— 56.5	4.7	BS-200E ⁴	nmol/L	111	74	— 148	12
	BS-230 ⁵	nmol/L	42.8	28.5	— 57.1	4.8	BS-230 ⁵	nmol/L	113	75	— 151	13
	BS-240E ⁶	nmol/L	42.7	28.5	— 56.9	4.7	BS-240E ⁶	nmol/L	112	75	— 149	12
	BS-300 ⁷	nmol/L	50.7	33.8	— 67.6	5.6	BS-300 ⁷	nmol/L	118	79	— 157	13
	BS-330 ⁸	nmol/L	44.0	29.3	— 58.7	4.9	BS-330 ⁸	nmol/L	112	75	— 149	12
	BS-330E ⁹	nmol/L	42.4	28.3	— 56.5	4.7	BS-330E ⁹	nmol/L	111	74	— 148	12
	BS-360E ¹⁰	nmol/L	40.8	27.2	— 54.4	4.5	BS-360E ¹⁰	nmol/L	112	75	— 149	12
	BS-380 ¹¹	nmol/L	42.0	28.0	— 56.0	4.7	BS-380 ¹¹	nmol/L	118	79	— 157	13
	BS-400 ¹²	nmol/L	47.6	31.7	— 63.5	5.3	BS-400 ¹²	nmol/L	126	84	— 168	14
	BS-430 ¹³	nmol/L	41.7	27.8	— 55.6	4.6	BS-430 ¹³	nmol/L	112	75	— 149	12
	BS-480 ¹⁴	nmol/L	41.7	27.8	— 55.6	4.6	BS-480 ¹⁴	nmol/L	112	75	— 149	12
	BS-600 ¹⁵	nmol/L	41.7	27.8	— 55.6	4.6	BS-600 ¹⁵	nmol/L	112	75	— 149	12
	BS-800 ¹⁶	nmol/L	41.7	27.8	— 55.6	4.6	BS-800 ¹⁶	nmol/L	112	75	— 149	12
	BS-2000 ¹⁷	nmol/L	41.4	27.6	— 55.2	4.6	BS-2000 ¹⁷	nmol/L	115	77	— 153	13
	BS-2800M ¹⁸	nmol/L	42.4	28.3	— 56.5	4.7	BS-2800M ¹⁸	nmol/L	113	75	— 151	13

English	Abbreviated name	Lp(a) II Control	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)
Русский	сокращенное наименование	Lp(a) II Контроль	модель	Прибор	Результат анализа	Диапазон(результат анализа \pm 3CO)
Português	Nome abreviado	Controle da Lp(a) II	Modelo	Unidade	Valores da análise	Faixa(Valores da análise \pm 3SD)
Español	nombre abreviado	Control Lp(a) II	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo \pm 3SD)
Italiano	abbreviazione	Controllo Lp(a) II	modelli	Unità	Valori di dosaggio	Intervallo(valore di concentrazione \pm 3SDS)
Türkçe	kısaltılmış ad	Lp(a) II Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri \pm 3SD)

Lp(a) II Control

English	Lipoprotein (a) Control
Русский	Липопротеин (а) Контроль
Português	Controle da Lipoproteína (a)
Español	Control de Lipoproteína (a)
Italiano	Controllo della Lipoproteina (a)
Türkçe	Lipoprotein (a) Kontrolü