

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

1.BS-120: BS-120, BS-130;

2.BS-180: BS-180, BS-190;

3.BS-200: BS-200, BS-220;

4.BS-200E: BS-200E, BS-220E;

5.BS-230: BS-230, BS-240;

6.BS-240E: BS240E, BS240Pro;

7.BS-300: BS-300, BS-320;

8.BS-330: BS-330, BS-350;

9.BS-330E:BS-330E(Serial Number starts with "XQ-"),

BS-350E(Serial Number starts with "XS-");

10.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00");

11.BS-380: BS-380, BS-390;

12.BS-400: BS-400, BS-420;

13.BS-430: BS-430, BS-450, BS-460;

14.BS-480: BS-480, BS-490;

15.BS-600: BS-600, BS-620;

16.BS-600M: BS-600M, BS-620M;

17.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

18.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

19.BS-2800M:BS-2600M, BS-2800M.

LOT : 061123003

📅 : 2024-10-06

Abbreviated name	Lp(a) II Control (L)					Lp(a) II Control (H)				
	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)	1 SD	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)	1 SD
Lp(a) II	BS-120 ¹	nmol/L	34.4	22.9 — 45.9	3.8	BS-120 ¹	nmol/L	121	81 — 161	13
	BS-180 ²	nmol/L	34.4	22.9 — 45.9	3.8	BS-180 ²	nmol/L	121	81 — 161	13
	BS-200 ³	nmol/L	44.3	29.5 — 59.1	4.9	BS-200 ³	nmol/L	127	85 — 169	14
	BS-200E ⁴	nmol/L	41.7	27.8 — 55.6	4.6	BS-200E ⁴	nmol/L	124	83 — 165	14
	BS-230 ⁵	nmol/L	45.0	30.0 — 60.0	5.0	BS-230 ⁵	nmol/L	128	85 — 171	14
	BS-240E ⁶	nmol/L	42.4	28.3 — 56.5	4.7	BS-240E ⁶	nmol/L	124	83 — 165	14
	BS-300 ⁷	nmol/L	51.2	34.2 — 68.2	5.7	BS-300 ⁷	nmol/L	131	87 — 175	15
	BS-330 ⁸	nmol/L	44.3	29.5 — 59.1	4.9	BS-330 ⁸	nmol/L	127	85 — 169	14
	BS-330E ⁹	nmol/L	41.7	27.8 — 55.6	4.6	BS-330E ⁹	nmol/L	124	83 — 165	14
	BS-360E ¹⁰	nmol/L	41.9	27.9 — 55.9	4.7	BS-360E ¹⁰	nmol/L	123	82 — 164	14
	BS-380 ¹¹	nmol/L	42.6	28.4 — 56.8	4.7	BS-380 ¹¹	nmol/L	124	83 — 165	14
	BS-400 ¹²	nmol/L	44.2	29.5 — 58.9	4.9	BS-400 ¹²	nmol/L	128	85 — 171	14
	BS-430 ¹³	nmol/L	42.2	28.1 — 56.3	4.7	BS-430 ¹³	nmol/L	125	83 — 167	14
	BS-480 ¹⁴	nmol/L	42.2	28.1 — 56.3	4.7	BS-480 ¹⁴	nmol/L	125	83 — 167	14
	BS-600 ¹⁵	nmol/L	42.2	28.1 — 56.3	4.7	BS-600 ¹⁵	nmol/L	125	83 — 167	14
	BS-600M ¹⁶	nmol/L	45.0	30.0 — 60.0	5.0	BS-600M ¹⁶	nmol/L	129	86 — 172	14
	BS-800 ¹⁷	nmol/L	42.2	28.1 — 56.3	4.7	BS-800 ¹⁷	nmol/L	125	83 — 167	14
	BS-2000 ¹⁸	nmol/L	45.8	30.5 — 61.1	5.1	BS-2000 ¹⁸	nmol/L	129	86 — 172	14
	BS-2800M ¹⁹	nmol/L	43.4	28.9 — 57.9	4.8	BS-2800M ¹⁹	nmol/L	128	85 — 171	14

English	Abbreviated name	Lp(a) II Control	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)
Русский	сокращенное наименование	Lp(a) II Контроль	модель	Прибор	Результат анализа	Диапазон(результат анализа \pm 3CO)
Português	Nome abreviado	Controle da Lp(a) II	Modelo	Unidade	Valores da análise	Faixa(Valores da análise \pm 3SD)
Español	nombre abreviado	Control Lp(a) II	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo \pm 3SD)
Italiano	abbreviazione	Controllo Lp(a) II	modelli	Unità	Valori di dosaggio	Intervallo(valore di concentrazione \pm 3SD)
Türkçe	kısaltılmış ad	Lp(a) II Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri \pm 3SD)

Lp(a) II Control	
English	Lipoprotein (a) Control
Русский	Липопротеин (а) Контроль
Português	Controle da Lipoproteína (a)
Español	Control de Lipoproteína (a)
Italiano	Controllo della Lipoproteina (a)
Türkçe	Lipoprotein (a) Kontrolü