

## Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

1.BS-120: BS-120, BS-130;

2.BS-180: BS-180, BS-190;

3.BS-200: BS-200, BS-220;

4.BS-200E: BS-200E, BS-220E;

5.BS-230: BS-230, BS-240;

6.BS-240E: BS240E, BS240Pro;

7.BS-300: BS-300, BS-320;

8.BS-330: BS-330, BS-350;

9.BS-330E:BS-330E(Serial Number starts with "XQ-"),

BS-350E(Serial Number starts with "XS-");

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

10.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00");

11.BS-380: BS-380, BS-390;

12.BS-400: BS-400, BS-420;

13.BS-430: BS-430, BS-450, BS-460;

14.BS-480: BS-480, BS-490;

15.BS-600: BS-600, BS-620;

16.BS-600M: BS-600M, BS-620M;

17.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

18.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

19.BS-2800M:BS-2600M, BS-2800M.

**LOT** : 061123006

**EXP** : 2025-04-13

Abbreviated name	Lp(a) II Control (L)					Lp(a) II Control (H)				
	Model	Unit	Assay Value	Range (Assay Value $\pm$ 3SD)	1 SD	Model	Unit	Assay Value	Range (Assay Value $\pm$ 3SD)	1 SD
Lp(a) II	BS-120 <sup>1</sup>	nmol/L	37.1	24.8 — 49.4	4.1	BS-120 <sup>1</sup>	nmol/L	109	73 — 145	12
	BS-180 <sup>2</sup>	nmol/L	37.1	24.8 — 49.4	4.1	BS-180 <sup>2</sup>	nmol/L	109	73 — 145	12
	BS-200 <sup>3</sup>	nmol/L	39.6	26.4 — 52.8	4.4	BS-200 <sup>3</sup>	nmol/L	118	79 — 157	13
	BS-200E <sup>4</sup>	nmol/L	40.3	26.8 — 53.8	4.5	BS-200E <sup>4</sup>	nmol/L	116	77 — 155	13
	BS-230 <sup>5</sup>	nmol/L	43.0	28.6 — 57.4	4.8	BS-230 <sup>5</sup>	nmol/L	121	82 — 160	13
	BS-240E <sup>6</sup>	nmol/L	40.8	27.3 — 54.3	4.5	BS-240E <sup>6</sup>	nmol/L	118	79 — 157	13
	BS-300 <sup>7</sup>	nmol/L	46.4	30.8 — 62.0	5.2	BS-300 <sup>7</sup>	nmol/L	120	81 — 159	13
	BS-330 <sup>8</sup>	nmol/L	39.6	26.4 — 52.8	4.4	BS-330 <sup>8</sup>	nmol/L	118	79 — 157	13
	BS-330E <sup>9</sup>	nmol/L	40.3	26.8 — 53.8	4.5	BS-330E <sup>9</sup>	nmol/L	116	77 — 155	13
	BS-360E <sup>10</sup>	nmol/L	38.4	25.5 — 51.3	4.3	BS-360E <sup>10</sup>	nmol/L	113	74 — 152	13
	BS-380 <sup>11</sup>	nmol/L	41.0	27.2 — 54.8	4.6	BS-380 <sup>11</sup>	nmol/L	117	78 — 156	13
	BS-400 <sup>12</sup>	nmol/L	40.5	27.0 — 54.0	4.5	BS-400 <sup>12</sup>	nmol/L	118	79 — 157	13
	BS-430 <sup>13</sup>	nmol/L	40.9	27.4 — 54.4	4.5	BS-430 <sup>13</sup>	nmol/L	120	81 — 159	13
	BS-480 <sup>14</sup>	nmol/L	40.9	27.4 — 54.4	4.5	BS-480 <sup>14</sup>	nmol/L	120	81 — 159	13
	BS-600 <sup>15</sup>	nmol/L	40.9	27.4 — 54.4	4.5	BS-600 <sup>15</sup>	nmol/L	120	81 — 159	13
	BS-600M <sup>16</sup>	nmol/L	41.0	27.2 — 54.8	4.6	BS-600M <sup>16</sup>	nmol/L	122	80 — 164	14
	BS-800 <sup>17</sup>	nmol/L	40.9	27.4 — 54.4	4.5	BS-800 <sup>17</sup>	nmol/L	120	81 — 159	13
	BS-2000 <sup>18</sup>	nmol/L	41.1	27.3 — 54.9	4.6	BS-2000 <sup>18</sup>	nmol/L	117	78 — 156	13
	BS-2800M <sup>19</sup>	nmol/L	43.1	28.7 — 57.5	4.8	BS-2800M <sup>19</sup>	nmol/L	121	82 — 160	13

English	Abbreviated name	Lp(a) II Control	Model	Unit	Assay Value	Range (Assay Value $\pm$ 3SD)
Русский	сокращенное наименование	Lp(a) II Контроль	модель	Прибор	Результат анализа	Диапазон(результат анализа $\pm$ 3CO)
Português	Nome abreviado	Controle da Lp(a) II	Modelo	Unidade	Valores da análise	Faixa(Valores da análise $\pm$ 3SD)
Español	nombre abreviado	Control Lp(a) II	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo $\pm$ 3SD)
Italiano	abbreviazione	Controllo Lp(a) II	modelli	Unità	Valori di dosaggio	Intervallo(valore di concentrazione $\pm$ 3DS)

# Lp(a) II Control

**mindray**

Türkçe	kısaltılmış ad	Lp(a) II Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)
<b>Lp(a) II Control</b>						
<b>English</b>				Lipoprotein (a) Control		
<b>Русский</b>				Липопротеин (а) Контроль		
<b>Português</b>				Controle da Lipoproteína (a)		
<b>Español</b>				Control de Lipoproteína (a)		
<b>Italiano</b>				Controllo della Lipoproteina (a)		
<b>Türkçe</b>				Lipoprotein (a) Kontrolü		