

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

1.BS-120: BS-120, BS-130;

2.BS-180: BS-180, BS-190;

3.BS-200: BS-200, BS-220;

4.BS-200E: BS-200E, BS-220E;

5.BS-230: BS-230, BS-240;

6.BS-240E: BS240E, BS240Pro;

7.BS-300: BS-300, BS-320;

8.BS-330: BS-330, BS-350;

9.BS-330E:BS-330E(Serial Number starts with "XQ-"),

BS-350E(Serial Number starts with "XS-");

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

10.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S,
BS-330E(V35.00)(Serial Number starts with "W8-"
and software version starts with "35.00"),
BS-350E(V35.00) (Serial Number starts with"W9-"
and software version starts with "35.00");

11.BS-380: BS-380, BS-390;

12.BS-400: BS-400, BS-420;

13.BS-430: BS-430, BS-450, BS-460;

14.BS-480: BS-480, BS-490;

15.BS-600: BS-600, BS-620;

16.BS-600M: BS-600M, BS-620M;

17.BS-800: BS-800, BS-820, BS-800M, BS-820M,
BS-1800, BS-1800plus;

18.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

19.BS-2800M:BS-2600M, BS-2800M.

LOT : 061123007

EXP : 2025-05-15

Abbreviated name	Lp(a) II Control (L)					Lp(a) II Control (H)				
	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)	1 SD	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)	1 SD
Lp(a) II	BS-120 ¹	nmol/L	37.1	24.8 — 49.4	4.1	BS-120 ¹	nmol/L	109	73 — 145	12
	BS-180 ²	nmol/L	37.1	24.8 — 49.4	4.1	BS-180 ²	nmol/L	109	73 — 145	12
	BS-200 ³	nmol/L	39.6	26.4 — 52.8	4.4	BS-200 ³	nmol/L	118	79 — 157	13
	BS-200E ⁴	nmol/L	40.3	26.8 — 53.8	4.5	BS-200E ⁴	nmol/L	116	77 — 155	13
	BS-230 ⁵	nmol/L	43.0	28.6 — 57.4	4.8	BS-230 ⁵	nmol/L	121	82 — 160	13
	BS-240E ⁶	nmol/L	40.8	27.3 — 54.3	4.5	BS-240E ⁶	nmol/L	118	79 — 157	13
	BS-300 ⁷	nmol/L	46.4	30.8 — 62.0	5.2	BS-300 ⁷	nmol/L	120	81 — 159	13
	BS-330 ⁸	nmol/L	39.6	26.4 — 52.8	4.4	BS-330 ⁸	nmol/L	118	79 — 157	13
	BS-330E ⁹	nmol/L	40.3	26.8 — 53.8	4.5	BS-330E ⁹	nmol/L	116	77 — 155	13
	BS-360E ¹⁰	nmol/L	38.4	25.5 — 51.3	4.3	BS-360E ¹⁰	nmol/L	113	74 — 152	13
	BS-380 ¹¹	nmol/L	41.0	27.2 — 54.8	4.6	BS-380 ¹¹	nmol/L	117	78 — 156	13
	BS-400 ¹²	nmol/L	40.5	27.0 — 54.0	4.5	BS-400 ¹²	nmol/L	118	79 — 157	13
	BS-430 ¹³	nmol/L	40.9	27.4 — 54.4	4.5	BS-430 ¹³	nmol/L	120	81 — 159	13
	BS-480 ¹⁴	nmol/L	40.9	27.4 — 54.4	4.5	BS-480 ¹⁴	nmol/L	120	81 — 159	13
	BS-600 ¹⁵	nmol/L	40.9	27.4 — 54.4	4.5	BS-600 ¹⁵	nmol/L	120	81 — 159	13
	BS-600M ¹⁶	nmol/L	41.0	27.2 — 54.8	4.6	BS-600M ¹⁶	nmol/L	122	80 — 164	14
	BS-800 ¹⁷	nmol/L	40.9	27.4 — 54.4	4.5	BS-800 ¹⁷	nmol/L	120	81 — 159	13
	BS-2000 ¹⁸	nmol/L	41.1	27.3 — 54.9	4.6	BS-2000 ¹⁸	nmol/L	117	78 — 156	13
	BS-2800M ¹⁹	nmol/L	43.1	28.7 — 57.5	4.8	BS-2800M ¹⁹	nmol/L	121	82 — 160	13

English	Abbreviated name	Lp(a) II Control	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)
Русский	сокращенное наименование	Lp(a) II Контроль	модель	Прибор	Результат анализа	Диапазон(результат анализа \pm 3CO)
Português	Nome abreviado	Controle da Lp(a) II	Modelo	Unidade	Valores da análise	Faixa(Valores da análise \pm 3SD)
Español	nombre abreviado	Control Lp(a) II	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo \pm 3SD)
Italiano	abbreviazione	Controllo Lp(a) II	modelli	Unità	Valori di dosaggio	Intervallo(valore di concentrazione \pm 3DS)

Lp(a) II Control

mindray

Türkçe	kısaltılmış ad	Lp(a) II Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)
Lp(a) II Control						
English				Lipoprotein (a) Control		
Русский				Липопротеин (а) Контроль		
Português				Controle da Lipoproteína (a)		
Español				Control de Lipoproteína (a)		
Italiano				Controllo della Lipoproteina (a)		
Türkçe				Lipoprotein (a) Kontrolü		