

Value sheet of Mindray BS Measurement System

Русский : Таблица результатов для системы BS компании Mindray

Português : Planilha de valores do Sistema de Medição BS da Mindray

Español : Hoja de valores del sistema de medición Mindray BS

Italiano : Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe : Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский : Данные совпадают во всех группах.

Español : la datos de cada grupo es la misma.

Türkçe : her grubun veri aynıdır.

Português : A dados de cada grupo é a mesma.

Italiano : la dati di ogni gruppo è la stessa.

1.**BS-120**: BS-120, BS-130;

2.**BS-180**: BS-180, BS-190;

3.**BS-200**: BS-200, BS-220;

4.**BS-200E**: BS-200E, BS-220E;

5.**BS-230**: BS-230, BS-240, BS-280;

6.**BS-240E**: BS-240E, BS-240Pro;

7.**BS-300**: BS-300, BS-320;

8.**BS-330**: BS-330, BS-350;

9.**BS-330E**:**BS-330E**(Serial Number starts with "XQ-"),

BS-350E(Serial Number starts with "XS-");

10.**BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S, **BS-330E(V35.00)**(Serial Number starts with "W8-" and software version starts with "35.00"), **BS-350E(V35.00)** (Serial Number starts with "W9-" and software version starts with "35.00");

11.**BS-380**: BS-380, BS-390;

12.**BS-400**: BS-400, BS-420;

13.**BS-430**: BS-410, BS-430, BS-450, BS-460, BS-470;

14.**BS-480**: BS-480, BS-490;

15.**BS-600**: BS-600, BS-620;

16.**BS-600M**: BS-600M, BS-620M;

17.**BS-800**: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

18.**BS-830**: BS-830, BS-830S, BS-840, BS-850, BS-860, BS-870;

19.**BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M.

20.**BS-2800M**: BS-2600M, BS-2800M

LOT : 061523005

: 2024-10-23

Abbreviated name	ADA Control(L)					ADA Control(H)								
	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD		
ADA	BS-120 ¹	U/L	26.2	21.4	—	31.0	1.6	BS-120 ¹	U/L	116	95	—	137	7
	BS-180 ²	U/L	26.2	21.4	—	31.0	1.6	BS-180 ²	U/L	116	95	—	137	7
	BS-200 ³	U/L	25.9	21.1	—	30.7	1.6	BS-200 ³	U/L	119	98	—	140	7
	BS-200E ⁴	U/L	26.2	21.4	—	31.0	1.6	BS-200E ⁴	U/L	119	98	—	140	7
	BS-230 ⁵	U/L	27.1	22.3	—	31.9	1.6	BS-230 ⁵	U/L	118	97	—	139	7
	BS-240E ⁶	U/L	25.9	21.1	—	30.7	1.6	BS-240E ⁶	U/L	117	96	—	138	7
	BS-300 ⁷	U/L	26.0	21.2	—	30.8	1.6	BS-300 ⁷	U/L	117	96	—	138	7
	BS-330 ⁸	U/L	25.9	21.1	—	30.7	1.6	BS-330 ⁸	U/L	119	98	—	140	7
	BS-330E ⁹	U/L	26.2	21.4	—	31.0	1.6	BS-330E ⁹	U/L	119	98	—	140	7
	BS-360E ¹⁰	U/L	25.8	21.3	—	30.3	1.5	BS-360E ¹⁰	U/L	117	96	—	138	7
	BS-380 ¹¹	U/L	25.3	20.8	—	29.8	1.5	BS-380 ¹¹	U/L	116	95	—	137	7
	BS-400 ¹²	U/L	26.0	21.2	—	30.8	1.6	BS-400 ¹²	U/L	118	97	—	139	7
	BS-430 ¹³	U/L	25.8	21.3	—	30.3	1.5	BS-430 ¹³	U/L	119	98	—	140	7
	BS-480 ¹⁴	U/L	25.9	21.1	—	30.7	1.6	BS-480 ¹⁴	U/L	118	97	—	139	7
	BS-600 ¹⁵	U/L	25.9	21.1	—	30.7	1.6	BS-600 ¹⁵	U/L	116	95	—	137	7
	BS-600M ¹⁶	U/L	26.5	21.7	—	31.3	1.6	BS-600M ¹⁶	U/L	121	100	—	142	7
	BS-800 ¹⁷	U/L	25.9	21.1	—	30.7	1.6	BS-800 ¹⁷	U/L	118	97	—	139	7
	BS-830 ¹⁸	U/L	25.5	21.0	—	30.0	1.5	BS-830 ¹⁸	U/L	118	97	—	139	7
	BS-2000 ¹⁹	U/L	26.4	21.6	—	31.2	1.6	BS-2000 ¹⁹	U/L	120	99	—	141	7
	BS-2800M ²⁰	U/L	26.5	21.7	—	31.3	1.6	BS-2800M ²⁰	U/L	121	100	—	142	7

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа \pm 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise \pm 3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo \pm 3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo(valore diconcentrazione \pm 3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri \pm 3SD)

ADA	
English	Adenosine Deaminase
Русский	является аденозиндеаминаза
Português	Adenosina deaminase
Español	adenosina desaminasa
Italiano	adenosina deaminasi
Türkçe	Adenozin Deaminazdır