

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

- | | |
|---|--|
| <p>1.BS-120: BS-120, BS-130;</p> <p>2.BS-180: BS-180, BS-190;</p> <p>3.BS-200: BS-200, BS-220;</p> <p>4.BS-200E: BS-200E, BS-220E;</p> <p>5.BS-230: BS-230, BS-240, BS-280;</p> <p>6.BS-240E: BS-240E, BS-240Pro;</p> <p>7.BS-300: BS-300, BS-320;</p> <p>8.BS-330: BS-330, BS-350;</p> <p>9.BS-330E:BS-330E(Serial Number starts with "XQ-"),
BS-350E(Serial Number starts with "XS-");</p> | <p>10.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S,
BS-330E(V35.00)(Serial Number starts with "W8-"
and software version starts with "35.00"),
BS-350E(V35.00) (Serial Number starts with "W9-"
and software version starts with "35.00");</p> <p>11.BS-380: BS-380, BS-390;</p> <p>12.BS-400: BS-400, BS-420;</p> <p>13.BS-430: BS-410, BS-430, BS-450, BS-460, BS-470;</p> <p>14.BS-480: BS-480, BS-490;</p> <p>15.BS-600: BS-600, BS-620;</p> <p>16.BS-600M: BS-600M, BS-620M;</p> <p>17.BS-800: BS-800, BS-820, BS-800M, BS-820M,
BS-1800, BS-1800plus;</p> <p>18.BS-830: BS-830, BS-830S, BS-840, BS-850,
BS-860, BS-870;</p> <p>19.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M.</p> <p>20.BS-2800M: BS-2600M, BS-2800M</p> |
|---|--|

LOT : 061523007

EXP : 2024-12-06

Abbreviated name	ADA Control(L)					ADA Control(H)								
	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD		
ADA	BS-120 ¹	U/L	25.7	21.2	—	30.2	1.5	BS-120 ¹	U/L	124	103	—	145	7
	BS-180 ²	U/L	25.7	21.2	—	30.2	1.5	BS-180 ²	U/L	124	103	—	145	7
	BS-200 ³	U/L	26.1	21.3	—	30.9	1.6	BS-200 ³	U/L	127	103	—	151	8
	BS-200E ⁴	U/L	25.5	21.0	—	30.0	1.5	BS-200E ⁴	U/L	126	102	—	150	8
	BS-230 ⁵	U/L	26.0	21.2	—	30.8	1.6	BS-230 ⁵	U/L	125	101	—	149	8
	BS-240E ⁶	U/L	25.7	21.2	—	30.2	1.5	BS-240E ⁶	U/L	126	102	—	150	8
	BS-300 ⁷	U/L	26.3	21.5	—	31.1	1.6	BS-300 ⁷	U/L	126	102	—	150	8
	BS-330 ⁸	U/L	26.1	21.3	—	30.9	1.6	BS-330 ⁸	U/L	127	103	—	151	8
	BS-330E ⁹	U/L	25.5	21.0	—	30.0	1.5	BS-330E ⁹	U/L	126	102	—	150	8
	BS-360E ¹⁰	U/L	25.5	21.0	—	30.0	1.5	BS-360E ¹⁰	U/L	126	102	—	150	8
	BS-380 ¹¹	U/L	25.4	20.9	—	29.9	1.5	BS-380 ¹¹	U/L	125	101	—	149	8
	BS-400 ¹²	U/L	25.6	21.1	—	30.1	1.5	BS-400 ¹²	U/L	126	102	—	150	8
	BS-430 ¹³	U/L	25.8	21.3	—	30.3	1.5	BS-430 ¹³	U/L	127	103	—	151	8
	BS-480 ¹⁴	U/L	25.6	21.1	—	30.1	1.5	BS-480 ¹⁴	U/L	127	103	—	151	8
	BS-600 ¹⁵	U/L	25.6	21.1	—	30.1	1.5	BS-600 ¹⁵	U/L	126	102	—	150	8
	BS-600M ¹⁶	U/L	26.2	21.4	—	31.0	1.6	BS-600M ¹⁶	U/L	129	105	—	153	8
	BS-800 ¹⁷	U/L	25.6	21.1	—	30.1	1.5	BS-800 ¹⁷	U/L	126	102	—	150	8
	BS-830 ¹⁸	U/L	25.9	21.1	—	30.7	1.6	BS-830 ¹⁸	U/L	126	102	—	150	8
	BS-2000 ¹⁹	U/L	26.0	21.2	—	30.8	1.6	BS-2000 ¹⁹	U/L	128	104	—	152	8
	BS-2800M ²⁰	U/L	26.2	21.4	—	31.0	1.6	BS-2800M ²⁰	U/L	129	105	—	153	8

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа \pm 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise \pm 3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo \pm 3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo(valore diconcentrazione \pm 3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri \pm 3SD)

ADA	
English	Adenosine Deaminase
Русский	является аденозиндезаминаза
Português	Adenosina deaminase
Español	adenosina desaminasa
Italiano	adenosina deamminasi
Türkçe	Adenozin Deaminazdır