

# Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

- 1.BS-120: BS-120, BS-130, BS-180, BS-190;
- 2.BS-200: BS-200, BS-220, BS-330, BS-350;
- 3.BS-200E: BS-200E, BS-220E,;
- 4.BS-240: BS-230, BS-240;
- 5.BS-240E: BS240E, BS240Pro;
- 6.BS-300: BS-300, BS-320;
- 7.BS-330E:BS-330E(Serial Number starts with "XQ-"), BS-350E(Serial Number starts with "XS-")
- 8.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00");
- 9.BS-380: BS-380, BS-390;

- 10.BS-400: BS-400, BS-420;
- 11.BS-430: BS-430, BS-450,BS-460;
- 12.BS-480: BS-480, BS-490;
- 13.BS-600: BS-600, BS-620;
- 14.BS-600M: BS-600M, BS-620M;
- 15.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;
- 16.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;
- 17.BS-2800M:BS-2600M, BS-2800M;
- 18.Applicable models of the chemistry shall be subject to the parameter list and instructions.

**LOT** : 061524001

**🕒**: 2025-04-15

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

**ADA**

English	Adenosine Deaminase
Русский	является аденозиндеаминаза
Português	Adenosina deaminase
Español	adenosina desaminasa
Italiano	adenosina deamminasi
Türkçe	Adenozin Deaminazdır

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
ADA	U/L	BS-120 <sup>1</sup>	27.8	22.7	—	32.9	1.7	BS-120 <sup>1</sup>	123	102	—	144	7
		BS-200 <sup>2</sup>	27.8	22.7	—	32.9	1.7	BS-200 <sup>2</sup>	128	104	—	152	8
		BS-200E <sup>3</sup>	27.0	22.2	—	31.8	1.6	BS-200E <sup>3</sup>	125	101	—	149	8
		BS-240 <sup>4</sup>	27.9	22.8	—	33.0	1.7	BS-240 <sup>4</sup>	126	102	—	150	8
		BS-240E <sup>5</sup>	27.0	22.2	—	31.8	1.6	BS-240E <sup>5</sup>	125	101	—	149	8
		BS-300 <sup>6</sup>	27.8	22.7	—	32.9	1.7	BS-300 <sup>6</sup>	126	102	—	150	8
		BS-330E <sup>7</sup>	27.0	22.2	—	31.8	1.6	BS-330E <sup>7</sup>	125	101	—	149	8
		BS-360E <sup>8</sup>	27.2	22.4	—	32.0	1.6	BS-360E <sup>8</sup>	126	102	—	150	8
		BS-380 <sup>9</sup>	26.9	22.1	—	31.7	1.6	BS-380 <sup>9</sup>	125	101	—	149	8
		BS-400 <sup>10</sup>	27.2	22.4	—	32.0	1.6	BS-400 <sup>10</sup>	126	102	—	150	8
		BS-430 <sup>11</sup>	27.3	22.5	—	32.1	1.6	BS-430 <sup>11</sup>	127	103	—	151	8
		BS-480 <sup>12</sup>	27.4	22.6	—	32.2	1.6	BS-480 <sup>12</sup>	128	104	—	152	8
		BS-600 <sup>13</sup>	27.2	22.4	—	32.0	1.6	BS-600 <sup>13</sup>	126	102	—	150	8
		BS-600M <sup>14</sup>	27.5	22.4	—	32.6	1.7	BS-600M <sup>14</sup>	128	104	—	152	8
		BS-800 <sup>15</sup>	27.4	22.6	—	32.2	1.6	BS-800 <sup>15</sup>	128	104	—	152	8

---

<b>BS-2000<sup>16</sup></b>	27.5	22.4	—	32.6	1.7	<b>BS-2000<sup>16</sup></b>	127	103	—	151	8
<b>BS-2800M<sup>17</sup></b>	27.5	22.4	—	32.6	1.7	<b>BS-2800M<sup>17</sup></b>	128	104	—	152	8


---

# mindray

## ADA Control

For use on: BS-2000, BS-2800M,  
BS-600M, BS-800

**LOT** 061524001

 2025-04-15

